

Minncare® Cold Sterilant

Application Brief: Sporicide Concentrate For Surface Disinfection in Pharmaceutical Cleanroom Environments

Minncare Cold Sterilant is an effective disinfectant and sterilant for hard surfaces commonly found in pharmaceutical cleanrooms. Surfaces made from non-porous materials such as plastics, stainless steel or glass can be disinfected using Minncare Cold Sterilant solutions. Examples include counter tops, walls, floors, etc.

Minncare Cold Sterilant Effectiveness

Minncare Cold Sterilant is a powerful disinfectant containing a proprietary mixture of 4.5% peracetic acid and 22.0% hydrogen peroxide. It is effective against a broad spectrum of organisms, including bacterial spores, fungal spores, bacteria, mycobacteria, yeast, molds and viruses. Minncare Cold Sterilant is not a systemic bactericide. Rather it is a powerful oxidizer that attacks and breaks down the organic components of microorganisms, not just killing them, but destroying their structure. This makes Minncare Cold Sterilant one of the fastest acting and most effective disinfectants available. Minncare Cold Sterilant is registered with the US EPA as a disinfectant and sterilant. It has also been successfully tested according to European Norms for bactericidal, mycobactericidal, virucidal and sporicidal activity.*



Procedure for use as a Germicidal Spray

1. Using water or mechanical action, remove heavy soil and gross filth from the surfaces to be disinfected.
2. Dilute Minncare Cold Sterilant with purified water to a 1% solution (100X dilution).
3. Apply the solution using a plastic spray bottle.
4. Allow the solution to remain on the surface for 10 minutes and let air dry or rinse with purified water, drain off any excess water if possible, and then allow to air dry.

Procedure for use as a Cleaner/Sanitizer

1. Using water or mechanical action, remove heavy soil and gross filth from the surfaces to be disinfected.
2. Dilute Minncare Cold Sterilant with purified water to a 3% solution (32X dilution).
3. Apply the solution to surfaces using an appropriate method such as wiping onto the surface with a cloth.
4. Allow the solution to remain on the surface for 10 minutes and let air dry or rinse with purified water, drain off any excess water if possible, and then allow to air dry.

*Tested for sporicidal activity according to prEN14347 and EN13704, for bactericidal activity according to EN1040 and EN1276, mycobactericidal activity by prEN14348 and fungicidal activity by EN1275 and EN1650.

Minncare® Cold Sterilant

For Surface Disinfection in Pharmaceutical Cleanroom Environments

Storage and Disposal

Once diluted, the solution must be used within seven (7) days. Discard any unused diluted solution. Fresh solution should be prepared daily or more often if the solution becomes diluted or soiled. Always follow local regulations for the disposal of any unused Minncare Cold Sterilant solution.

Cautions

1. DO NOT USE AFTER EXPIRATION DATE.
2. Do not allow Minncare Cold Sterilant to mix with alkaline substances such as bleach (sodium hypochlorite).
3. If spilled, flush away with large quantities of water.
4. Use purified water when diluting Minncare Cold Sterilant.
5. Minncare Cold Sterilant is a single-use product not intended for reuse.
6. Fresh sanitizing solution should be prepared daily or more often if the solution becomes diluted or soiled.
7. Always follow the directions for use as shown on the Minncare Cold Sterilant label.

Other Information

Other documents are available for your information. These include the “Minncare Cold Sterilant Research Data Report” which contains information on efficacy and stability, the “Actril Cold Sterilant for Surface Disinfection” Application Note, which has information similar to that of this note for that disinfectant product and the Minncare Cold Sterilant and Actril Cold Sterilant brochures.

Contact your Mar Cor Purification representative for these documents and other information on Mar Cor Purification products.

Minncare Cold Sterilant is registered with the US Environmental Protection Agency (EPA Registration Number 52252-4).

Minncare Ordering Information

U.S. Part No.	Canadian Part No.	European Part No.	Description	Packaging
176-01-001	78399-647	n/a	Minncare Cold Sterilant 4 x 1 quart	4 x 1 quarts
176-01-002	78399-646	n/a	Minncare Cold Sterilant 4 US gallons	4 x 1 US gallons
176-01-003	78399-648	n/a	Minncare Cold Sterilant 40 US gallons	40 gallons in 55 gal. drum
n/a	n/a	176-01-125	Minncare Cold Sterilant 2 x 5 liters	2 x 5 liters
n/a	n/a	176-01-121	Minncare Cold Sterilant 2 x 10 liters	2 x 10 liters
n/a	n/a	176-01-129	Minncare Cold Sterilant 6 x 1 kg	6 x 1 kg
185-40-004	78338-000	185-40-004	Minncare Residual Test Strips	100 per vial
185-40-005	78339-000	185-40-005	Minncare 1% Test Strips	100 per vial

Minncare® is a registered trademark of Minntech Corporation, a Cantel Medical Company



Mar Cor Purification 4450 Township Line Road Skipack, PA 19474-1429 Tel: (484) 991-0220 Toll Free: (800) 346-0365 Fax: (484) 991-0230	Mar Cor Purification 14550 28th Avenue North Plymouth, MN 55447 Tel: (484) 991-0220 Toll Free: (800) 633-3080 Fax: (763) 210-3868	Mar Cor Purification 3250 Harvester Road - Unit 6 Burlington, ON L7N 3W9 Tel: (905) 639-7025 Toll Free: (800) 268-5035 Fax: (905) 639-0425	Mar Cor Purification Sourethweg 11 6422 PC Heerlen The Netherlands Tel: (+31) 45 5471 471 Fax: (+31) 45 5429 695	Mar Cor Purification 1A International Business Park, #05-01 Singapore 609933 Tel: (+65) 6227 9698 Fax: (+65) 6225 6848
--	--	---	---	---