

### USEFUL FORMULAS

<b>1 gallon</b>	=	3.785 liters
<b>1 grain per gallon</b>	=	17.1 mg/liter
<b>1 mg/liter</b>	=	1 ppm (parts per million)
<b>1 micromho/cm</b>	=	.4 mg/liter as NaCl .65 mg/liter as ion .5 mg/liter as CaCO <sub>3</sub>
<b>TDS in ppm as CaCO<sub>3</sub></b>	=	Conductivity (micromhos) x .5
<b>Conductivity (micromho/cm)</b>	=	$\frac{1}{\text{Resistivity (megohm-cm)}}$
<b>Resistivity (megohm-cm)</b>	=	$\frac{1}{\text{Conductivity (micromho/cm)}}$
<b>Concentration Factor</b>	=	Percent Recovery
<b>Percent of Recovery</b>	=	$\frac{\text{flow rate of product}}{\text{flow rate of feed}} \times 100$

$$\text{Percent of Rejection} = \left( \frac{\text{Conductivity (micromho/cm) of feed} - \text{Conductivity of Product}}{\text{Conductivity of Feed}} \right) \times 100$$

### HARDNESS EQUATION:

$$(\text{mg CaCO}_3/\text{l} = 2.497 \times \text{Ca}(\text{mg/l}) + 4.118 \times \text{mg}(\text{mg/l}) = \text{Hardness as CaCO}_3 \div 17.1 = \text{GPG}.$$

### CARBON SIZING EQUATION:

$$\text{Empty bed contact time: } V = \frac{Q (\text{EBCT})}{7.48}$$

V = Volume (cu. ft)

Q = Flow

EBCT = Empty bed contact time needed.

**NOTE:** Recommended contact time is 6 minutes to remove free chlorine.  
Recommended contact time is 10 minutes to remove chloramine.

### RO SIZING EQUATION:

Total the flow needed for hemodialysis equipment, re-use equipment and other water uses. For a recirculating loop add 1 GPM to the total flow requirements to meet the minimum standard of 1.5 ft/sec through the loop.