

MINNCARE[®] HD Disinfectant

MinnCare HD is the first chemical disinfectant to hold a FDA 510(k) cleared claim for the disinfection of water purification systems for hemodialysis. MinnCare HD is a peracetic acid solution developed specifically for use on hemodialysis reverse osmosis (RO) membranes and their associated distribution systems. Peracetic acid formulations were developed by Minntech Corporation (Medivators Inc.) in 1976 and are currently being used in a wide variety of applications. Mar Cor Purification distributes MinnCare HD Disinfectant throughout the renal community.

Protect Both You and Your System

MinnCare HD Disinfectant ensures that your dialysis water purification equipment disinfection procedures are fully validated. The disinfectant was specifically developed and tested for conditions present within water purification systems licensed for use in hemodialysis applications.

Elimination of Risk

MinnCare HD eliminates the risk of using unregulated industrial chemicals, such as chlorine bleach, in patient treatment areas. The disinfection protocol is validated, the chemistry has been proven efficacious, and accurate test strips are available to insure dosage, as well as to verify complete disinfectant rinse out of the disinfectant. The use of MinnCare HD in a dialysis clinic will help maintain the water system in optimum condition and make it easier for the water system operator to complete an effective disinfection.

Biodegradable and Ease of Disposal

MinnCare HD Disinfectant is biodegradable and will decompose into oxygen, water and acetic acid. It poses no threat to municipal waste treatment plants.* MinnCare HD Residual test strips are available that accurately verify rinse-out levels of MinnCare HD Disinfectant as low as one part per million.

510(k) Medical Device



Superior Biocidal Activity

Protect reverse osmosis membranes and their associated distribution systems against short- and long-term bacterial problems. MinnCare HD Disinfectant is a powerful oxidant that stops organism growth by oxidizing microbial cell proteins and enzyme systems and has been shown to be effective against biofilms.

MinnCare HD Disinfectant provides a 6 log reduction of bacteria population at a 100X dilution (1%) within a 36-minute soak time. Elimination of bacterial contamination will help improve product quality and contribute to increased RO membrane life.

Compatible w/ Thin Film Composite Membranes

MinnCare HD Disinfectant gently and completely cleans membranes. Original MinnCare Cold Sterilant has been successfully used at hundreds of locations on thin film composite membranes and cellulose acetate membranes. Many of the applications have included extended soak times.

*Spent MinnCare HD Disinfectant is acidic and requires normal neutralization as specified by your state and local regulations. A 1% solution of MinnCare HD Disinfectant has a pH of 3.5.

Minnicare® HD Test Strips

Medical Device

510 (k) Peracetic Acid Test Strips

Minnicare HD Test Strips provide quick results with easy-to-read indicators with the assurance of accurate results associated with a 510(k) medical device. Use the Minnicare HD 1% TS Test Strips as a pass/fail measurement for adequate concentration of sterilant after dilution. Then verify residual levels after rinse-out using Minnicare HD Residual Test Strips.

Minnicare HD Peracetic Acid Test Strips represent the leading brand of test strip technology. Designed for hemodialysis applications, these test strips combine ease-of-use with precise results.

Both the Minnicare HD 1% Indicator Test Strips and Minnicare HD Residual Test Strips may be used up to and including the last day of the month indicated by the expiration printed on the container.



Minnicare® HD 1% Indicator Test Strips

Minnicare HD 1% Test Strips verify adequate concentration of a 1% Minnicare solution during sanitization of your hemodialysis water purification system. These convenient dip-and-read strips can be used at any test port in your system to provide an accurate pass/fail measurement of Minnicare HD Disinfectant.

Minnicare® HD Residual Test Strips

Minnicare HD Residual Test Strips are easy to use and are conveniently packaged 100 test strips per container. Each strip has a “reaction zone” pad located at one end of the strip. When stored and used properly, the test strip pad can indicate residual Minnicare HD levels as low as 1 part per million (ppm). The reaction zones are set to provide six different reading levels at 100, 30, 10, 3, 1 and 0 (zero) parts per million.

Minnicare HD Ordering Information

Part No.	Description	Packaging
3029765	Minnicare HD Cold Disinfectant 4 x 1 quart	4 x 1 quarts
3029764	Minnicare HD Cold Disinfectant 4 US gallons	4 x 1 US gallons
3029795	Minnicare HD Residual Test Strips	100 per vial
3029796	Minnicare HD 1% Test Strips	100 per vial

Minnicare® is a registered trademark of Medivators Inc., a Cantel Medical Company



Mar Cor Purification
4450 Township Line Road
Skipack, PA 19474-1429
Tel: (484) 991-0220
Toll Free: (800) 346-0365
Fax: (484) 991-0230

Mar Cor Purification
14550 28th Avenue North
Plymouth, MN 55447
Tel: (484) 991-0220
Toll Free: (800) 633-3080
Fax: (763) 210-3868

Mar Cor Purification
3250 Harvester Road - Unit 6
Burlington, ON L7N 3W9
Tel: (905) 639-7025
Toll Free: (800) 268-5035
Fax: (905) 639-0425

Mar Cor Purification
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: (+31) 45 5471 471
Fax: (+31) 45 5429 695

Mar Cor Purification
1A International
Business Park, #05-01
Singapore 609933
Tel: (+65) 6227 9698
Fax: (+65) 6225 6848