

Biolab Equipment Ltd. o/a Mar Cor Purification

BioPure HX2[®] SYSTEM

System Operation Manual



Unit may not appear exactly as shown

INTENDED USE

CANADA

This device is intended to remove organic, inorganic, and microbial contaminants from water that is used to dilute dialysis concentrate to form dialysate, and to produce purified water for dialyzer reprocessing and equipment rinse and sterilization.

CAUTION: *When used as a Medical Device, federal law restricts this device to sale by and or on the order of a physician.*

UNITED STATES OF AMERICA

The BioPure HX2 Reverse Osmosis water purification system is designed to purify pre-treated water using reverse osmosis for making dialysate for hemodialysis applications. The device is intended to be a component in a complete water purification system. It must be preceded by pre-treatment devices, and may need to be followed by post-treatment devices as well to meet current AAMI and Federal (U.S.) standards. BioPure HX2 is intended for use in water rooms in a hospital, clinic, or dialysis center. The device includes an integrated heat disinfection process.

The BioPure HX2 is available two different double pass configurations that supply 5 gallons per minute (gpm) and 11 gpm. Version HX2P-05 is a double pass triple membrane RO that produces up to 5 gpm of product water. Version HX2P-11 is a double pass six membrane RO that produces up to 11 gpm of product water.

CAUTION: *When used as a Medical Device, federal law restricts this device to sale by and or on the order of a physician.*

This manual is intended to provide instructions to operate and provide basic maintenance instructions for the BioPure HX2® System.

BIOLAB SERIAL #:

DEVICE IDENTIFIER:

DATE MANUFACTURED:

MANUAL REVISION HISTORY

Revision	Author	Date	Changes
A	Peter Hang	June 19, 2013	Initial
B	Janet Parkin	July 30, 2013	Update Intended Use with US Specific Wording, 4.5.3-4.5.4 addition of manual bypass option
C	Janet Parkin	September 9, 2014	50 Hz option added, - Updates: Logo / Format / BioLab Serial # / 3.6 Define Security levels / 5.9.1 Remove test strip storage info and replace with Caution
D	Janet Parkin	July 19, 2016	Correction French Language Labels Equipment Weight Updated in sections: SPECIFICATION, 2.2 and 2.3. Update Spare Parts List, Screen Shots 4.3 rapid recovery mode 8 addition 5.1.4 High Voltage Connection Checks
E	Jeno Takacs / Janet Parkin	August 16, 2016	SPECIFICATIONS – Clarification Wet Weight 3.5.22 Temperature correction 3.1 Security Levels corrected 3.5.22 Temperature/Flow Rate correction 3.6 Footnotes added 4.1 Clarification re filters 5.15.4 Test solution given in °F and °C 8 Descriptions clarified Minor punctuation, grammar, spelling corrections

MANUAL PREPARATION

Richard Marmen – Engineering Manager
 Biolab Equipment Ltd o/a Mar Cor Purification

August 20, 2016
 Date

TABLE OF CONTENTS

INTENDED USE	2
SAFETY	8
DEFINITIONS	8
LABELS	9
SPECIFICATIONS	10
GENERAL PRECAUTIONS	12
1.0 BIOPURE HX2® SYSTEM OVERVIEW	15
2.0 INSTALLATION REQUIREMENTS	17
3.0 SYSTEM CONTROL OVERVIEW	18
3.1 SECURITY	18
3.2 ALARM MANAGEMENT	19
3.2.1 Main Control Panel Status and Alarm Indicator Function	19
3.3 REMOTE NURSE'S STATION MONITOR AND ALARM INDICATOR	20
3.4 USB DRIVE AND PRINTER (OPTIONAL)	20
3.5 THE HMI DISPLAY NAVIGATION	21
3.5.1 Valve Indicator Positions	22
3.5.2 Screen Layout.....	23
3.5.3 Screen Selection	25
3.5.4 Main Screen	27
3.5.5 System Mode Control.....	28
3.5.6 Pre-treatment	29
3.5.7 Break Tank	30
3.5.8 Heater	31
3.5.9 RO Feed.....	32
3.5.10 1 st Pass Pump	33
3.5.11 1 st Pass RO.....	34
3.5.12 2 nd Pass Pump	35
3.5.13 2 nd Pass RO.....	36
3.5.14 Overview	37
3.5.15 Status.....	38
3.5.16 Loop Filter (Optional)	40
3.5.17 Active Alarm.....	41
3.5.18 Alarm History	42
3.5.19 Heat Sanitization History	43
3.5.20 Chemical Cleaning History	44
3.5.21 Sequence History.....	45
3.5.22 Setpoint 1	46
3.5.23 Setpoint 2	49
3.5.24 Setpoint 3	52
3.5.25 Setpoint 4	55
3.5.26 Setpoint 5	57
3.5.27 Production and Heat Sanitization Scheduler	61
3.5.28 Program Revision	62

3.5.29	Time and Date.....	63
3.5.30	Users Password.....	64
3.5.31	Reset User Password.....	65
3.5.32	Misc. Devices.....	66
3.5.33	Data Logging.....	67
3.5.34	Remote Nurses Station.....	68
3.6	SETPPOINT REFERENCE TABLE.....	69
4.0	SYSTEM OPERATION	73
4.1	SYSTEM VERIFICATION BEFORE STARTING UP.....	73
4.2	PANEL START-UP AND SHUT DOWN	73
4.2.1	Panel Start-up.....	73
4.2.2	Panel Shut Down.....	73
4.3	MANUAL START-UP / SHUT DOWN PROCEDURE.....	74
4.4	HEAT SANITIZATION	74
4.4.1	Auto Heat Sanitization.....	76
4.4.2	Manual Heat Sanitization	76
4.4.3	Loop and Points of use Heat Sanitization	76
4.5	EMERGENCY BY-PASS OPERATION	77
4.5.1	1st Pass By-Pass Operation.....	77
4.5.2	2nd Pass By-Pass Operation	78
4.5.3	Manual By-pass Operation	78
4.5.4	Conductivity Transmitter Manual By-pass Setpoint Change	79
4.6	RO PERFORMANCE OVERVIEW.....	80
4.6.1	System Performance.....	80
4.6.2	Effects of Temperature on System Performance	82
4.6.3	High Recovery Operation	82
5.0	MAINTENANCE	84
5.1	MAINTENANCE REQUIREMENTS	84
5.1.1	RO System Shutdown Procedure.....	85
5.1.2	RO Pre-filter Cartridge Replacement Procedure	85
5.1.3	Log Sheet	86
5.1.4	Service Schedule.....	88
5.1.5	Alarm / Alert Log Sheet.....	89
5.2	CLEANING OVERVIEW	90
5.3	SANITIZATION OVERVIEW	91
5.4	CHEMICAL SELECTION AND USE OF CLEANING CHEMICALS	93
5.5	GENERAL PRECAUTIONS IN CLEANING CHEMICAL SELECTION AND USAGE	93
5.6	CLEANING/SANITIZATION FREQUENCY	95
5.7	RO CLEANING/SANITIZATION CHEMICALS	95
5.8	RO CLEANING/SANITIZATION CHEMICALS CONCENTRATIONS	95
5.9	NEUTRALIZING CLEANING/SANITIZATION SOLUTIONS	96
5.10	MINNCARE® DISINFECTANT USE.....	96
5.10.1	MINNCARE® Test Strip Overview	99
5.10.2	Test Strip Instruction	99
5.11	SEMI-AUTOMATIC CHEMICAL CLEANING AND SANITIZATION	100
5.11.1	Chemical Cleaning and Sanitization Preparation.....	101
5.12	1ST PASS CHEMICAL CLEANING / SANITIZATION PROCEDURE.....	101
5.13	2ND PASS CHEMICAL CLEANING / SANITIZATION PROCEDURE	103
5.14	LOOP CHEMICAL CLEANING / SANITIZATION PROCEDURE	105

5.15	SERVICE / REPLACEMENT PROCEDURES FOR SYSTEM COMPONENTS.....	107
5.15.1	RO Pre-filter cartridge replacement Procedure.....	107
5.15.2	Loop filter Cartridge Replacement (Optional Component)	108
5.15.3	Compressor	109
5.15.4	Conductivity Transmitter Calibration Procedure.....	112
6.0	SYSTEM TROUBLESHOOTING	117
7.0	SYSTEM STATUS MESSAGES.....	134
8.0	SPARE PARTS LIST	137
9.0	GENERAL LIMITED WARRANTY	140
9.1	CONTACT INFORMATION	140

SAFETY

This Manual describes procedures necessary to operate and maintain the BioPure HX2® system. It may make references to the Component Manuals of other manufacturers.

Please read and follow the instructions in this manual carefully to operate the RO system. Pay particular attention to all danger and caution statements to avoid serious injury to the operator and damage to the system.

This equipment is approved by regulatory agencies. Changes to the design and substitution of components are performed only under design controls and shall be documented. Changes without design controls are in contravention to regulation(s).

Operation of this system outside of its intended use and/or without adequate pre-treatment/post-treatment that meets ISO 13959-11 standard guidelines is not allowed and is in contravention to regulation(s).

DEFINITIONS



WARNING *Indicates a potentially hazardous situation which, if not avoided, could result in death or serious personal injury.*

CAUTION *Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate personal injury or possible damage to the system.*

NOTE *Information that requires special emphasis.*

LABELS

Read all labels and tags attached to the equipment. Personal injury or damage to the equipment could occur if not observed.

	
<p style="text-align: center;">WARNING</p> <p style="text-align: center;">DISCONNECT POWER BEFORE ENTERING</p>	<p style="text-align: center;">AVERTISSEMENT</p> <p style="text-align: center;">COUPER L'ALIMENTATION ELECTRIQUE AVANT D'OUVRIRE CETTE PORTE</p>
<p style="text-align: center;">WARNING</p> <p style="text-align: center;">SYSTEM MAY START AUTOMATICALLY. ISOLATE ALL ENERGY SOURCES PRIOR TO MAINTENANCE</p>	<p style="text-align: center;">AVERTISSEMENT</p> <p style="text-align: center;">SYSTÈME PEUT DEMARRER AUTOMATIQUEMENT. ISOLER TOUTES SOURCES ELECTRIQUES AVANT TOUT ENTRETIEN</p>
<p style="text-align: center;">WARNING</p> <p style="text-align: center;">SYSTEM USES HEAT FOR SANITIZATION. SURFACES MAY BE HOT ENOUGH TO CAUSE PERSONAL INJURY</p>	<p style="text-align: center;">AVERTISSEMENT</p> <p style="text-align: center;">SYSTÈME UTILISE L'EAU CHAUDE POUR DESINFECTION. LES SURFACES PEUVENT ETRE ASSEZ CHAUDE POUR PROVOQUER DES BLESSURES</p>
<p style="text-align: center;">WARNING</p> <p style="text-align: center;">UNINTERRUPTABLE POWER SUPPLY LOCATED INSIDE CONTROL PANEL. DISCONNECT DOES NOT DE- ENERGIZE ALL CIRCUITS</p>	<p style="text-align: center;">AVERTISSEMENT</p> <p style="text-align: center;">BLOC D'ALIMENTATION SANS INTERRUPTION A L'INTERIEUR DU PANNEAU DE CONTRÔLE. DISJONCTEUR N'OUVRE PAS TOUS LES CIRCUITS ELECTRIQUES</p>
<p style="text-align: center;">CAUTION</p> <p style="text-align: center;">VOLTAGE SUPPLIED FROM OTHER SOURCES. DISCONNECT POWER BEFORE ENTERING</p>	<p style="text-align: center;">ATTENTION</p> <p style="text-align: center;">LES SOURCES D'ALIMENTATION ELECTRIQUES FOURNIES PAR D'AUTRES DOIVENT ETRE COUPE AVANT D'OUVRIRE CETTE PORTE</p>

SPECIFICATIONS

BioPure HX2® SYSTEM	
Product Flow	HX2P-05 3 US gpm (11.3 lpm) to 7 US gpm (26.5 lpm) ¹ HX2P-11 7 US gpm (26.5 lpm) to 15 US gpm (56.7 lpm) ¹
System Operation	Loop Direct Feed
Quantity of Pumps	2
Total Motor Horsepower	12.5 HP (60 Hz Models) / 15.5 HP (50 Hz Models)
Heater Power	20 kW
Pre-Filter	Yes
Loop Filter	Optional
Operational Weight	HX2P-05 System: Dry: 2,504 lbs (1,138 kg), Wet: 2,671 lbs (1,214 kg)
	HX2P-11 System: Dry: 3,234 lbs (1,470 kg), Wet: 3,684 lbs (1,675 kg)
Total Power	41 kW
Electrical Requirement	Standard voltages: 208VAC/3Ø/60Hz, 460VAC/3Ø/60Hz, 575VAC/3Ø/60Hz, 380VAC/3Ø/50Hz, 415VAC/3Ø/50Hz
Frame and Panel Colour	RAL 7035 Pebble Grey

TANK	
Tank Volume	100 US gal (378 litres)
Material	304SS Stainless Steel & 304SS Cladding
Finish	Interior: ≈25 Ra , Exterior: Ground & Polish to #4 Finish
Operating Weight	Dry: 381 lbs (173 kg), Wet: 1,214 lbs (551 kg)
Tank Rating	Atmospheric

DIMENSIONS – INCH (MM)			
RO Skid Including Tank	166 (4217) Width	37 (940) Depth (with Drain Header)	82 (2083) Height

¹ Product Flow based on new RO Composite Polyamide membranes and feed water temperature of 25°C (77°F).
 A 3% loss in product flow from the RO system can be expected per degree below 25°C (77°F).

GENERAL MATERIAL OF CONSTRUCTION	
Component	Standard Material
Piping	316 / 316L SS
Tank	304 SS / EPDM
Pumps	316 SS / Cast Iron / Tungsten carbide / Carbon / EPDM / Viton or Silicon carbide (Q) / Silicon carbide (Q) complete with EPDM (E) seals
Solenoid Valves	316 SS / EPDM
Heater Element	316 SS
Valves	316 SS / PTFE
Sensors	316 SS / PTFE / EPDM
Fittings	316 SS
Gaskets	EPDM or Tuf-steel
Pipe Thread Seal	PTFE
Panel	Powder Coated Carbon Steel
Frame	Powder Coated Carbon Steel

PRODUCT WATER COMPONENTS MATERIAL OF CONSTRUCTION	
Component	Standard Material
Piping	316 / 316L SS
Valves	316 SS / PTFE
Conductivity/Temperature Sensors	316 SS / PTFE / EPDM
Fittings	316 SS
Gaskets	EPDM or Tuf-steel
Pressure Sensors	316 SS
Flow Sensors	316 SS

GENERAL PRECAUTIONS

The following are general precautions that must be taken with the BioPure HX2® System.

Please read these precautions before operating the system.

WARNING: *This system is a registered medical device. Operation of this system outside of its intended use and without adequate pre-treatment/post-treatment that meets ISO 13959-11 standards will void the system's medical device license.*

WARNING: *Changes to the design, substitution of components or other changes will void the system's medical device license.*

WARNING: *The input water quality to the RO must meet the National Primary Drinking Water, Standards of the United States Environmental Protection Agency (EPA) or Canadian Federal and Provincial guidelines governing drinking water in Canada.*

WARNING: *During the heat sanitization process the water temperature can reach 98°C (208°F). Extreme caution must be taken when handling components on the system. Only appropriate heat resistant materials can be used for fluid pathways that are sanitized with hot water. Exposed surfaces that cannot be insulated must be properly labeled.*

WARNING: *Only chemicals listed in the manual may be used with the BioPure HX2® system. Other chemicals may not be compatible with the material of construction. When using chemical cleaners or sanitizers, ensure the system is properly rinsed prior to use for dialysis.*

WARNING: *Ensure that the system is connected to a power source in compliance with local and national electrical codes. Failure to comply may create a shock or fire hazard.*

WARNING: *After installation and subsequent use, when any component of the RO system is changed or replaced, the user should conduct appropriate tests to ensure that the revised system meets the initial design criteria and water quality requirements.*

WARNING: *A suitable carbon filter should always be present as pre-treatment for the RO unit in order to prevent chlorine/chloramine damage to the Composite Polyamide RO membranes and harm to the patient. Mar Cor strongly recommends the use of two carbon tanks, used in a series configuration with a minimum empty bed contact time of 10 minutes.*

WARNING: *Some settings are user adjustable via the HMI (Human-Machine-Interface). Only personnel that understand the operation of the BioPure HX2® should change or update any password protected settings.*

- WARNING:** *Do not attempt to tamper with or perform any unauthorized modifications to the device including programming changes. Any unauthorized work on the device could cause damage to the equipment or harm to personnel and patients.*
- WARNING:** *While a water treatment system may produce water of sufficient quality to meet the requirements of ISO 13959-11, distribution of the water may degrade its quality to the point where it no longer meets the requirements of the above mentioned standards if the distribution system is not designed, installed and/or maintained appropriately. Construction and materials of the distribution loop must be compatible with the RO system and avoid dead legs in the pipe routing, and ID of pipe must be passivated. 0.035" wall 316ss tube or 0.065" wall 316LSS tubing with sanitary connection is recommended to withstand hot water temperature up to 98°C (208°F).*
- WARNING:** *Loop piping and any other piping that can be exposed to hot water from the system must be insulated with a minimum thickness of ½". Insulation must be chloride free to prevent degradation of the stainless steel surfaces.*
- WARNING:** *Connect this device to a proper ground connection in accordance with the National Electrical Code if applicable. The system skid, panel and piping shall be properly grounded and bonded for safety.*
- WARNING:** *Surge suppression shall be used on the incoming main power feed to the panel if frequent voltage spikes (e.g. lightning, utility variations) are common in the geological area.*
- WARNING:** *Do not, under any circumstances, remove any Caution, Warning or any other descriptive labels from the devices until the conditions warranting the label are eliminated.*
- WARNING:** *To avoid fire or explosion, do not operate this device in an explosive environment or near flammable materials or anesthetics.*
- WARNING:** *To prevent electrical shock, disconnect the electrical power to the system before servicing.*
- WARNING:** *Ensure that all piping connections are tight to avoid leakage.*
- WARNING:** *Always relieve the pressure in the process lines on the system before disassembly to prevent personal injuries against water spray.*
- WARNING:** *Ensure that there is adequate ventilation around the system to avoid the build-up of chemical fumes.*
- WARNING:** *Follow carefully the manufacturer's safety instructions in their individual manuals and labels on chemical containers.*

- CAUTION:** *To avoid physical and/or equipment harm, be cautious to use appropriate tools to move any components on the RO system (i.e. Use a mobile hoist to remove pump motor).*
- NOTE:** *The term 'RO system', "water system or "water purification system" used throughout this manual equally identify BioPure HX2® System.*
- NOTE:** *Use proper wiring and connection methods to satisfy hospital electrical codes if applicable.*
- NOTE:** *Use proper water and drain connections to prevent contamination of the RO and danger to the patient.*
- NOTE:** *The percent rejection is calculated based on feed water after mixing with some recycled reject water.*
- NOTE:** *Use a clean soft cloth with a mild soap or detergent to clean the HMI display. Dry the display with a chamois or moist cellulose sponge to avoid water spots. Remove fresh paint splashes and grease before drying by rubbing lightly with isopropyl alcohol (70% concentration). Afterward, wash using a mild soap or detergent. Rinse with clean water.*

1.0 BIOPURE HX2® SYSTEM OVERVIEW

The BioPure HX2® is a stand-alone RO system in which the final product water does not come into contact with the patient.

The main core of the RO system consists of two trains (1st & 2nd pass) of membrane elements. By pumping pretreated potable water through the 1st pass membrane elements, product water is produced and is pumped directly through the 2nd pass membrane elements to further improve water quality as the final product water.

Final product water is then fed through an optional loop filter or delivered directly to the points of use where it is consumed by the dialysis machines. Any unused water within the distribution loop is returned back to the feed of the 1st pass membrane elements to be reprocessed through the RO system. The final product water of the RO system shall meet ISO 13959-11 standard for hæmodialysis applications.

NOTE: *Under emergency circumstances, either the 1st or 2nd pass membrane elements can be bypassed to allow the RO system to continue producing purified water.*

A Programmable Logic Controller (PLC) shall monitor and control the operation of the water purification system. The system includes monitoring for pressure, flow, water conductivity, percent rejection, percent recovery, differential pressure, pump speed, Break Tank level and valve position (% open).

The system status and performance can be viewed via the Human-Machine-Interface (HMI). The HMI has a colour touch screen display with virtual push buttons to provide a means of controlling the system and managing alarms and settings. The system has user-adjustable, password protected setpoints for controlling the RO system.

A remote indicator panel is available to monitor critical alarms (Nurses Remote). As an option, additional capabilities can be provided to enable remote monitoring via Ethernet within the dialysis facility or from outside.

The RO system is interlocked with the pre-treatment and the PLC can inhibit a regeneration / backwash during a heat sanitization sequence.

An optional final submicron filter can be installed on the RO product exiting the system and is monitored by the PLC.

System chemical clean sequence can be initiated manually by a Mar Cor Service Technician or operator with equivalent skill at any time it is convenient. The sequence is semi-automatic requiring the user to manually open, close valves, and add chemicals. The technician will be prompted to perform these tasks one step at a time via the HMI. The operator can move back and forth between each step but cannot modify fixed settings. This sequence is password protected.

Heat sanitization can be initiated manually by the user at any time it is convenient, or automatically as predetermined on the system scheduler. The scheduler will allow the user to select what part of the system to sanitize and the day and time to initiate a sanitization (typically in the off hours or weekends).




FIGURE 1-1 BIOPURE HX2® SYSTEM

The BioPure HX2® system is equipped with the following main components:

1. Stainless Steel Break Tank (Used during Hot Water and Chemical Sanitization only)
2. Electric Heater
3. Programmable Logic Controller
4. Stainless Steel VFD Booster Pumps
5. Conductivity and Temperature Sensors
6. Flow Sensors
7. NEMA 4/12 Control Enclosures with HMI colour display
8. All Stainless Steel Piping and Housings

NOTE: *The BioPure HX2® system does not include external pre-treatment and post-treatment devices, piping and fittings to the point of use.*

2.0 INSTALLATION REQUIREMENTS

Refer to the BioPure HX2® “System Installation Manual” for instructions of the following:

- UNPACKING THE SYSTEM
- SPACE / ROOM REQUIREMENTS
- LOADING AND UNLOADING SYSTEM SKID
- PIPING REQUIREMENTS
- PIPING INSULATION
- ELECTRICAL HOOK-UP
- PRE-TREATED FEED WATER REQUIREMENTS

3.0 SYSTEM CONTROL OVERVIEW

The system is controlled by a Programmable Logic Controller (PLC) that sends and receives signals to and from various devices and sensors on the system. The status of selected devices and sensors, as well as the overall system performance can be viewed at various screens of the Human-Machine-Interface (HMI) display located on front of the control panel. Screens containing setpoint values are adjustable with password protection.

NOTE: *There is a red emergency stop push button (E-Stop) on the front of the control panel that allows the users to terminate the system operation upon an emergency situation. Once E-Stop is activated, the system will transition to 'Off' mode and all output devices will stop.*

NOTE: *The remote monitoring HMI at the nursing station is strictly designed for monitoring purposes only. It has absolutely no access to make setpoints adjustment to alter the system performance.*

3.1 SECURITY

Various screens can be viewed and accessed from the HMI display. Each screen serves a specific function to allow users to monitor system status, performance and/or adjust the system setpoints. Some of the screens and setpoints have security restrictions and require a password.

NOTE: *Once the user has logged in and there is no HMI activity detected a 10 minute timer will count down and automatically log out the user. The user can also log out by pressing the icon if desired.*

There are three levels of security and they can be accessed by pressing the login icon on the HMI and must enter the login user name followed by the password. Access levels are as follows:

- Login Name: Operator and Default Password: 1 (Level 1)
 - Mode selection ("OFF", "STANDBY", "AUTOMATIC")
 - Acknowledge alarms
 - Reset alarms
 - No access to setpoints
- Login Name: Supervisor and Default Password: 2 (Level 2)
 - All of the above
 - Modify scheduler
 - Adjust setpoints with User Access Level = 2 or less.
 - Access MAINTENANCE Mode

- Login Name: Administration and Default Password: 3 (Level 3)
 - All of the above
 - Adjust setpoints with User Access Level = 3 or less.
 - Security, change passwords
 - Time and date

3.2 ALARM MANAGEMENT

When an alarm has occurred the user may press the alarm silence button on the HMI panel to silence the alarm.

The user shall determine the corrective action by looking at the alarm banner on the HMI and correct the alarm condition(s). Once corrected, the alarm banner can be cleared and the reset button pressed on the HMI to restart the system.

When the alarm silence is pressed and no action is taken the horn will reactivate after 120 seconds. This cycle will repeat every time the alarm silence button is pressed and no action is taken.

The user will be able to view alarms, both current and historical. The user will be able to do the following;

- Silence the alarm horn by pressing the “SILENCE ALARMS” button.
- Reset all alarms/alerts allowing the system to be restarted by pressing the “ALARM RESET” button.
- Acknowledge and date-stamp the current active alarm by pressing the “ACK” button.
- Acknowledge and date-stamp all active alarms by pressing the “ACK ALL” button

3.2.1 Main Control Panel Status and Alarm Indicator Function

The system status is displayed on the local HMI screen.

The description and function is as follows:

Automatic Green Indicator (AU): The green indicator shall turn on when the system is in AUTOMATIC.

Maintenance Amber Indicator (M): The amber indicator shall turn on when the system is in MAINTENANCE. The indicator shall turn on for the following.

- Pre-treatment Lockout (only if single filter pre-treatment devices are used – non-redundant systems)
- RO Flush
- Loop Flush
- Heat Sanitization
- Chemical Clean

- When any device (control module) is not in auto mode on the HMI.

Alarm Red Indicator (AL): The indicator shall turn on when an alarm has occurred. An alarm banner will appear on the HMI, and will be logged in the alarm history.

Loss of communication between any of the HMI devices (Main or Remote Nurse's Station or Remote Monitoring system) and the PLC will initial a local alert on the device with lost communication. Alert only.

NOTE: *An optional 24VDC remote stack light can be installed if system status and alarm indicator is required.*

3.3 REMOTE NURSE'S STATION MONITOR AND ALARM INDICATOR

The system status is displayed on a remotely mounted HMI display (Nurses Station). It shall be limited on view only information except for silencing the audible alarm. The description and function are as follows:

The remote HMI mimics the main screen on the central system. It will display the system status, operational information, alerts and alarms in real time. If a critical alarm occurs the remote audible alarm will activate. The user can silence the audible alarm from the remote HMI or from the main control panel. The user can only reset the alarm from the central system HMI. The remote nurse's station monitor shall indicate loss of communication with the PLC as an alert. Alert only.

During a heat sanitization or chemical cleaning / sanitization the audible alarm will sound intermittently (every 15 seconds for approx. ½ a second).

The remote audible alarm, if silenced, will re-alarm in less than 120 seconds.

The maximum length for the Ethernet cable segment used to connect the nurse's station to the control panel must not be longer than 300'. If the cable distance is greater than 300' the use of active hardware such as a repeater, switch or router shall be necessary to boost the signal.

3.4 USB DRIVE AND PRINTER (optional)

The system is equipped with two (2) USB ports on the left side of the panel. Description and use of each port is as follows:

The first port is used to connect a USB drive (optional) to collect real time raw data from the HMI in a csv format that can be viewed using Microsoft Excel. When viewed the columns have tags and or descriptions that indicate what the values correspond to. The tags can be referenced on the process flow drawing. The data can be filtered using Excel to create system performance reports. This device is set at the factory level security access.

RowNum	LocalDate	LocalTime	{[HX2]AI3635.Val}	{[HX2]AI3665.Val}	{[HX2]AI3695.Val}
1	09/06/2012	0:00:03	53.81	6.8	46.44
2	09/06/2012	0:00:08	54.75	6.8	46.2
3	09/06/2012	0:00:13	54.75	6.8	45.96
4	09/06/2012	0:00:18	54.75	6.8	45.78
5	09/06/2012	0:00:23	53.81	6.8	45.55
6	09/06/2012	0:00:28	53.81	6.8	45.26
7	09/06/2012	0:00:33	52.88	6.8	45.05
8	09/06/2012	0:00:38	52.88	6.86	44.86
9	09/06/2012	0:00:43	52.88	6.86	44.67
10	09/06/2012	0:00:48	52.88	6.86	44.47
11	09/06/2012	0:00:53	52.88	6.86	44.29
12	09/06/2012	0:00:58	52.88	6.86	44.14
13	09/06/2012	0:01:04	52.88	6.86	44
14	09/06/2012	0:01:09	51.94	6.86	43.9
15	09/06/2012	0:01:14	51.94	6.86	43.8

The second port is used to connect a USB printer (Optional). The system will print out messages to the printer. It will print based on the time interval and message count. This device is set at the factory level security access.

Example:

*09/14/2012 10:06:52 AM SYSTEM SEQUENCE STEP 1150 HEAT SANI 1ST PASS
MATRIX MODE 75 RO 1ST PASS HEAT HOLD*

*09/14/2012 10:07:00 AM SYSTEM SEQUENCE STEP 1150 HEAT SANI 1ST PASS
MATRIX MODE 75 RO 1ST PASS HEAT HOLD*

3.5 THE HMI DISPLAY NAVIGATION

The user can access different screens by pressing the appropriate button on the HMI display. The top of the screen contains basic information, such as the system screen title, the system clock time and date, and the alarm indicator. Each blue button will navigate to a specific screen when pressed. The name shown on the button is the associated screen name.



NOTE: Screens shown under section 3.5 and sub-sections are displaying all function keys available for logged in as *Factory* user. Same screens will show less function keys for different users (i.e. *Administrator, Supervisor, etc.*).

3.5.1 Valve Indicator Positions



Indicates the 2-way valve is closed




Indicates the 2-way valve is open



Indicates the 3-way valve open from the bottom port to the right port

3.5.2 Screen Layout



Multi-Stage Indicator displays one of the forty-seven (47) operation modes:

Operating Sequence Modes

- Mode 1: Off
- Mode 2: Standby
- Mode 3: Sleep
- Mode 4: 1st RO Flush
- Mode 5: 2nd Pass RO Flush
- Mode 6: RO Run To Drain
- Mode 7: Loop Flush
- Mode 8: RO Run To Loop
- Mode 10: 1st Pass Bypass Run To Drain
- Mode 11: 1st Pass Bypass Run To Loop
- Mode 12: 2nd Pass Bypass Run To Drain
- Mode 13: 2nd Pass Bypass Run To Loop
- Mode 14: Safe
- Mode 72: Heat Sanitization Tank Fill
- Mode 73: Heat Sanitization Line Rinse
- Mode 74: RO 1st Pass Heat Ramp
- Mode 75: RO 1st Pass Heat Hold
- Mode 76: RO 1st Pass Temperature Cool Down
- Mode 77: RO 1st Pass Tank Drain
- Mode 78: RO 1st Pass TOC Flush
- Mode 79: RO 2nd Pass Heat Ramp
- Mode 80: RO 2nd Pass Heat Hold
- Mode 81: RO 2nd Pass Temperature Cool Down
- Mode 82: RO 2nd Pass Tank Drain
- Mode 83: RO 2nd Pass TOX Flush
- Mode 84: Loop Heat Ramp
- Mode 85: Loop Heat Hold

Screen Navigation Icons (BLUE buttons):

Screen navigation icons enable HMI user to move to different monitoring screens or control/setpoint setting screens.

LOGIN/LOGOUT = Allow user to log in and gain access to password protected screens.

MODE CONTROL = Press to view and change system operation mode.

SCREEN SELECT = Press to view accessible screens and setpoints.

ACTIVE ALARM = Press to display only alarms currently active.

ALARM RESET = Press to temporarily reset alarms.

ALARM SILENCE = Press to silence audible alarms.

ALARM HISTORY = Press to display all old and active alarms.

ICONS



PRINT SCREEN = Press to print the current screen Icon to the optional printer.



CONTACT INFORMATION = Mar Cor Phone and email information.

Mode 86: Loop Tank Drain

Mode 87: RO Loop Temperature Cool Down

Mode 88: Loop Temperature Cool Down

Mode 89: Loop Heat Hold With Points Of Use

Mode 30: Chemical Clean Tank Fill

Mode 31: Chemical Clean Line Rinse

Mode 32: Chemical Clean Chemical Mixing

Mode 33: Chemical Clean 1st Pass Recycle

Mode 34: Chemical Clean 1st Pass Neutralize

Mode 35: Chemical Clean Tank Drain

Mode 36: Chemical Clean 1st Pass Flush

Mode 37: Chemical Clean 2nd Pass Recycle

Mode 38: Chemical Clean 2nd Pass Neutralize

Mode 39: Chemical Clean 2nd Pass Flush

Mode 40: Chemical Clean Loop Recycle


Mode 41: Chemical Clean Loop Neutralize


Mode 42: Chemical Clean Loop Flush

Mode 43: Chemical Clean 2nd Pass Line Rinse

Mode 45: Chemical Clean Pause For Operation

3.5.3 Screen Selection

Screen Title: Screen Select	Screen Description: Main Screen for Screen Navigation
The purpose of this screen is to allow the USER to navigate to different HMI screens.	
<p>Factory Access</p>  <p>User Screen (No Password)</p>	<p>Multi-Stage Indicator displays one of the forty-six (46) operation modes. See section 3.5.2 to review the list of modes and description on the HMI.</p> <p>MONITORING SCREENS</p> <p>MAIN - Displays system performance on 1st Pass RO water quality, 2nd Pass RO water quality and temperature, Loop Feed flow rate, and Loop Return flow rate.</p> <p>STATUS - Control interface for current operation status and maintenance functions.</p> <p>SCHEDULE - Displays the pre-set time scheduler for system operation in Normal Production and Heat Sanitization.</p> <p>OVERVIEW - Displays real-time system performance.</p> <p>QUALITY TREND - Displays Conductivity measurements versus Time on the system.</p> <p>MATRIX CHECK - Displays the position indicator on the valve and determines its physical position.</p> <p>MISC. DEVICES - Displays operation status of miscellaneous components: Remote Audible Alarm, Audible Alarm, System Alarm Relay, Water Leak, DI Alarm, DI Active, Auxiliary Input, Main Power, Main E-Stop, Alarm Silence, and Compressed Air.</p> <p>PRE-TREATMENT - Displays operation status of automatic pre-treatment systems: Softeners, Media Filters, and Carbon Filters.</p> <p>BREAK TANK - Displays component operation and system performance associated with the break tank. Break tank is strictly used for heat sanitization.</p> <p>HEATER - Displays component operation and system performance associated with the heater. Heater is strictly used for heat sanitization.</p> <p>RO FEED - Display water conditions feeding the RO system.</p> <p>TIME & DATE - Adjustment for HMI time and date.</p> <p>1ST PASS PUMP - Display water conditions feeding the RO system via 1st pass RO booster pump.</p> <p>1ST PASS RO - Display water conditions before and after feeding through the 1st pass RO membrane.</p>

Screen Title: Screen Select


NOTE: Display buttons that are grayed out have no function on the current screen or are password protected.

Screen Description: Main Screen for Screen Navigation

2nd PASS PUMP - Display water conditions feeding the RO system via 2nd pass RO booster pump.

2nd PASS RO - Display water conditions before and after feeding through the 2nd pass RO membrane.

FLOW TABLE - Displays access to the system flow tables to set up the system product output.

LOOP FILTER - Display water conditions associated with the loop filter (Optional).

SETPOINT 1 - Adjustment for various system operation settings.

SETPOINT 2 - Adjustment for various system operation settings.

SETPOINT 3 - Adjustment for various system operation settings.

SETPOINT 4 - Adjustment for various system operation settings.

SETPOINT 5 - Adjustment for various system operation settings.

FACTORY - Displays access to factory setup information.

ACTIVE ALARMS - Displays the active alerts / alarms in the system.

ALARM HISTORY - Displays all the system alerts / alarm in a history log.

HEAT SANITIZATION HISTORY - Displays all the heat sanitization history.

CHEMICAL CLEAN HISTORY - Displays all the chemical clean history.

SEQUENCE HISTORY - Displays the sequence history including modes and steps.

DATA LOGGING - User can disable or enable the data log file when removing the USB storage device located on the left side of the control panel.

PROGRAM REVISION - Displays current PLC and HMI program file name and revision date.

DIAG - Displays diagnostic information on the system operation for troubleshooting.


RESET PASSWORD - Update security passwords for new users.

CHANGE PASSWORD - Manage security passwords for various users.


SYSTEM RESET - Displays access to reset the system

SHUTDOWN HMI - Displays access to shut down and resets the HMI application.


3.5.4 Main Screen

Screen Titles: Main	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: Full Access Supervisor: Full Access</p> <p>Administrator: Full Access</p>  <p>System Performance Indicator</p>	<p>The purpose of this screen is to allow for system performance monitoring. In addition, the screen displays the following information:</p> <p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> 1ST PASS QUALITY - Shows current water quality after 1st pass RO filtration. 2ND PASS QUALITY - Shows current water quality after 2nd pass RO filtration. 2ND PASS TEMP - Shows current water temperature after 2nd pass RO filtration. LOOP FEED FLOW - Shows current water flow rate feeding the distribution loop. LOOP RETURN FLOW - Shows current water flow rate returning from the distribution loop to the break tank.


3.5.5 System Mode Control

Screen Title: Mode Control	Screen Description: View System Operation Mode.
Security: No Login: Partial Access Supervisor: Full Access Administrator: Full Access	The purpose of this screen is to allow the user to change system operating mode.
 <p>COMBO HEAT SANITIZATION</p> <ul style="list-style-type: none"> 2ND PASS & LOOP - Press to manually activate heat sanitizing 2nd Pass RO followed by heat sanitizing the Loop. 1ST & 2ND PASS - Press to manually activate heat sanitizing 1st Pass RO followed by heat sanitizing 2nd Pass RO. 1ST & 2ND PASS & LOOP - Press to manually activate heat sanitizing 1st Pass RO followed by heat sanitizing 2nd Pass RO. Then finally heat sanitizing the Loop. <p>SINGLE HEAT SANITIZATION</p> <ul style="list-style-type: none"> 1ST PASS - Press to manually activate heat sanitizing 1st Pass RO only. 2ND PASS - Press to manually activate heat sanitizing 2nd Pass RO only. LOOP - Press to manually activate heat sanitizing Loop only. LOOP POU - Press to manually activate heat sanitizing Loop with Point Of Use only. 	<p>NORMAL MODES</p> <ul style="list-style-type: none"> ? - Press to access help screen. OFF MODE - Press to shut down system. AUTO MODE - Press to place system in auto mode. SLEEP - Indicate system in non-production. PRODUCTION - Indicate system in production. STANDBY MODE - Press to place system in standby. Auto Flush as scheduled. <p>MAINTENANCE MODES</p> <ul style="list-style-type: none"> MAINT MODE - Press to place system in Maintenance mode. Disable all system activities in Auto mode. RO FLUSH - Press to activate a RO flush. LOOP FLUSH - Press to activate a Loop flush. RO RUN TO DRAIN - Press to divert product water to drain. RO RUN TO LOOP - Press to divert product water to Loop. 1ST PASS BYPASS - Press to run 2nd pass RO only. 2nd PASS BYPASS - Press to run 1st pass RO only. <p>CLEANING MODES</p> <ul style="list-style-type: none"> CLEAN MODE - Press to place system in Clean mode. 1ST PASS - Press to set all automated valve inline for chemical clean 1st Pass RO only. 2ND PASS - Press to set all automated valve inline for chemical clean 2nd Pass RO only. LOOP - Press to set all automated valve inline for chemical clean Loop only.


3.5.6 Pre-treatment

Screen Titles: Pre-treatment	Screen Description: System Monitoring
Security: No Login: View Only Supervisor: View Only Administrator: View Only	The purpose of this screen is to allow for system performance monitoring. In addition, the screen displays the following information:
 <p>The screenshot shows the 'PRETREATMENT' screen with the following sections:</p> <ul style="list-style-type: none"> SOFTENERS: <ul style="list-style-type: none"> UA1800: On, XS1800: Off UA1805: On, XS1805: Off MEDIA FILTERS: <ul style="list-style-type: none"> UA1835: On, XS1835: Off UA1840: On, XS1840: Off CARBON: <ul style="list-style-type: none"> UA1810: On, XS1810: Off UA1815: On, XS1815: Off UA1820: On, XS1820: Off GENERAL PRETREATMENT INTERLOCK: <ul style="list-style-type: none"> UA1850: On <p>A red arrow points from the 'System Performance Indicator' label to the 'UA1850' status indicator.</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> UA1800 & UA1805 = Softener #1 and or #2, Service= On, Regeneration = Off. UA1810, UA1815 & UA1820 = Carbon Filter #1, #2, and or #3, Service=On, Backwash=Off. UA1835 & UA1840 = Media Filter #1 and or #2, Service=On, Backwash=Off. UA1850 = General Pre-treatment Interlock is On when pre-treatment system stops supplying feed water to the RO system. The system operation is then automatically switched to standby until feed water supply is resumed. XS_____ = Indicate the pre-treatment device is lockout during a heat sanitization. This depends on type of timer used for the pre-treatment.


3.5.7 Break Tank

Screen Title: Break Tank	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Administrator: View Only</p> <p>Supervisor: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> YV3400 - Drain cooling water valve. YV3410 - Tank drain valve. YV3405 - Feed water divert valve. LI3615 - Tank level indicator. ZC3402 - Feed water control valve % open indicator. F3005 - Pre-filter Indicator. PI3600 - Filter inlet pressure indicator. PI3605 - Filter outlet pressure indicator. CPI9850 - Pre-filter calculated differential pressure indicator. CFI9866 - Calculated feed flow indicator. CFI9868 - Calculated feed flow totalizer.


3.5.8 Heater

Screen Title: Heater	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> TI3630 - Loop return temperature indicator. TI3640 - Feed water temperature. CF19874 - Calculated loop return flow velocity indicator. FI3620 - Loop return flow indicator. EC3010 - Heater output Indicator. YC3010 - Heater power contactor. TAHH3610 - Heater high temperature switch indicator. PID9400 - Heater control Indicator. YV3415 - Tank return fill valve. YV3420 - Tank return drain valve.


3.5.9 RO Feed

Screen Title: RO Feed	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> FI3655 - 1st pass reject flow indicator. FI3660 - 1st pass reject to drain flow indicator. FI3690 - 2nd pass recycle flow indicator. CFI9872 - Calculated 1st pass recycle flow indicator. AI3635 - Feed conductivity indicator. TI3640 - Feed temperature indicator. ZC3440 - 1st pass reject control valve % open indicator. PID9430 - 1st pass reject control indicator. ZC3445 - 1st pass reject drain control valve % open indicator. PID9440 - 1st pass reject drain control indicator. ZC3460 - 2nd pass recycle control valve % opened indicator. PID9470 - 2nd pass recycle control indicator. YV3430 - Feed water valve indicator. YV3450 - 1st pass reject drain divert valve indicator. YV3465 - 2nd pass recycle bypass valve.


3.5.10 1st Pass Pump

Screen Title: 1st Pass Pump	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>Pump Run Time Meter</p> <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> AI3635 - Feed conductivity indicator. TI3640 - Feed temperature indicator. SIC3015 - 1st pass pump speed indicator. PI3645 - 1st pass pump discharge pressure indicator. FI3710 - 2nd Pass product flow indicator. PID9410 - 1st pass pump pressure control indicator. PID9420 - 1st pass pump flow control indicator. YV3435 - Loop bypass divert valve indicator. YV3475 - 1st pass product bypass valve indicator.


3.5.11 1st Pass RO

Screen Title: 1 st Pass RO	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> AI3665 - 1st pass product conductivity indicator. TI3670 - 1st pass product temperature indicator. FI3655 - 1st pass reject flow indicator. CFI9870 - Calculate 1st pass product flow indicator. PI3650 - 1st Pass reject pressure indicator. CPI9854 - 1st pass calculated differential pressure Indicator. CFI9860 - 1st pass calculated % array recovery indicator. CFI9858 - 1st pass calculated % system recovery indicator. ZC3440 - 1st pass reject control valve % open indicator. PID9430 - 1st pass flow control indicator. YV3475 - 1st pass product bypass valve indicator. YV3455 - 2nd pass product bypass valve indicator. F3020 - 1st pass membrane indicator.


3.5.12 2nd Pass Pump

Screen Title: 2 nd Pass Pump	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>The screenshot shows the '2ND PASS PUMP' monitoring interface. It includes a schematic of the pump system with a bypass valve (YV3455) and a pump. Key indicators include:</p> <ul style="list-style-type: none"> Pump Run Time Meter: 1302.0 HOURS System Performance Indicator: 1.64 PSI PID9450 PRESS: PV 1.64 PRESS, SP 1.64 PRESS, CV 0.01 % FLOW PID9450 FLOW: PV 0.00 FLOW, SP 0.00 FLOW, CV 0.00 % PID9465 VEL: PV 0.00 VEL, SP 0.00 VEL, CV 0.01 % FLOW PID9465 FLOW: PV 0.00 FLOW, SP 0.00 FLOW, CV 0.00 % PID9460 FLOW: PV 0.00 FLOW, SP 0.00 FLOW, CV 0.33 % PRES PID9460 PRESS: PV 1.64 PRESS, SP 1.64 PRESS, CV 0.00 % <p>Buttons at the bottom include: LOGOUT, MODE, SCREEN, ACTIVE ALARM, ALARM, ALA, SILEN, ALARM, HISTORY.</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> SIC3025 - 2nd pass pump speed indicator. PID9440 - 2nd pass pump maximum flow control indicator. PID9450 - 2nd pass pump pressure control indicator. PID9460 - 2nd pass pump flow control indicator. PI3675 - 2nd pass pump discharge pressure indicator. FI3710 - 2nd pass product flow indicator. YV3455 - 2nd pass product bypass valve indicator.


3.5.13 2nd Pass RO

Screen Title: 2nd Pass RO	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> A3695 - 2nd pass product conductivity indicator. T3700 - 2nd pass product temperature indicator. F3710 - 2nd pass product flow indicator. F3690 - 2nd pass recycle flow indicator. ZC3460 - 2nd pass recycle flow control valve % open indicator. P3685 - 2nd pass recycle pressure indicator. P3705 - 2nd pass product pressure indicator. CPI9856 - 2nd pass calculated differential pressure Indicator. CFI9862 - 2nd pass calculated % system recovery indicator. CQI9864 - System percent rejection indicator. YV3455 - 2nd pass product bypass valve indicator. YV3435 - Loop bypass divert valve indicator. YV3470 - 2nd pass product divert valve indicator. (Loop Feed)

3.5.14 Overview


Screen Title: Overview	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>System Performance Indicator</p> <p><i>NOTE: Due to the amount of information displayed most of the units of measure are not shown. Refer to the individual screens for more details.</i></p>	<p>System Performance Indicator will show the following information:</p> <p>This is a general overview of the system piping and instrumentation. The previous screens provide details of each instrument and indicator as per section listed below:</p> <ul style="list-style-type: none"> • Section: 3.5.7 • Section: 3.5.8 • Section: 3.5.9 • Section: 3.5.10 • Section: 3.5.11 • Section: 3.5.12 • Section: 3.5.13

3.5.15 Status


Screen Title: Status	Screen Description: System Monitoring and Maintenance
<p>Security:</p> <p>No Login: View Only Supervisor: Partial Access Administrator: Partial Access</p>	<p>The purpose of this screen is to allow the user to monitor the status of the system at any time and perform chemical sanitization or cleanings. The user would first place the system in appropriate cleaning mode from the mode control screen in section 3.5.5.</p> <p>In the middle of the screen the following information will be based on what the system is doing.</p> <ul style="list-style-type: none"> • Current step the system is in. • What the system is waiting for to move on. • Time remaining in the step or a delay. • Instructions and or prompts to the user.
 <p>NOTE: If the indicator is greyed out it cannot be accessed.</p>	<p>System Performance Indicator will show the following information:</p> <p>1st Pass Recycle - Press this indicator to perform a 1st pass recycle of the cleaning or sanitization solution.</p> <p>2nd Pass Recycle - Press this indicator to perform a 2nd pass recycle of the cleaning or sanitization solution.</p> <p>Loop Recycle - Press this indicator to perform a loop recycle of the cleaning or sanitization solution.</p> <p>1st Pass Flush - Press this indicator to flush out the 1st pass of any cleaning or sanitization solution.</p> <p>2nd Pass Flush - Press this indicator to flush out the 2nd pass of any cleaning or sanitization solution.</p> <p>Loop Flush - Press this indicator to flush out the loop of any cleaning or sanitization solution.</p> <p>Tank Drain - Press this indicator to drain the tank.</p> <p>Tank Fill - Press this indicator to fill the tank with clean permeate water.</p> <p>Line Rinse - Press this indicator to perform a line.</p> <p>Chemical Mixing - Press this indicator to mix the chemicals added to the tank and verify concentration.</p> <p>RO Run To Drain - Press this indicator to place the system in RO run to drain to check system performance after cleaning and or sanitization.</p> <p>Neutralize - Press this indicator to perform a neutralization of the cleaning solution before dumping it to drain.</p> <p>2nd Pass Line Rinse - Press this indicator to perform a 2nd pass rinse the line.</p>

Screen Title: Status	Screen Description: System Monitoring and Maintenance
	<p>Next - Press this indicator to move to the next cleaning or sanitization step.</p> <p>Clean Complete - Press this indicator exit the cleaning sequence. This is logged in the history.</p> <p>Sanitization Complete - Press this indicator exit the sanitization sequence. This is logged in the history.</p> <p>Pump Start - Press this indicator to start the pump. The indicator will turn from red to green when the pump is on.</p> <p>Pump Stop - Press this indicator to stop the pump. The indicator will turn from red to green when the pump is stopped.</p> <p>Stop - Press this indicator to stop a cleaning or sanitization. This can be done to add more chemical. If in a recycle mode the clean duration timer will pause.</p> <p>Clean Dur 20 Min - Press this indicator to change the clean or sanitization recycle duration. This can be done when the pump is running. Once the timer has counted down the pump will stop automatically.</p> <p>Reset Timer - Press this indicator to reset the clean duration timer.</p> <p>Bump Timer - Press this indicator to bump the clean duration timer to zero.</p> <p>Abort Heat Sani - Press this indicator to abort a heat sanitization.</p>


3.5.16 Loop Filter (Optional)

Screen Title: Loop Filter	Screen Description: System Monitoring
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to allow for system performance monitoring of the RO skid. In addition, the screen displays the following information:</p>
 <p>The screenshot displays the 'LOOP FILTER' screen. At the top, it shows the date and time (JUN 28, 2016 15:09:47) and the status (OFF MODE, STEP 0 OFF, MATRIX 1 OFF, SCHEDULE IN PRODUCTION MODE). The main area shows a process flow diagram with a filter labeled F3035. A pressure indicator PI3715 is shown with a value of -24.99 PSI. A red callout box points to the indicator with the text 'System Performance Indicator'. The bottom of the screen has a blue bar with buttons: LOGOUT, MODE CONTROL, SCREEN SELECT, ACTIVE ALARM, ALARM RESET, ALARM SILENCE, and ALARM HISTORY.</p>	<p>System Performance Indicator will show the following information:</p> <ul style="list-style-type: none"> PI3715 - Loop filter F3035 outlet pressure / Loop feed pressure. CPI9852 - Loop filter calculated differential pressure indicator. F3035 - Loop filter indicator.


3.5.17 Active Alarm

Screen Title: Active Alarm	Screen Description: System Monitoring (Popup Screen)
Security: No Login: View Only Supervisor: View Only Administrator: View Only	The purpose of this screen is to display all active alarms /messages. In addition, the screen displays the following information:
 <div data-bbox="365 1182 669 1252" style="border: 1px solid red; padding: 5px; display: inline-block;"> Alarms/Messages Display Box </div>	Operator interfaces To view all active alarms/messages in the display box. The operator can acknowledge (Ack) one alarm/message at a time or all (Ack All). They are time and date stamped when they occur and when they are acknowledged. When they are acknowledged they move to the history log. The operator can scroll up and down one alarm/message at time using the single arrows or to the very top or bottom using the double arrows. !=Active The alarm/messages with an (!) beside them have not been acknowledged.


3.5.18 Alarm History

Screen Title: Alarm History	Screen Description: System Monitoring (Popup Screen)
Security: No Login: View Only Supervisor: View Only Administrator: View Only	The purpose of this screen is to display all the alarms/messages. In addition, the screen displays the following information:
 <p>NOTE: The recorded events cannot be erased manually.</p>	Operator interfaces To view all history alarm/message display box. The operator can view the history of alarms/messages. They are time and date stamped when they occurred. The operator can acknowledge (Ack) one alarm/message at a time or all (Ack All). They are time and date stamped when they occur and when they are acknowledged. The operator can scroll up and down one alarm/message at time using the single arrows or to the very top or bottom using the double arrows. !=Active The alarm/messages with an (!) beside them have not been acknowledged. *=Ack'ed The alarm/messages with a (*) beside them have been acknowledged. Based on the HMI internal memory capacity the history log can record approximately 4,000 events including the heat sanitization and cleaning history. Once full, the next event recorded will in turn erase the oldest event recorded.


3.5.19 Heat Sanitization History

Screen Title: Heat Sani History	Screen Description: System Monitoring (Popup Screen)
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to display the heat sanitization messages. In addition, the screen displays the following information:</p>
 <p>Heat Sani History Display Box</p>	<p>Operator interfaces</p> <p>To view all executed heat sanitization in the display box.</p> <p>The operator can view the history of messages. They are time and date stamped when they are recorded.</p> <p>The operator can scroll up and down one message at time using the single arrows or to the very top or bottom using the double arrows.</p> <p>Once a heat sanitization begins and the heater ramps up in temperature and is recorded in the history every 10 minutes. Once the sanitization is completed it also recorded.</p>


3.5.20 Chemical Cleaning History

Screen Title: Chem Clean History	Screen Description: System Monitoring (Popup Screen)
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to display all cleaning messages. In addition, the screen displays the following information:</p>
 <p>Chemical Clean Display Box</p>	<p>Operator interfaces</p> <p>To view all executed cleanings in the display box.</p> <p>The operator can view the history of messages. They are time and date stamped when they are recorded.</p> <p>The operator can scroll up and down one message at a time using the single arrows or to the very top or bottom using the double arrows.</p> <p>The start and completion of the chemical clean are recorded in the history.</p>

3.5.21 Sequence History

Screen Title: Sequence History	Screen Description: System Monitoring (Popup Screen)
<p>Security:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of this screen is to display all the sequence messages. In addition, the screen displays the following information:</p>
 <p>Sequence History Display Box</p>	<p>Operator interfaces</p> <p>To view all sequences in the display box.</p> <p>The operator can view the history of messages. They are time and date stamped when they are recorded.</p> <p>The operator can scroll up and down one message at time using the single arrows or to the very top or bottom using the double arrows. This is a good troubleshooting tool.</p>


3.5.22 Setpoint 1

Screen Title: Setpoint 1	Screen Description: System Monitoring & Setpoint Adjustment
Security: No Login: View Only Supervisor: View Only With Limited Access To Setpoint Administrator: View Only With Limited Access To Setpoint	The purpose of this screen is to allow for viewing or changing system setpoints. In addition, the screen displays the following information:
 <p>Setpoint Boxes</p>	Operator Interfaces To adjust any setpoints on this screen, press the desired “Setpoint Box” and then change the entry as required. A numerical entry keypad will appear and it will only accept numbers with working range indicated at the top of keypad. Softener Configuration Choice between: Duplex (Parallel) and Single Default: Single If Duplex is selected the system will remain in operation if one of the softeners regenerates via an open signal (dry contact). If Single is selected the system will go into off mode once the softener regeneration is complete it will go back into the mode it was in prior to the softener regenerating. NOTE: This is set up at the factory or 1st time start-up. UA1850 General Pre-treatment Interlock RO On/Off Choice between: RO on and RO off Default: RO Off If RO On is selected the system will remain in operation if the general pre-treatment device sends an open signal (dry contact) indicating it is out of service, general interlock or an alarm. If RO Off is selected the system will go into off mode if the general pre-treatment device sends an open signal (dry contact) and will go back into the mode it was in prior, once the signal closes. NOTE: This is set up at the factory or the 1st time start-up. General pre-treatment devices could be a chemical feed pump or multi-layer filter...

Screen Title: Setpoint 1	Screen Description: System Monitoring & Setpoint Adjustment
<p><u>Setpoint 1 Screen Navigation Shortcut (BLUE buttons):</u> Quick access to view different setpoint setting screens.</p> <div data-bbox="205 337 310 391" style="background-color: #0056b3; color: white; padding: 2px; text-align: center; margin-bottom: 5px;">SETPOINT 2</div> <div data-bbox="310 370 661 391" style="margin-left: 10px;">- Press to display Setpoint 2 screen.</div> <div data-bbox="205 407 310 461" style="background-color: #0056b3; color: white; padding: 2px; text-align: center; margin-bottom: 5px;">SETPOINT 3</div> <div data-bbox="310 435 661 456" style="margin-left: 10px;">- Press to display Setpoint 3 screen.</div> <div data-bbox="205 477 310 531" style="background-color: #0056b3; color: white; padding: 2px; text-align: center; margin-bottom: 5px;">SETPOINT 4</div> <div data-bbox="310 500 661 521" style="margin-left: 10px;">- Press to display Setpoint 4 screen.</div> <div data-bbox="205 535 310 589" style="background-color: #0056b3; color: white; padding: 2px; text-align: center; margin-bottom: 5px;">SETPOINT 5</div> <div data-bbox="310 565 661 586" style="margin-left: 10px;">- Press to display Setpoint 5 screen.</div>	<p><u>LAL3615 POU Sanitization Tank Level Low Alert</u> Setpoint Range: 6 to 80 gallons (23 to 302 litres) Default: 25 gallons (95 litres) This will provide an alert to the operator when drawing water off during point of use heat sanitization so they have time to close the points of use before the tank reaches the low level alarm. NOTE: This is set up based on the operator water volume requirement for a POU heat sanitization.</p> <p><u>CAL 9874 Loop Return Flow Velocity Low Alert</u> Setpoint Range: 1.0 to 10.0 fts (0.305 to 3.05 mps) Default: 1.5 fts (0.457 mps) This will provide an alert, after a 10 min. time delay, to the operator that the loop return flow velocity is below the setpoint indicating the water consumption from the loop is too high and should be reduced. NOTE: This is set up at the factory or the 1st time start-up and the calculation is based on the loop internal diameter setpoint entered in section 3.5.2.</p> <p><u>KAHH 3640 Heat Sanitization Ramp Too Long Alarm</u> Setpoint Range: 120 to 480 minutes Default: 240 minutes This timer will count down during the heat ramp up of a hot water sanitization. If the final temperature hold setpoint is not reached in the time set the system will abort the sanitization. This setpoint is set based on loop length and expected temperature losses. NOTE: This is set up at the 1st time start-up.</p> <p><u>AAH 3655 Feed Conductivity High Alert</u> Setpoint Range: 1000 to 3000 µS/cm Default: 2000 µS/cm This provides an alert to the operator that the feed conductivity is above the setpoint. This would indicate that there is a pre-treatment issue or an upset in the potable feed water source. NOTE: This is set up at 1st time start-up.</p>

Screen Title: Setpoint 1	Screen Description: System Monitoring & Setpoint Adjustment
	<p><u>CAH 9854 1st Pass Membrane Pressure Differential High Alert</u> Setpoint Range: 5 to 60 psig Default: 10 psig This provides an alert to the operator the 1st pass membrane pressure differential is above the setpoint. This could indicate that a cleaning of the first pass membrane is required. NOTE: This is set up at the 1st time start-up and after every cleaning. It is set to 15% above the normal operating 1st pass membrane pressure differential.</p> <p><u>AAH 3665 1st Product Conductivity High Alert</u> Setpoint Range: 1.0 to 100.0 µS/cm Default: 10.0 µS/cm This provides an alert to the operator that the 1st pass product conductivity is above the setpoint. This could indicate that a cleaning of the first pass membrane is required or the membranes need to be replaced. NOTE: This is set up at 1st time start-up, after a cleaning or as the membranes age.</p> <p><u>CAH 9856 2nd Pass Membrane Pressure Differential High Alert</u> Setpoint Range: 5 to 60 psig Default: 10 psig This provides an alert to the operator the 2nd pass membrane pressure differential is above the setpoint. This could indicate that a cleaning of the first pass membrane is required. NOTE: This is set up at 1st time start-up and after every cleaning. It is set to 15% above the normal operating 2nd pass membrane pressure differential</p>


3.5.23 Setpoint 2

Screen Title: Setpoint 2	Screen Description: System Monitoring & Setpoint Adjustment
Security: No Login: View Only Supervisor: View Only With Limited Access To Setpoint Administrator: View Only With Limited Access To Setpoint	The purpose of this screen is to allow for viewing or changing system setpoints. In addition, the screen displays the following information:
 <p>Setpoint Boxes</p>	Operator Interfaces To adjust any setpoints on this screen, press the desired “Setpoint Box” and then change the entry as required. A numerical entry keypad will appear and it will only accept numbers with working range indicated at the top of keypad. <u>AAH 3695 2nd Pass Product Conductivity High Alert</u> Setpoint Range: 1.0 to 100.0 μ S/cm Default: 10.0 μ S/cm This provides an alert to the operator that the 2nd pass product conductivity is above the setpoint. This could indicate that a cleaning of the second pass membrane is required or the membranes need to be replaced. NOTE: This is set up at 1st time start-up, after a cleaning or as the membranes age. <u>AAHH 3695 2nd Pass Product Conductivity High High Alarm</u> Setpoint Range: 1.0 to 100.0 μ S/cm Default: 20.0 μ S/cm This provides an alarm to the operator that the 2nd pass product conductivity is above the setpoint. The 2nd pass product will divert to drain and the system will continue to operate. If the quality does not improve after 10 minutes the system will shut down. This could indicate that a cleaning of the second pass membrane is required or the membranes need to be replaced. NOTE: This is set up at 1st time start-up, after a cleaning or as the membranes age. <u>CAL 9864 System Percent Reject Alert</u> Setpoint Range: 70.0 to 99.0 % Default: 95.0%

Screen Title: Setpoint 2	Screen Description: System Monitoring & Setpoint Adjustment
<p><u>Setpoint 2 Screen Navigation Shortcut (BLUE buttons):</u> Quick access to view different setpoint setting screens.</p> <div data-bbox="205 337 310 386">SETPOINT 1</div> - Press to display Setpoint 1 screen. <div data-bbox="205 402 310 451">SETPOINT 3</div> - Press to display Setpoint 3 screen. <div data-bbox="205 467 310 516">SETPOINT 4</div> - Press to display Setpoint 4 screen. <div data-bbox="205 532 310 581">SETPOINT 5</div> - Press to display Setpoint 5 screen.	<p>This provides an alert to the operator that the system percent salt rejection is below the setpoint. This calculation based on your feed conductivity versus you 2nd pass conductivity.</p> <p>NOTE: This is set up at 1st time start-up, after a cleaning or as the membranes age.</p> <p><u>SP1 3615 Loop And Point Of Use Heat Sanitization Tank Full</u> Setpoint Range: 30 to 80 gallons (113 to 302 litres) Default: 80 gallons (302 litres) This setpoint is used to set the chemical clean tank fill level. This is based on the amount of purified water the operator requires during the chemical cleaning or sanitization. This setpoint can also be access from the status screen in section 3.5.15.</p> <p>NOTE: This is set by the operator. The system will perform a line rinse after the tank is filled. Approximately 10 gallons (37.9 litres) is sent to drain during the sequence. This needs to be taken into account when determining the desired water level.</p> <p><u>SP3 3615 Chemical Clean Tank Full</u> Setpoint Range: 40 to 80 gallons (151 to 302 litres) Default: 80 gallons (302 litres) This setpoint is used to set the chemical clean tank fill level. This is based on the amount of purified water the operator requires during the chemical cleaning or sanitization. This setpoint can also be access from the status screen in section 3.5.15.</p> <p>NOTE: This is set by the operator. The system will perform a line rinse after the tank is filled. Approximately 10 gallons (37.9 litres) is sent to drain during the sequence. This needs to be taken into account when determining the desired water level.</p> <p><u>SP4 3615 Break Tank Gallons / Inch Of WC</u> Setpoint Range: 0.1 to 20.0 gallons Default: 2.5823 This setpoint is a multiplier used for converting inches of water to gallons.</p> <p>NOTE: This is set by the Mar Cor Tech 2.</p>

Screen Title: Setpoint 2	Screen Description: System Monitoring & Setpoint Adjustment
	<p><u>SP1 3402 Heat Sanitization Cool Down Feed Water Valve % Open</u></p> <p>Setpoint Range: 0 to 100 % valve opened Default: 20 %</p> <p>This setpoint is used to set the valve opening to control cool. It is an upper limit setting.</p> <p><i>NOTE: This is set by the supervisor.</i></p>


3.5.24 Setpoint 3

Screen Title: Setpoint 3	Screen Description: System Monitoring & Setpoint Adjustment
Security: No Login: View Only Supervisor: View Only With Limited Access To Setpoint Administrator: View Only With Limited Access To Setpoint	The purpose of this screen is to allow for viewing or changing system setpoints. In addition, the screen displays the following information:
 <p>Setpoint Boxes</p>	Operator Interfaces To adjust any setpoints on this screen, press the desired “Setpoint Box” and then change the entry as required. A numerical entry keypad will appear and it will only accept numbers with working range indicated at the top of keypad. <u>SP1 9880 Heat Sanitization Temperature Hold</u> Setpoint Range: 65°C to 85°C (149°F to 185°F). Default: 80°C (176°F) This setpoint is used to set the hot water sanitization temperature hold setpoint and is measured at the 2nd pass product during a 1st and 2nd pass heat sanitization and at the loop return during a POU and loop heat sanitization. Once the setpoint is reached the hold timer will start to count down. NOTE: This is set at 1st time start-up. <u>SP2 9885 Heat Sanitization Hold Time</u> Setpoint Range: 20 to 120 minutes Default: 60 minutes This setpoint is used to set the hot water sanitization temperature hold time. The timer will count down once the temperature hold setpoint is reach. The timer will pause if the temperature should fall below the setpoint but no more than the low temperature setpoint. Once the temperature ramps back up to the setpoint the timer will restart. NOTE: This is set at 1st time start-up. <u>SP2 9880 POU Loop Heat Sanitization Hold Time</u> Setpoint Range: 10 to 480 minutes Default: 480 minutes

Screen Title: Setpoint 3	Screen Description: System Monitoring & Setpoint Adjustment
<p><u>Setpoint Screen Navigation Shortcut (BLUE buttons):</u></p> <p>Quick access to view different setpoint setting screens.</p> <div data-bbox="205 337 310 394">SETPOINT 1</div> - Press to display Setpoint 1 screen. <div data-bbox="205 402 310 459">SETPOINT 2</div> - Press to display Setpoint 2 screen. <div data-bbox="205 467 310 524">SETPOINT 4</div> - Press to display Setpoint 4 screen. <div data-bbox="205 532 310 589">SETPOINT 5</div> - Press to display Setpoint 5 screen.	<p>This setpoint is used to set the hot water temperature hold time for sanitizing the POU Loop. The timer will count down once the temperature hold setpoint is reach. The timer will pause if the temperature should fall below the setpoint but no more than the low temperature setpoint. Once the temperature ramps back up to the setpoint the timer will restart.</p> <p><i>NOTE: This is set at 1st time start-up.</i></p> <p><u>SP5 9880 Membrane Clean/Sanitization Pressure</u></p> <p>Setpoint Range: 30 to 60 psig Default: 50 psig</p> <p>This setpoint is used to set the pump operating pressure during a chemical cleaning / sanitization.</p> <p><i>NOTE: This is set by the operator.</i></p> <p><u>SP7 9880 Line Rinse Pump Speed</u></p> <p>Setpoint Range: 0 to 50% Default: 10 %</p> <p>This setpoint is used to set the pump speed to rinse the pump inlet and outlet piping or pre-treated water during a heat sanitization for cleaning.</p> <p><i>NOTE: This is set up at the factory or 1st time start-up.</i></p> <p><u>SP1 3645 1st Pass Pump Pressure</u></p> <p>Setpoint Range: 30 to 250 psig Default: 150 psig</p> <p>This setpoint is used to set the pump speed at start-up during normal operation. Once at pressure and in production the system will adjust the pressure and flows automatically to maintain the required 2nd pass product flow.</p> <p><i>NOTE: This is set up at the factory or 1st time start-up.</i></p> <p><u>SP2 3645 Loop Heat Sanitization Pressure</u></p> <p>Setpoint Range: 10 to 99 psig Default: 60 psig</p> <p>This setpoint is used to set the pump speed during a loop heat sanitization.</p> <p><i>NOTE: This is set up at 1st time start-up.</i></p>


Screen Title: Setpoint 3	Screen Description: System Monitoring & Setpoint Adjustment
	<p>WARNING: <i>Verify the loop material maximum operating pressure versus temperature. Set the heat sanitization pump speed 10 % below that pressure. Possible hot water leaks and or flooding can occur if operating above the loop material specifications</i></p> <p><u>SP1 3675 2nd Pass Pump Pressure</u> Setpoint Range: 30 to 250 psig Default: 150 psig This setpoint is used to set the pump speed at start-up during normal operation. Once at pressure and in production the system will adjust the pressure and flows automatically to maintain the required 2nd pass product flow. NOTE: <i>This is set up at the factory or 1st time start-up.</i></p> <p><u>SP2 3655 1st Pass Reject Flow</u> Setpoint Range: 10 to 40 gpm (38 to 151 lpm) Default: Select entry from Flow Table This setpoint is used to set the 1st pass reject flow. This setpoint is automatically entered from a predefined flow chart. NOTE: <i>This is set up at the factory.</i></p>

3.5.25 Setpoint 4

Screen Title: Setpoint 4	Screen Description: System Monitoring & Setpoint Adjustment
Security: No Login: View Only Supervisor: View Only With Limited Access To Setpoint Administrator: View Only With Limited Access To Setpoint	The purpose of this screen is to allow for monitoring of system setpoints. In addition, the screen displays the following information:
 <p style="text-align: center;">Setpoint Boxes</p>	Operator Interfaces To adjust any setpoints on this screen, press the desired “Setpoint Box” and then change the entry as required. A numerical entry keypad will appear and it will only accept numbers with working range indicated at the top of keypad. <u>SP1 3660 1st Pass Reject Drain Flow</u> Setpoint Range: 1 to 20 gpm (4 to 76 lpm) Default: Select entry from Flow Table This setpoint is used to set the 1st pass reject drain flow at initial start-up. This setpoint is automatically entered from a predefined flow chart. NOTE: This is set up at the factory. <u>SP1 3690 2nd Pass Recycle Flow</u> Setpoint Range: 1 to 20 gpm (4 to 76 lpm) Default: Select entry from Flow Table This setpoint is used to set the 2nd pass recycle flow at initial start-up. This setpoint is automatically entered from a predefined flow chart. NOTE: This is set up at the factory. <u>SP1 3710 2nd Pass Maximum Product Flow</u> This setpoint is used to set the 1st pass maximum product flow. This setpoint is automatically entered from a predefined flow chart. NOTE: This is set up at the factory. <u>Feed Low Pressure Auto Retry</u> Setpoint Range: 0 to 5 Default: 0

Screen Title: Setpoint 4	Screen Description: System Monitoring & Setpoint Adjustment
<p><u>Setpoint Screen Navigation Shortcut (BLUE buttons):</u> Quick access to view different setpoint setting screens.</p> <div data-bbox="203 370 308 423"> <p>SETPOINT 1</p> </div> <p>- Press to display Setpoint 1 screen.</p> <div data-bbox="203 435 308 488"> <p>SETPOINT 2</p> </div> <p>- Press to display Setpoint 2 screen.</p> <div data-bbox="203 500 308 553"> <p>SETPOINT 3</p> </div> <p>- Press to display Setpoint 3 screen.</p> <div data-bbox="203 565 308 618"> <p>SETPOINT 5</p> </div> <p>- Press to display Setpoint 5 screen.</p>	<p>This setpoint is used to set how many times the system will attempt an auto restart after a low feed pressure alarm shut down has occurred. NOTE: This is set at 1st time start-up.</p> <p><u>Feed Low Pressure Retry Delay Time</u> Setpoint Range: 10 to 60 minutes Default: 10 minutes This setpoint is used to set the time before the system attempts to perform an auto restart after a low feed pressure alarm shut down has occurred. NOTE: This is set at 1st time start-up or if a feed water low pressure issue at the facility is common.</p> <p><u>Loop Enabled</u> Choice between Yes and No. Default: Yes If yes is selected the system will perform its normal auto flush sequence through the loop during non -production time. If no is selected the system will perform its normal auto flush sequence through the RO only during non-production times. NOTE: This is set up at the factory or if maintenance is required in the loop.</p>

3.5.26 Setpoint 5


Screen Title: Setpoint 5	Screen Description: System Monitoring & Setpoint Adjustment
Security: No Login: View Only Supervisor: View Only With Limited Access To Setpoint Administrator: View Only With Limited Access To Setpoint	The purpose of this screen is to allow for monitoring of system setpoints. In addition, the screen displays the following information:
 <p>The screenshot shows the 'SET POINT 5' screen. At the top, it displays the date 'JUN 28, 2010', time '14 : 58 : 11', and 'COM' status with three red lights. Below this, it says 'OFF MODE', 'STEP 0 OFF', 'MATRIX 1 OFF', and 'SCHEDULE IN PRODUCTION MODE'. The main area contains a list of parameters with their current values and units: AUTO FLUSH CYCLE TIME (60 MIN), AUTO FLUSH DURATION (90 SEC), RO FLUSH WITH PUMP (YES), CLEANING ELAPSE TIMER (30 MIN), TEMPERATURE DISPLAY UNITS (°C, °F, °C), VOLUME DISPLAY UNITS (LITRES, GALLONS, GALLONS), VELOCITY DISPLAY UNITS (MPS, FPS, FPS), STOP ALARM - ENABLE/DISABLE (DISABLE, ENABLE, ENABLED), STOP ALARM DISABLE DURATION (REMAINING - 02 : 00 : 00) (120 MIN), LOOP PIPING INTERNAL DIAMETER (1.50 INCH), and 15 MINUTE WATER SAMPLE ALERT (YES, NO, NO). At the bottom, there are buttons for 'SETPOINT 1' through 'SETPOINT 4', 'POINT 5', 'LOGOUT Administrator', 'MODE CONTROL', 'SCREEN SELECT', 'ACTIVE ALARM', 'ALARM RESET', 'ALARM SILENCE', and 'ALARM HISTORY'. A red arrow points to the 'SETPOINT 1' through 'SETPOINT 4' buttons, which are labeled 'Setpoint Boxes'.</p>	Operator Interfaces To adjust any setpoints on this screen, press the desired “Setpoint Box” and then change the entry as required. A numerical entry keypad will appear and it will only accept numbers with working range indicated at the top of keypad. Auto Flush Cycle Timer Setpoint Range: 30 to 120 minutes Default: 60 minutes This setpoint is used to set the time between auto flushes of the system including the loop (if loop enabled) during non-production times. NOTE: This is set up at 1st time start-up. Auto Flush Duration Setpoint Range: 60 to 900 seconds Default: 90 seconds This setpoint is used to set how long the auto flush will be active. NOTE: This is set up at 1st time start-up. RO Flush With Pump Choice between Yes and No. Default: Yes If yes is selected the pumps will turn on during a flush if no is selected the pump will not turn on. NOTE: This is set up at 1st time start-up. CAUTION: If no is select the system will perform a flush at normal feed pressure. The flush will not be as effective. This option should only be used temporarily if there are feed water pressure issues.

Screen Title: Setpoint 5	Screen Description: System Monitoring & Setpoint Adjustment
<p><u>Setpoint Screen Navigation Shortcut (BLUE buttons):</u> Quick access to view different setpoint setting screens.</p> <div data-bbox="205 337 310 386">SETPOINT 1</div> - Press to display Setpoint 1 screen. <div data-bbox="205 402 310 451">SETPOINT 2</div> - Press to display Setpoint 2 screen. <div data-bbox="205 467 310 516">SETPOINT 3</div> - Press to display Setpoint 3 screen. <div data-bbox="205 532 310 581">SETPOINT 4</div> - Press to display Setpoint 4 screen.	<p><u>Clean Elapse Timer</u> Setpoint Range: 0 to 120 minutes Default: 30 minutes This setpoint is used to set the duration of a chemical cleaning or sanitization recycle sequence. The timer will start when the pump is turns on and stop after it is complete. This setpoint can be accessed from the status screen during when performing a chemical cleaning or sanitization. See section 0. NOTE: <i>This is set up when performing a chemical cleaning and or sanitization.</i></p> <p><u>Temperature Display Units</u> Setpoint Range: °C or °F Default: °F This setpoint is used to change the temperature units of measure. NOTE: <i>1st time start-up.</i> CAUTION: <i>System must be in Off mode when making this change. This will prevent the system from shutting down unexpectedly when converting the settings.</i></p> <p><u>Volume Display Units</u> Setpoint Range: Litres or Gallons (US) Default: Gallons This setpoint is used to change the volume units of measure. NOTE: <i>1st time start-up.</i> CAUTION: <i>System must be in Off mode when making this change. This will prevent the system from shutting down unexpectedly when converting the settings.</i></p> <p><u>Velocity Display Units</u> Setpoint Range: MPS or FPS Default: FPS This setpoint is used to change the loop velocity units of measure. MPS (Meters Per Second), FPS (Feet Per Second)</p>


Screen Title: Setpoint 5	Screen Description: System Monitoring & Setpoint Adjustment
	<p>NOTE: 1st time start-up.</p> <p>CAUTION: System must be in Off mode when making this change. This will prevent the system from shutting down unexpectedly when converting the settings.</p> <p><u>Stop Alarm Enable/Disable</u> Setpoint Range: Enable or Disable Default: Enable This setpoint is used to enable or disable the stop function when an alarm occurs.</p> <p>NOTE: This is set up at the factory and if set to disabled shall only be used for troubleshooting or during emergency situations. The heater safety switch and poor quality divert are still active.</p> <p>CAUTION: If set to disabled the system will not stop on certain alarms. The operator shall monitor the system manually during this time to prevent damage to system.</p> <p><u>Stop Alarm Disable Duration (Remaining - *:*)</u> Setpoint Range: 10 to 120 minutes Default: 10 minutes This setpoint is used to the disable time counter. When the stop alarm is disabled the counter will begin to count down. Once it has elapsed the stop alarm setpoint will automatically go back to enabled.</p> <p>NOTE: This is set up at the factory.</p> <p><u>Loop Piping Internal Diameter</u> Setpoint Range: 0.00 to 2.00 inches Default: 0.625 inches This setpoint is used to enter the loop internal diameter. This used to calculate the loop return flow velocity.</p> <p>NOTE: This is set up at the 1st time start-up. Make sure it is the installed loop internal diameter that is entered.</p> <p><u>15 Minute Water Sample Alert</u> Setpoint Range: Yes or No</p>

Screen Title: Setpoint 5	Screen Description: System Monitoring & Setpoint Adjustment
	<p>Default: No</p> <p>When set to yes the system will provide an alert to the operator that daily samples can be taken. The alert will occur 15 minutes after the 2nd pass product water has diverted to the loop. This will occur when the system transition from standby / sleep to production.</p> <p>NOTE: This is set up at the 1st time start-up.</p>


3.5.27 Production and Heat Sanitization Scheduler

Screen Title: Schedule	Screen Description: Production and Heat Sanitization Schedule Setting
Security: No Login: View Only Supervisor: Full Access Administrator: Full Access	The purpose of this screen is to allow user to pre-set automatic system production and heat sanitization schedules.
 <p>Multi-selection Button</p> <p>Schedule Indicator</p> <p>CAUTION: System must be in Off mode when making changes to the schedule. This will prevent the system from shutting down unexpectedly.</p> <p>NOTE: If the stop production time is changed while the system is running and the time set is before the original stop time the system will continue to operate until the next stop time in the next day.</p>	Scheduling Procedure Place the system in Off mode when making schedule changes. Press the desired Schedule Indicator to: <ul style="list-style-type: none"> Select DAY of the week to ENABLED or DISABLED production schedule and heat sanitization schedule. Set time of the day to start production and heat sanitization (i.e. Enter hour and minute). Set time of the day to stop production and heat sanitization (i.e. Enter hour and minute). Only ONE sanitization sequence can be executed in a day. Available sequences are as follow: <ul style="list-style-type: none"> 1ST PASS = Heat sanitize 1st pass RO only 2ND PASS = Heat sanitize 2nd pass RO only LOOP = Heat sanitize loop only 1ST & 2ND PASS = Heat sanitize 1st pass RO, 2nd pass RO LOOP POU = Heat sanitize loop and point of use. 1ST & 2ND & LOOP = Heat sanitize 1st pass RO, 2nd pass RO, loop. If multiple sequences are set, only one sequence will heat sanitize at a time followed by the next one. System must be in Sleep (Non-Production) or Standby for a heat sanitization to occur.


3.5.28 Program Revision

Screen Title: Program Revision	Screen Description: Display Software Information																		
Security: No Login: View Only Supervisor: View Only Administrator: View Only	This screen displays information about the current revision of HMI and PLC code.																		
 <p> Software Version: 2.00.06 </p> <table border="1"> <thead> <tr> <th>Field</th> <th>Value</th> </tr> </thead> <tbody> <tr> <td>HMI IP:</td> <td>0.0.0.0</td> </tr> <tr> <td>HMI MASK:</td> <td>0.0.0.0</td> </tr> <tr> <td>HMI GW:</td> <td>0.0.0.0</td> </tr> <tr> <td>HMI MAC:</td> <td>00-ff-8e-86-01-5b</td> </tr> <tr> <td>PLC IP:</td> <td>192.168.53.100</td> </tr> <tr> <td>PLC MASK:</td> <td>255.255.255.0</td> </tr> <tr> <td>PLC GW:</td> <td>192.168.53.53</td> </tr> <tr> <td>PLC MAC:</td> <td>f4-f4-f3-f4-f5-f0</td> </tr> </tbody> </table> <p> Navigation Bar: LOGOUT, MODE, SCREEN, ACTIVE ALARM, ALARM, ALARM, ALARM </p>	Field	Value	HMI IP:	0.0.0.0	HMI MASK:	0.0.0.0	HMI GW:	0.0.0.0	HMI MAC:	00-ff-8e-86-01-5b	PLC IP:	192.168.53.100	PLC MASK:	255.255.255.0	PLC GW:	192.168.53.53	PLC MAC:	f4-f4-f3-f4-f5-f0	
Field	Value																		
HMI IP:	0.0.0.0																		
HMI MASK:	0.0.0.0																		
HMI GW:	0.0.0.0																		
HMI MAC:	00-ff-8e-86-01-5b																		
PLC IP:	192.168.53.100																		
PLC MASK:	255.255.255.0																		
PLC GW:	192.168.53.53																		
PLC MAC:	f4-f4-f3-f4-f5-f0																		


3.5.29 Time and Date

Screen Title: Time And Date	Screen Description: Date and Time Adjustment
Security: No Login: View Only Supervisor: Full Access Administrator: Full Access	The purpose of this screen is to allow for adjusting the system time/date/year in the HMI.
 <p>Current Time & Date Adjustment Boxes</p> <p>DST Fall Indicator Boxes</p> <p>DST Spring Indicator Boxes</p>	Date And Time Adjustment <ul style="list-style-type: none"> To change current date and time, enter desired year/date/time at the Current Time & Date Adjustment Boxes. Once entry is complete, press the SET TIME AND DATE icon. To set the DST Spring time, enter the desired date/time at the indicator boxes under the SET DST SPRING heading. To set the DST Fall time, enter the desired date/time at the indicator boxes under the SET DST FALL heading. <p>NOTE: The AUTOMATIC DAYLIGHT SAVING TIME ENABLED icon must be press to YES for automatic time change. Otherwise use the + PLUS 1 HOUR icon or the - MINUS 1 HOUR icon to manually adjust current date and time to DST.</p> <p>NOTE: If a heat sanitization is scheduled within the one hour time switch over the sanitization will not occur.</p>


3.5.30 Users Password

Screen Title: CHANGE PASSWORD	Screen Description: Change Password For An Existing User
Security: No Login: No Access Supervisor: No Access Administrator: Full Access	The purpose of this screen is to allow for modifying user accounts.
	Operator interfaces <p>Step 1: Press the LOGIN icon and enter the desired user name and then enter the associated password. Login Name: Supervisor Default Password: 2 Login Name: Administrator Default Password: 3</p> <p>Step 2: Press the CHANGE PASSWORD icon to bring up the pop-up screen. First type in the User Account. Then enter the old and new password according to the labeled boxes. Then enter new password again in the confirmation box and press the ENTER symbol.</p>


3.5.31 Reset User Password

Screen Title: RESET PASSWORD	Screen Description: Reset Password For A New User
Security: No Login: No Access Supervisor: No Access Administrator: Full Access	The purpose of this screen is to allow for resetting password for new user accounts.
	Operator interfaces Step 1: Enter the desired user name in the User Account box. User Names: Supervisor; Administrator. Step 2: Enter new password and re-enter to confirm new password entry. Step 3: Press the CHANGE PASSWORD icon. The RESULT box will report the results of the reset operation as one of the following: <ul style="list-style-type: none"> • Success. • Internal error, cannot access security component. • Internal error, cannot update MER File with new password. • User %s credentials are incorrect. • %s is replaced with specified user account. • User %s account has been disabled. • %s is replaced with specified user account. • Password cannot exceed %d characters. • %d is read from FactoryTalk Security policy. • "The password cannot be changed. • FactoryTalk Security policy is preventing the operation from succeeding. • Password must be active for over %d days before it can be changed. • %d is read from FactoryTalk Security policy settings. • The password has already been used. • FactoryTalk Security policy is preventing the operation from succeeding. • Password must be at more than %d characters long. • %d is read from FactoryTalk Security policy settings. • Password does not meet complexity requirements. • FactoryTalk Security policy is preventing the operation from succeeding. • Windows-Linked user account password can only be changed in Windows. • Password cannot be changed. Access denied. • FactoryTalk Security is preventing the operation from succeeding.


3.5.32 Misc. Devices

Screen Title: MISC. DEVICES	Screen Description: Misc. Device Current State
<p>Security:</p> <p>No Login: View Only Administrator: View Only</p> <p>Supervisor: View Only</p>  <p>The screenshot shows the 'MISC DEVICES' screen with a status bar at the top indicating 'OFF MODE', 'STEP 0', 'OFF', 'MATRIX 1 OFF', and 'SCHEDULE IN PRODUCTION MODE'. The screen is divided into several sections with toggle buttons for various devices: REMOTE STACK LIGHT GREEN (XL9830, Off), REMOTE AUDIBLE ALARM (YC9945, Off), WATER LEAK (YC9957, On), MAIN POWER (YC9800, On), REMOTE STACK LIGHT YELLOW (XL9835, On), AUDIBLE ALARM (YC9820, Off), DI ALARM (YC9980, On), MAIN E-STOP (YC9805, On), REMOTE STACK LIGHT RED (XL9840, Off), SYSTEM ALARM RELAY (YC9815, Off), DI ACTIVE (YC9985, On), ALARM SILENCE (YC9810, Off), AUXILIARY INPUT (YC9990, On), and COMPRESSED AIR (PI6640, 95.54 PSI). At the bottom, there is a navigation bar with buttons for LOGOUT, MODE CONTROL, SCREEN SELECT, ACTIVE ALARM, ALARM RESET, ALARM SILENCE, and ALARM HISTORY.</p>	<p>The purpose of this screen is to allow the user to view the current state and or reading of the miscellaneous devices.</p> <p>System Performance Indicator will show the following:</p> <p>Remote Audible Alarm YC9945 – [Off] Audible remote nurses alarm is not active; [On] Audible remote nurses alarm is active.</p> <p>Audible Alarm YC9820 - [Off] Audible panel alarm is not active; [On] Audible panel alarm is active.</p> <p>System Alarm Relay YC9815 - [Off] The panel general alarm relay is active; [On] The panel general alarm relay is not active.</p> <p>Water Leak YV9957 - [On] External leak monitor is not in alarm, [Off] External leak monitor is in alarm. System will go to off mode.</p> <p>DI Alarm YC9980 - [On] External DI monitor is not in alarm, [Off] External DI monitor is in alarm.</p> <p>DI Active YV9985 - [On] External DI monitor is not active, [Off] External DI monitor is active. System will go to off mode.</p> <p>Auxiliary Input YC9990 - [On] Auxiliary input is not active, [Off] Auxiliary is active. System will go to off mode.</p> <p>Main Power YV9800 - [On] Main power detection monitor is not in alarm. [Off] Main power detection monitor is in alarm.</p> <p>Main E-Stop YC9805 - [On] Main emergency stop button is not in alarm, [Off] Main emergency stop button is in alarm. System will go to off mode.</p> <p>Alarm Silence YC9810 - [On] Alarm silence is not active; [Off] Alarm silence is active. This is a momentary button.</p> <p>Compressed Air PI6640 - Indicates the current air pressure on the system mounted air compressor.</p>

3.5.33 Data Logging

Screen Title: DATA LOGGING	Screen Description: Data Logging
<p>Security:</p> <p>No Login: No Access Supervisor: No Access</p> <p>Administrator: Full Access</p>	<p>The purpose of this screen is to control the data logging file.</p>
	<p>From this screen the user can disable or enable the data log file when removing the USB storage device located on the left side of the control panel.</p> <p>Press enabled when a new storage device is connected.</p> <p>Press disabled before removing a storage device or when changing the file destination.</p> <p>CAUTION: <i>Failure to performing the above instruction can lead to possible corruption of the data log file and or data logging application stopping.</i></p>

3.5.34 Remote Nurses Station

Screen Description: Nurses Remote Display	
<p>No Security Password Required:</p> <p>No Login: View Only Supervisor: View Only</p> <p>Administrator: View Only</p>	<p>The purpose of these screen is provide information to the operator from the remote nurse's station.</p>
 <p>Mode Indicator</p> <p>Sequence in Step</p> <p>Step</p> <p>Current Mode State</p> <p>Software Version</p> <p>Alarm Banner</p> <p>System Performance Indicator</p>	<p>System Performance Indicator will show the following:</p> <ul style="list-style-type: none"> OFF - When lit up it indicates the system is off. STANDBY - When lit up it indicates the system is in standby. AUTO - When lit up it indicates the system is on auto. 1ST PASS QUALITY - Shows current water quality after 1st pass RO filtration. 2ND PASS QUALITY - Shows current water quality after 2nd pass RO filtration. 2ND PASS TEMP - Shows current water temperature after 2nd pass RO filtration. LOOP FEED FLOW - Shows current water flow rate feeding the distribution loop. LOOP RETURN FLOW - Shows current water flow rate returning from the distribution loop to the break tank. SILENCE ALARM - Push button to silence the audible alarm. ALARM HISTORY - Push button to view history. CLOSE - Push button to close history and return to main screen.

3.6 SETPOINT REFERENCE TABLE

Tag #	Parameter Name	Setting Range					
		Security level*	Lower Limit	Upper Limit	Default	Customer Setting	Unit
No Tag	Softener Configuration	4	Duplex or Single		Single	No Unit	No Tag
UA 1850	General Pre-treatment Interlock RO On/Off	n/a	RO On or RO Off		RO Off	No Unit	UA 1850
LAL 3615	POU Sanitization Tank Level Low Alert	5	6	80	25		gallons
			23	302	189		litres
LALL 3615	Tank Level Low Low Alarm	4	Fixed		5		gallons
			Fixed		19		litres
LAHH 3615	Tank Level High High Alarm	4	Fixed		82		gallons
			Fixed		310		litres
CAL 9874	Loop Return Flow Velocity Low Alert	4	1.0	10.0	1.5		fts
			0.305	3.05	0.915		mps
CAH 9850	Pre-filter Pressure Differential High Alert	4	Fixed		≥10		psig
PAL 3605	Feed Pressure Low Low Alarm	4	Fixed		≤ 15		psig
TAHH 3610	Heater Element Temperature High High Alarm	4	Fixed		212		°F
			Fixed		100		°C
TALL 3630	Heat San Temperature Low Low Alarm	4	Fixed		CTC1 9800		°F
			Fixed		CTC1 9800		°C
AAH 3635	Feed Conductivity High Alert	2	1000	3000	≥ 2000		μS/cm
TAHH1 3640	Feed Temperature High High Alarm	4	Fixed		≥ 104		°F
			Fixed		≥ 40		°C
TAHH2 3640	Heat Sanitization Temperature High High Alarm	n/a	Fixed		≥ 208		°F
			Fixed		≥ 98		°C
KAHH 3640	Heat Sanitization Ramp Too Long Alarm	2	120	480	240		min
PAHH1 3645	1 st Pass Pump Pressure High High Alarm	4	Fixed		≥ 250		psig
PAHH2 3645	1 st Pass Heat Sani Pressure High High Alarm	4	Fixed		≥ 30		psig
AAH 3665	1 st Pass Product Conductivity High Alert	2	1.0	100.0	≥ 10		μS/cm
FALL 3660	1 st Pass Reject Drain Flow Low Low Alarm	4	Fixed		1.00		gpm
					3.78		litres

Tag #	Parameter Name	Setting Range					
		Security level*	Lower Limit	Upper Limit	Default	Customer Setting	Unit
CAH 9854	1st Pass Membrane Pressure Differential High Alert	n/a	5	60	10		psig
AAH 3665	1st Pass Product Conductivity High Alert	2	1.0	100.0	≥ 10.0		μS/cm
PAHH1 3675	2 nd Pass Pump Pressure High High Alarm	4	Fixed		≥ 280		psig
PAHH2 3675	2 nd Pass Heat Sani Pressure High High Alarm	4	Fixed		≥ 30		psig
CAH 9856	2nd Pass Membrane Pressure Differential High Alert	2	5	60	10		psig
AAH 3695	2nd Pass Product Conductivity High Alert	2	1.0	100.0	≥ 10.0		μS/cm
AAHH 3695	2nd Pass Product Conductivity High High Alarm	2	1.0	100.0	≥ 20.0		μS/cm
TALL 3700	Heat San Temperature Low Low Alarm	n/a	Fixed		CTC1 9880		°F
			Fixed		CTC1 9880		°C
PAHH 3705	2nd Pass Product Pressure High High Alarm	4	Fixed		≥ 95		psig
FALL 3710	2nd Pass Product Flow Low Low Alarm	4	Fixed		1.00		gpm
			Fixed		3.78		litres
CAH 9852	Loop Filter Pressure Differential High Alert	4	Fixed		≥ 15		psig
PAL 6640	Compressed Air Pressure Low Alert	4	Fixed		85		psig
PALL1 6440	Compressed Air Pressure Low Low Alarm	4	Fixed		75		psig
PALL2 6440	Compressed Air Pressure Low Low Alarm	n/a	Fixed		65		psig
SP1 3615	Loop And Point Of Use Heat Sanitization Tank Full	2	30	80	80		gallons
			113	302	302		litres
SP2 3615	RO And Loop Heat San Tank Full	n/a	Fixed		30		gallons
			Fixed		113		litres
SP3 3615	Chemical Clean Tank Full	2	30	80	80		gallons
			113	302	302		litres
SP1 3402	Heat San Cool Down Feed Water Valve % Open	2	0	100	20		%

Tag #	Parameter Name	Setting Range					
		Security level*	Lower Limit	Upper Limit	Default	Customer Setting	Unit
CTC1 9880	Heat San Temperature Low Setpoint	n/a	Fixed		SP19880 minus 9		°F
			Fixed		SP19880 minus 5		°C
SP1 9880	Heat San Temperature Hold	5	149	185	176		°F
			65	85	80		°C
SP2 9880	Heat San Hold Time	2	20	120	60		min
SP3 9880	Heat San Cool Down Temperature	n/a	Fixed		113		°F
			Fixed		45		°C
SP4 9880	Membrane Heat San Pressure	n/a	Fixed		25		psig
SP5 9880	Membrane Clean/San Pressure	2	30	60	50		psig
SP6 9880	System High Recovery	2	Fixed		>300=75 300=80 150=85 100=90 50=95		µS/cm = %
SP7 9880	Line Rinse Pump Speed	4	0	50	10		%
SP8 9880	Tank Empty Level	n/a	Fixed		10		gallons
					37.8		litres
SP1 3645	1st Pass Pump Pressure	2	30	250	150		psig
SP2 3645	Loop Heat San Pressure	4	10	99	60		psig
SP1 3675	2nd Pass Pump Pressure	2	30	250	150		psig
SP1 3655	1st Pass Reject Flow	2	Factory Set				gpm
			Factory Set				lpm
SP1 3660	1st Pass Reject Drain Flow	2	Factory Set				gpm
			Factory Set				lpm
CFC1 3660	1 st Pass Reject To Drain Flush Flow	n/a	Factory Set		1.25 multiply SP1 3660		gpm
							lpm
CFC2 3660	1 st Pass Reject To Drain Bypass Flow	n/a	Factory Set		SP1 3710		gpm
							lpm
SP1 3690	2nd Pass Recycle Flow	2	Factory Set				gpm
			Factory Set				lpm
SP1 3710	2nd Pass Maximum Product Flow	5	Factory Set				gpm
			Factory Set				lpm
SP1 3620	Loop Return Minimum Flow	4	1	20	3		gpm
			4	76	11		lpm

Tag #	Parameter Name	Setting Range					
		Security level*	Lower Limit	Upper Limit	Default	Customer Setting	Unit
No Tag	Feed Low Pressure Auto Retry	2	0	5	0		No Unit
No Tag	Feed Low Pressure Retry Delay Time	2	10	60	10		min
No Tag	Auto Flush Cycle Time	2	30	120	60		No Unit
No Tag	Auto Flush Duration	2	60	900	90		sec
No Tag	Loop Enabled?	n/a	Yes or No		Yes		No Unit
No Tag	RO Flush With Pump	4	Yes or No		Yes		No Unit
No Tag	Cleaning Elapse Timer	2	0	120	30		min
No Tag	Temperature Display Units	2	°C or °F		°F		°C or °F
No Tag	Volume Display Units	2	Litres or gallons		gallons		Litres or gallons
No Tag	Velocity Display Units	2	fps or mps		fps		fps or mps
No Tag	Stop Alarm – Enable/Disable	2	Enable or Disable		Enable		No Unit
No Tag	Stop Alarm Disable Duration (Remaining – * : *)	2	10	120	10		min
No Tag	Loop Pipe Internal Diameter	4	0.00	2.00	0.625		Inch
No Tag	15 Minute Water Sample Alert	2	Yes or No		No		No Unit
No Tag	Loop Filter Option	n/a	Enable or Disable		Disable		No Unit

Security Level* Defined:

- 1 Operator
- 2 Supervisor
- 3 Administrator
- 4 Mar Cor Technician 1 Access²
- 5 Mar Cor Technician 2 Access³
- n/a Fixed Set Points

² Security Level 4 is only accessible by a trained Mar Cor Technician

³ Security Level 5 is only accessible by a trained Mar Cor Technician

4.0 SYSTEM OPERATION

4.1 SYSTEM VERIFICATION BEFORE STARTING UP

Prior to starting the BioPure HX2® system, verify the following is complete. Once the system has been flooded with pretreated water it must be placed into operation within a day.

1. Temperature blending valve is installed, flushed and temperature is set to 77°F (25°C).
2. Softener(s) are installed, programmed, flushed and outlet water hardness tested.
3. A minimum of two carbon filters are installed in series, programmed (if automatic), flushed and outlet water tested for chlorine.
4. Pressure regulator after the last pre-treatment step is installed and set to between 16 and 30 psig.
5. Pre-treatment electrical interlocks are installed and verified.
6. Electrical supply is connected and power is adequate to support the RO system.
7. Interconnect and loop piping including drains are complete and verified.
8. All system pipelines are pressure tested for leaks and flushed to remove debris.
9. RO pre-filter cartridges and/or optional loop filter cartridges are installed.
10. All other accessories are installed and ready to receive water.

NOTE: *Other pre-treatment devices (Media filter, Chemical Injection, Iron Filter.....) may be required depending on the potable feed water analysis. Please consult with a Mar Cor sales representative.*

4.2 PANEL START-UP AND SHUT DOWN

The panel is equipped with a 120VAC uninterruptible power supply (UPS) located in the upper control cabinet. It will provide approximately 20 minutes of control power to the system. You must turn off the main power and UPS to remove power from the system.

4.2.1 Panel Start-up

1. Open the top control panel and press and hold the UPS power button for 3 seconds to turn it on. The UPS is located in the top left hand corner.
2. Turn the main power disconnect to the ON position. It is located on the top right hand corner of the lower panel.
3. Close the top panel and secure door.
4. Allow the panel to initialize for 5 minutes.

4.2.2 Panel Shut Down

1. Turn the main power disconnect to the OFF position. It is located on the top right hand corner of the lower panel.
2. Open the top control panel and press and hold the UPS power button for 3 seconds to turn it off. The UPS is located in the top left hand corner.

WARNING: *Uninterruptible power supply located inside control panel. Main disconnect does not de-energize all circuits.*

WARNING: *Use proper safety procedures when working in the electrical panel. Disconnect and/or lock-out the power supply to the system and tag it out. The system must be properly grounded to prevent personal injury or damage to the system.*

4.3 MANUAL START-UP / SHUT DOWN PROCEDURE

After the system has been placed in operation for the 1st time it will not require any user input to operate normally with the condition that the system is controlled based on the scheduled time programmed in the calendar. If desired the nurse / operator (no password) can start the system when required and place the system in standby at the end of the day:

1. From the SCHEDULE screen (security level Supervisor or higher), disable the production schedule for each day.
2. From the MAIN screen press the AUTO icon to place the system into production when required.
3. From the MAIN screen press the STANDBY icon to place the system into standby when required. Placing the system in standby will allow the scheduled heat sanitization to initiate as programmed.

WARNING: *Do not place the system in OFF unless maintenance is required for facility issues. Placing the system in OFF mode for more than 24 hours can result in excessive bacterial growth. Should this occur; a full system heat sanitization and water quality sampling must be performed prior to use.*

NOTE: *When the system is in full production (Mode #8) and a power failure should occur the system will attempt to perform a rapid recovery but only if the main power to the panel is restored in less than 5 minutes. If the power is restored the system will restart in Mode #8 in less than 5 seconds and ramp up to the last production flow rate in less than 15 seconds. If the power is restored after 5 minutes the system will perform the normal start-up flushes to assure the water and microbiological quality is met.*

4.4 HEAT SANITIZATION

Heat sanitization can be initiated manually by the user at any time in Standby or automatically based on the system scheduler when in Standby, or in Automatic Mode, during Sleep.

The system scheduler will allow the user to select what part of the system to sanitize, the day and hour to initiate a sanitization, typically in the off hours or weekends.

The sequence is password protected requiring a security level Supervisor or higher. The scheduler will allow a single heat sanitization cycle per day, and can be one of the following choices:

Heat Sanitization Selection	Type	Approximate Duration
1 st Pass	Manual or Auto	≈ 4 hrs
2 nd Pass	Manual or Auto	≈ 4 hrs
Loop	Manual or Auto	≈ 4 hrs
Loop And Points of Use	Manual	User Aborted or Programmable Timer
1 st Pass → 2 nd Pass	Manual or Auto	≈ 8 hrs
2 nd Pass → Loop	Manual or Auto	≈ 8 hrs
1 st Pass → 2 nd Pass → Loop	Manual or Auto	≈ 12 hrs

NOTE: *Duration time will vary depending on loop length, hold temperature, hold time, ambient temperature, insulation type and site installation. Times indicated above include a hold temperature of 80°C (176°F), hold time of 60 minutes and loop length of 305 meters (1,000 ft.) partially insulated.*

The typical RO membrane sanitization frequency is once a month to a maximum of once a week. A loop sanitization can be performed on a daily basis. If more frequent sanitizations are required it will affect the normal life span or calibration frequency of critical components on the system.

The schedulers within the system are monitored by the PLC. The PLC will close a digital output to delay the regeneration and/or backwash initiation sequence of the pre-treatment (if equipped) during heat sanitization modes. Once the sanitization is completed the PLC will open the digital output.

During a heat sanitization the horn will sound every 15 seconds for 0.5 seconds. Every 10 minutes the following message will appear on the HMI and recorded in the history.

XXXXX Heat Sanitization In Progress.
Temperature XX.X °C

Once a heat sanitization is complete the following message will appear on the HMI and recorded in the history.

XXXXX Heat Sanitization Completed.
XXX Minutes at XX.X °C
Started HH:MM:SS YYYY/MM/DD
Completed HH:MM:SS YYYY/MM/DD

If for some reason the sanitization should abort the following message will appear on the HMI and recorded in the history.

XXXXX Heat Sanitization Aborted.
Started HH:MM:SS YYYY/MM/DD

Aborted HH:MM:SS YYYY/MM/DD

WARNING: *During the heat sanitization process the water temperature can reach 98°C (208°F). Extreme caution must be taken when handling components on the system.*

WARNING: *Verify the loop material maximum operating pressure versus temperature. Set the heat sanitization pump speed at least 10 % below that pressure. Possible hot water leaks and or flooding can occur if operating above the loop material specifications.*

4.4.1 Auto Heat Sanitization

An auto heat sanitization can be scheduled during non-usage periods. It is important to leave enough time for the heat sanitization to complete. If pre-treatment devices used have a regeneration or backwash schedule, ensure they do not conflict with the auto heat sanitization schedule unless they are equipped with a lock-out input connected to the panel. An auto heat sanitization will only occur if the system is in Standby or Auto (Sleep).

Refer to section 3.5.29 for instructions on how to set the time schedule for auto heat sanitizations.

4.4.2 Manual Heat Sanitization

A manual heat sanitization can be scheduled during non-usage periods. It is important to leave enough time for the heat sanitization to complete. If pre-treatment devices used have a regeneration or backwash schedule, ensure they do not conflict with the auto heat sanitization schedule unless they are equipped with a lock-out input connected to the panel.

1. From the Mode Control screen, log in with username and password.
2. Press the Standby mode icon if the system is operating.
3. Press the desired heat sani icon in the single heat sanitization headline and allow the system to complete its cycle. The sequence of heat sani process is displayed in the Status screen for monitoring purpose.
4. Once sani is complete the system will return to STANDBY mode.
 - The heat sanitization is logged in the HMI as complete.

4.4.3 Loop and Points of use Heat Sanitization

A manual heat sanitization of the loop and points of use can be scheduled during non-usage periods. It is important to leave enough time for the heat sanitization to complete. If pre-treatment devices used have a regeneration or backwash schedule, ensure they do not conflict with the auto heat sanitization schedule unless they are equipped with a lock-out input connected to the panel.

1. From the Mode Control screen, log in with username and password.
2. Press the Standby Mode icon if the system is normally operating.
3. Press the Loop & POU icon under the heading for single heat sanitization and allow the system water temperature to ramp to the set temperature.

4. Once the temperature is reached the PLC will provide an HMI indication that the temperature is reached and hot water can be drawn from the points of use.
 - The system will maintain the temperature and display the tank level in real time throughout the hold period.
 - An alert on the HMI and horn will activate once the tank reaches the adjustable low level setpoint.
 - At the beginning of the mode an adjustable counter will start. Should the counter reach zero before any other alarm, the sanitization sequence will abort and begin a cool down sequence.
5. If the user does not press Sanitization Complete button and the Break Tank level reaches the low low level alarm, the 1st pass RO pump will stop along with the heater.
6. Once the low level alert setpoint is reached the operator shall close all the points of use and press the Sanitization Complete icon. The system will transition to the cool down sequence and the heat sanitization is logged in the HMI as complete.
7. Once the cool down sequence is complete the system will return to Standby mode.
8. The sequence of heat sani process is displayed in the Status screen for monitoring purpose.

4.5 EMERGENCY BY-PASS OPERATION

In an emergency situation the system can be operated in either first pass by-pass or second pass by-pass. When placed in either mode the 2nd pass quality alert setpoint will automatically change to 10 $\mu\text{S}/\text{cm}$ and the 2nd pass quality divert setpoint will automatically change to 20 $\mu\text{S}/\text{cm}$. The reason for this is because in by-pass mode the system operates as a single pass RO producing lower quality water.

NOTE: *An optional by-pass circuit can be added to the panel to allow for manual operation of the system if for some reason the PLC or HMI should fail.*

4.5.1 1st Pass By-Pass Operation

1. Determine and record the issue.
2. Place the system in OFF mode from the Mode Control screen.
3. Close valve V3230 on the inlet to the 1st pass pump P3015.
4. Log in with username and password.
5. Press the Maintenance Mode icon and then press the 1st Pass Bypass icon.
6. Press the Screen Select icon.
7. Press the Status icon.
8. Read the Message and then press the Next icon. The system will automatically restart and operate at a programmed product flow rate. The recovery is automatically adjusted to 50%.
9. Contact Mar Cor service with the issue to arrange a service call.

NOTE: *No heat sanitization can be performed when in emergency by-pass operation.*

CAUTION: *The system should only operate in either by-pass modes for no more than a week.*

WARNING: *The system must be shut down when performing service work on the system in by-pass operation.*

4.5.2 2nd Pass By-Pass Operation

1. Determine and record the issue.
2. Place the system in OFF mode from the Mode Control screen.
3. Close valve V3260 on the inlet to the 2nd pass pump P3025.
4. Log in with username and password.
5. Press the Maintenance Mode icon and then press the 2nd Pass Bypass icon. The system will automatically restart and operate at a programmed product flow rate. The recovery is automatically adjusted to 50%.
6. Press the Screen Select icon.
7. Press the Status icon.
8. Read the Message and then press the Next icon. The system will automatically restart and operate at a programmed product flow rate. The recovery is automatically adjusted to 50%.
9. Contact Mar Cor service with the issue to arrange a service call.

NOTE: *No heat sanitization can be performed when in emergency by-pass operation.*

CAUTION: *The system should only operate in either by-pass modes for no more than a week.*

WARNING: *The system must be shut down when performing service work on the system in by-pass operation.*

4.5.3 Manual By-pass Operation

If the PLC and or HMI should fail the system can be place into manual operation.

1. Push the emergency stop red button on front of panel and wait for 10 seconds then reset the emergency stop push button.
2. Open the top control panel and locate the manual by pass key switch. With the key installed turn the switch to manual. The following will occur.
 - Feed water valve YV3430 will open and 1st pass pump P3015 will start and ramp to the pre-set speed programmed in the variable speed drive.
 - After a time delay the 2nd pass pump P3025 will start and ramp to the pre-set speed programmed in the variable speed drive.
 - Once the conductivity quality setpoint (10µS/cm) is met on the second pass product outlet AE3695 the 2nd pass product divert valve YV3470 will divert to the loop feed. The conductivity setpoint is programmed in the quality meter AIT9600 located inside the top control panel.
3. Monitor feed water pressure and pre-treatment devices during this time to prevent running the pump without water.

4. At the end of the day turn the key switch back to auto to shut the system down.

CAUTION: *The system must not be left in operation during non-dialysis hours without local monitoring. Pre-treatment devices are not interlocked with the system when in by-pass. If any pre-treatment device should perform a backwash or regeneration there is a risk of the pumps running dry causing possible damage.*

Once the issue has been resolved make sure the key switch is in the auto position.

4.5.4 Conductivity Transmitter Manual By-pass Setpoint Change

CAUTION: *The BioPure HX2 system is equipped with a conductivity transmitter and sensors to measure the feed water, 1st pass product and 2nd pass product temperature and conductivity. The transmitter is located inside the top control panel. The procedure below describes how to change the 2nd pass product temperature and conductivity poor quality setpoint if the system is placed in manual by-pass due to a PLC and or HMI system failure.*

NOTE: *The transmitter must be calibrated by an operator that is familiar with accessing and navigating the settings on the transmitter. General navigation and access instruction are not covered in this procedure. Refer to the transmitter operator manual.*


CAUTION: *Do not alter any other settings in the transmitter. System may become inoperable if changes are made.*

Factory password: 3658


Default 2nd Pass Conductivity Manual By-pass Setpoint: 10 μ S/cm.

If the user needs to change the default setpoint the follow procedure will assist them in making the change.

1. Access the parameter and the function menu.



- Navigate to the setpoint value and make the change.



- Save all changes before exiting the menu.

4.6 RO PERFORMANCE OVERVIEW

4.6.1 System Performance

The system flows and pressures are controlled by the PLC and are programmed at the factory level. Depending on the product flow rate and model installed the system product flow output can be increased or decreased. If it is required the pre-treatment capacity will need to be verified prior to making any flow changes. A Mar Cor Technician will need to make the necessary setpoint changes.

The approximate flows are listed in gpm / lpm:

	2nd Pass Product	1st Pass Recovery	1st Pass Reject SP1 3655	1st Pass Reject Drain SP1 3660	2nd Pass Recycle SP1 3690	2nd Pass Product SP1 3710
3 Membrane Nominal 5 gpm/18.9 lpm System	3.0 / 11.4	75%	30.0 / 113.6	1.6 / 6.1	1.8 / 6.8	3.0 / 11.4
		50%	30.0 / 113.6	4.8 / 18.2	1.8 / 6.8	3.0 / 11.4
	3.5 / 13.2	75%	30.0 / 113.6	1.9 / 7.2	2.1 / 7.9	3.5 / 13.2
		50%	25.0 / 94.6	5.6 / 21.2	2.1 / 7.9	3.5 / 13.2
	4.0 / 15.1	75%	30.0 / 113.6	2.1 / 7.9	2.4 / 9.1	4.0 / 15.1
		50%	25.0 / 94.6	6.4 / 24.2	2.4 / 9.1	4.0 / 15.1
	4.5 / 17.0	75%	30.0 / 113.6	2.4 / 9.1	2.7 / 10.2	4.5 / 17.0
		50%	25.0 / 94.6	7.2 / 27.3	2.7 / 10.2	4.5 / 17.0
	5.0 / 18.9	75%	30.0 / 113.6	2.7 / 10.2	3.0 / 11.3	5.0 / 18.9
		50%	25.0 / 94.6	8.0 / 30.3	3.0 / 11.3	5.0 / 18.9
	5.5 / 20.8	75%	30.0 / 113.6	2.9 / 11.0	3.3 / 12.5	5.5 / 20.8
		50%	25.0 / 94.6	8.8 / 33.3	3.3 / 12.5	5.5 / 20.8
	6.0 / 22.7	75%	30.0 / 113.6	3.2 / 12.1	3.6 / 13.6	6.0 / 22.7
		50%	25.0 / 94.6	9.6 / 36.3	3.6 / 13.6	6.0 / 22.7
	6.5 / 24.6	75%	30.0 / 113.6	3.5 / 13.2	3.9 / 14.8	6.5 / 24.6
		50%	25.0 / 94.6	10.4 / 39.4	3.9 / 14.8	6.5 / 24.6

	7.0 / 26.5	75%	30.0 / 113.6	3.3 / 12.5	3.0 / 11.4	7.0 / 26.5
		50%	25.0 / 94.6	10.0 / 37.9	3.0 / 11.4	7.0 / 26.5

	2nd Pass Product	1st Pass Recovery	1st Pass Reject SP1 3655	1st Pass Reject Drain SP1 3660	2nd Pass Recycle SP1 3690	2nd Pass Product SP1 3710
6 Membrane Nominal 11 gpm/14.6 lpm System	7.0 / 26.5	75%	30.0 / 113.6	3.7 / 14.0	3.0 / 11.4	7.0 / 26.5
		50%	20.0 / 75.7	11.2 / 42.4	3.0 / 11.4	7.0 / 26.5
	7.5 / 28.4	75%	30.0 / 113.6	4.0 / 15.1	4.5 / 17.0	7.5 / 28.4
		50%	20.0 / 75.7	12.0 / 45.4	4.5 / 17.0	7.5 / 28.4
	8.0 / 30.3	75%	30.0 / 113.6	4.3 / 16.3	4.8 / 18.2	8.0 / 30.3
		50%	20.0 / 75.7	12.8 / 48.5	4.8 / 18.2	8.0 / 30.3
	8.5 / 32.2	75%	30.0 / 113.6	4.5 / 17.0	5.1 / 19.3	8.5 / 32.2
		50%	20.0 / 75.7	13.6 / 51.5	5.1 / 19.3	8.5 / 32.2
	9.0 / 34.1	75%	30.0 / 113.6	4.8 / 18.2	5.4 / 20.4	9.0 / 34.1
		50%	20.0 / 75.7	14.4 / 54.5	5.4 / 20.4	9.0 / 34.1
	9.5 / 35.9	75%	30.0 / 113.6	5.1 / 19.3	5.7 / 21.6	9.5 / 35.9
		50%	15.0 / 56.7	15.2 / 57.5	5.7 / 21.6	9.5 / 35.9
	10.0 / 37.9	75%	25.0 / 94.6	5.3 / 20.1	6.0 / 22.7	10.0 / 37.9
		50%	15.0 / 56.7	16.0 / 60.6	6.0 / 22.7	10.0 / 37.9
	10.5 / 39.7	75%	25.0 / 94.6	5.6 / 21.2	6.3 / 23.8	10.5 / 39.7
		50%	15.0 / 56.7	16.8 / 63.6	6.3 / 23.8	10.5 / 39.7
	11.0 / 41.6	75%	25.0 / 94.6	5.9 / 22.3	6.6 / 25.0	11.0 / 41.6
		50%	15.0 / 56.7	17.6 / 66.6	6.6 / 25.0	11.0 / 41.6
	11.5 / 43.5	75%	25.0 / 94.6	6.1 / 23.1	6.9 / 26.1	11.5 / 43.5
		50%	15.0 / 56.7	18.4 / 69.7	6.9 / 26.1	11.5 / 43.5
	12.0 / 45.4	75%	25.0 / 94.6	6.4 / 24.2	7.2 / 27.3	12.0 / 45.4
		50%	15.0 / 56.7	19.1 / 72.3	7.2 / 27.3	12.0 / 45.4
	12.5 / 47.3	75%	25.0 / 94.6	6.7 / 25.4	7.5 / 28.4	12.5 / 47.3
		50%	15.0 / 56.7	20.0 / 75.7	7.5 / 28.4	12.5 / 47.3
	13.0 / 49.2	75%	25.0 / 94.6	6.2 / 23.5	5.6 / 21.2	13.0 / 49.2
		50%	15.0 / 56.7	18.6 / 70.4	5.6 / 21.2	13.0 / 49.2
	13.5 / 51.1	75%	25.0 / 94.6	6.4 / 24.2	5.8 / 22.0	13.5 / 51.1
		50%	15.0 / 56.7	19.3 / 73.1	5.8 / 22.0	13.5 / 51.1
	14.0 / 53.0	75%	25.0 / 94.6	6.7 / 25.4	6.0 / 22.7	14.0 / 53.0
		50%	10.0 / 37.9	20.0 / 75.7	6.0 / 22.7	14.0 / 53.0
	14.5 / 54.9	75%	25.0 / 94.6	6.4 / 24.2	4.8 / 18.2	14.5 / 54.9
		50%	15.0 / 56.7	19.3 / 73.1	4.8 / 18.2	14.5 / 54.9
	15.0 / 56.8	75%	25.0 / 94.6	6.7 / 25.4	5.0 / 18.9	15.0 / 56.8
		50%	30.0 / 113.6	3.7 / 14.0	3.0 / 11.4	7.0 / 26.5

4.6.2 Effects of Temperature on System Performance

The design product flow rate for the RO system is based on a feed water temperature of 25°C (77°F). The product flow rate will be reduced at temperatures below 25°C (77°F) as the viscosity of the water increases, making it more difficult to force water through the reverse osmosis membrane.

CAUTION: *At temperatures higher than 25°C (77°F), the operating pressure should be reduced to maintain the design flow rate. Operating at a higher product flow will invalidate the RO membrane warranty.*

If the temperature is not 25°C (77°F), a correction factor must be used to calculate the expected product flow rate from the system. The table below gives temperature correction factors to be applied for feed temperatures other than 25°C (77°F).

Temperature correction factors are listed for all reverse osmosis elements. The reference temperature is 25°C (77°F).	
Temperature °F (°C)	Factor
40 (4)	0.48
50 (10)	0.60
60 (16)	0.73
70 (21)	0.88
77 (25)	1.00
80 (27)	1.06
90 (32)	1.26

CAUTION: *Water quality can be affected when operating at higher temperatures.*

4.6.3 High Recovery Operation

The system will automatically adjust the reject flow to drain. The 1st pass RO pump VFD will decrease or increase the speed of the pump to maintain the flow setpoint.

The PLC will monitor the feed water conductivity and increase or decrease the 1st pass system recovery as the feed water conductivity improves by closing or opening the 1st Pass Reject Drain Control Valve per the chart below. The minimum flow measured at the reject drain flow sensor must be greater than 1.0 gpm.

Feed Water Conductivity (µS/cm)	System Recovery (%)
>300	75%
300	80%
150	85%
100	90%
50	95%

NOTE: *If the system is set to operate at 50% recovery the high recovery operation is disabled.*

5.0 MAINTENANCE

5.1 MAINTENANCE REQUIREMENTS

After the BioPure HX2® system has been commissioned; it will require scheduled maintenance at a defined frequency. Regular maintenance is necessary in order to ensure:

- The product water quality meets or exceeds the ISO 13959-11 water quality for haemodialysis applications.
- Efficiency of RO performance.
- Maximum membrane life.

Section 5.1.3 (log sheet) lists the operating data that needs to be monitored and recorded on the RO system at a specified frequency. Alerts and Alarms are automatically recorded in the PLC history to provide indication when maintenance is required.

All readings and samples must be taken when the system is in normal operation (RO Run To Loop).

Section 5.1.4 includes a list of the recommend service schedule of the system.

A 15 minute water sample alert can be activated on the HMI to provide an audible alert to the nurse / technicians that readings and water samples can be taken. The timer will begin to count once the system has transition to RO Run To Loop mode.

An alarm / alert log in section 5.1.5 is provided to record any alarms / alerts that could occur during maintenance.

All alarms / alerts / messages are recorded in the history on the HMI.

CAUTION: *Any changes in feed water chemistry must be dealt with immediately to prevent system poor performance and or system failure.*

CAUTION: *The log sheet does not cover any parameters that need to be monitored on the pre-treatment or post treatment of the system. Please consult with you Mar Cor representative.*

NOTE: *Maintenance of the BioPure HX2® system following its installation is the responsibility of the end users to ensure that product water is tested periodically. Appropriate hot water sanitization at scheduled frequency and/or use of chemical sanitization must be performed to maintain acceptable product water quality.*

NOTE: *Use this system only with the specified feed water and pre-treatment chemicals determined by a Mar Cor representative. Use all sanitizing and cleaning agents according to the instructions in this Manual.*

5.1.1 RO System Shutdown Procedure

If the RO system is not being used for less than one week, it is usually sufficient to rely on the Auto Flush feature for flushing the system automatically at a pre-set time interval to minimize biological growth. Place the system in Standby to allow the system to perform its auto flush sequence.

If the RO system is not being used for more than one week, it is recommend that the heat sanitization schedule be programmed to perform a heat cycle at a minimum of once a week along with its normal Auto Flush cycles. Place the system in Standby to allow the system to perform its scheduled heat cycles and flush sequences.

If the system must be shut down for an extended time (greater than one week) without water and or power, a preservative solution (1% Sodium Bisulfate) should be added to the system to reduce biological growth. The solution will have to be pumped through the system manually and must be performed by a Mar Cor representative. Membranes shall be flushed with pre-treated water and re-preserved every three months.

WARNING: Placing the system in 'Off' mode for greater than 24 hours can result in excessive bacterial growth. A full system heat sanitization and water quality sampling must be performed prior to use.

5.1.2 RO Pre-filter Cartridge Replacement Procedure

NOTE: *As with any water system, water spillage on the ground or on the equipment when servicing the equipment is expected. Make sure proper equipment in on hand to dry up excess water during maintenance.*

1. Press the Mode Control icon followed by the Off Mode icon.
2. Fully close manual valve (customer installed) to isolate feed water supply.
3. Open vent valve V3210 to depressurize filter chamber.
4. Open sample valve V3205 to drain off water in the housing.
5. Close valve V3205.
6. Unlock the V-clamp to remove filter top cover.
7. Remove the hold-down plate after turning the knob loose. The knob keeps the plate firmly against the filter cartridges.
8. Remove and discard old filter cartridges.
9. Inspect and replace o-ring if signs of cracking or degrading are found on the housing.
10. Wear clean gloves and install new filter cartridges.
11. Return the hold-down plate and make sure it is set flat on the cartridges.
12. Reinstall knob to keep hold-down plate firmly against the cartridges.
13. Install filter housing top cover and V-clamp.
14. Slowly open the manual feed water valve (customer installed).
15. Once a steady stream of water is seen exiting the vent valve V3210, close the valve.
16. Check for leaks and repair if any are found.

5.1.3 Log Sheet

CAUTION: All readings and samples must be taken when the system is in normal operation (RO Run To Loop).
Any changes in feed water chemistry must be dealt with immediately to prevent poor performance and/or system failure.

Parameters	Sample / Reading Location	Test D/M/Y	Units	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
Date			D/M/Y							
Time			HH:MM							
Temperature	TI3640	D	°C / °F							
pH	V3225	D								
Total Chlorine	V3205	D	ppm							
Feed Iron	V3205	M	ppm							
Total Hardness	V3205	D	ppm							
Feed SDI	V3225	Y								
Pre-Filter Inlet Pressure	PI3600	D	psig							
Pre-Filter Outlet Pressure	PI3605	D	psig							
Concentrated Feed Quality	AI3635	D	µS/cm							
Feed Flow	CFI9868	D	lpm/gpm							
1st Pass Pump Pressure	PI3645	D	psig							
1st Pass Reject Pressure	PI3650	D	psig							
1st Pass Reject Flow	FI3655	D	lpm/gpm							
1st Pass Reject Drain Flow	FI3660	D	lpm/gpm							
1st Pass Recycle Flow	CFI9872	D	lpm/gpm							

CAUTION: All readings and samples must be taken when the system is in normal operation (RO Run To Loop).

Any changes in feed water chemistry must be dealt with immediately to prevent poor performance and/or system failure.

Parameters	Sample / Reading Location	Test D/M/Y	Units	Mon	Tue	Wed	Thurs	Fri	Sat	Sun
1st Pass Product Quality	AI3665	D	µS/cm							
2nd Pass Pump Pressure	PI3675	D	psig							
2nd Pass Reject Pressure	PI3685	D	psig							
2nd Pass Recycle Flow	FI3690	D	lpm/gpm							
2nd Pass Product Flow	FI3710	D	lpm/gpm							
2nd Pass Product Quality	AI3695	D	µS/cm							
2nd Pass Product Pressure	PI3705	D	psig							
Loop filter Outlet Pressure	PI3715	D	psig							
Loop Return Flow	FI3620	D	lpm/gpm							
Loop Return Temperature	TI3630	D	°C / °F							
Operator Initials:										
Notes:										

5.1.4 Service Schedule

Item	Tag Reference	Service Interval
Pre-filter *	F3005	Replace every 3 months or Delta P Alert.
		Bleed air from pre-filter using valve V3210 once a week.
Loop filter * (Optional)	F3035	Replace every 6 months or Delta P Alert.
1 st Pass RO Membranes **	F3020	Heat Sanitize once a month to a maximum of once a week. Clean and Chemical Sanitize every 6 months or Delta P Alert or replace as required.
2 nd Pass RO Membranes **	F3030	Heat Sanitize once a month to a maximum of once a week. Clean and Chemical Sanitize every 6 months or Delta P Alert or replace as required.
Solenoid Valves **	YV3400,YV3410	Replace Diaphragm Every 2 years.
Compressor Inlet Filter *	F6000	Inspect every 6 months and replace as required.
Compressor Outlet Filter *	F6030	Purge filter monthly, inspect every 6 months and replace yearly.
Compressor Oil Level *	P6010	Check oil level weekly and replace yearly by a Mar Cor Technician.
Compressor Air Tank *	T6020	Purge tank monthly.
Compressor Relief Valve *	V6230	Test monthly.
Feed Conductivity Sensor **	AE3635	Verify reading every 6 month calibrate once a year.
1 st Pass Conductivity Sensor **	AE3665	Verify reading every 6 month calibrate once a year.
2 nd Pass Conductivity Sensor **	AE3695	Verify reading every 6 month calibrate once a year.
PLC Battery **		Replace when low battery alert shows on HMI banner.
Spring Locking Ring on RO Housing *		Check locking ring is secured by visual inspection weekly. If the locking ring is loose, depressurize the RO housing and make proper adjustment to secure the ring in place.
High Voltage Terminals **		Inspect every 6 months and tighten as required.

* Indicates components can be serviced by End User.

** Indicates components shall be serviced by a Mar Cor Technician or technician with equivalent skills.

5.1.5 Alarm / Alert Log Sheet

Alarm Log		
Time	Alarm Condition	Corrective Measure

5.2 CLEANING OVERVIEW

This overview provides general information about the usual foulants affecting the performance of Composite Polyamide Reverse Osmosis (RO) membrane elements and the removal of these foulants.

NOTE: *The Composite Polyamide type of RO membrane elements may not be exposed to chlorinated water under any circumstances. Any such exposure will cause irreparable damage to the membrane. Absolute care must be taken following any disinfection of piping or equipment or the preparation of cleaning or storage solutions to ensure that no trace of chlorine is present in the feed water to the RO membrane elements.*

NOTE: *It is recommended that all RO membrane cleaning operations should be closely coordinated with Mar Cor during the RO membrane element warranty period. Mar Cor field service personnel are available to be on site for cleaning assistance, should the need arise. Please contact Mar Cor for current charges for this service.*

NOTE: *The use of cationic surfactant should be avoided in cleaning solutions, since irreversible fouling of the membrane elements may occur.*

During normal operation over a period of time, RO membrane elements are subject to fouling by suspended or sparingly soluble materials that may be present in the feed water. Common examples of foulants are:

1. Calcium carbonate scale
2. Sulfate scale of calcium, barium or strontium
3. Metal oxides (iron, manganese, copper, nickel, aluminum, etc.)
4. Polymerized silica scale
5. Inorganic colloidal deposits
6. Mixed inorganic/organic colloidal deposits
7. NOM (Natural Organic Matter)
8. Man-made organic material (e.g. antiscalant/dispersants, cationic polyelectrolytes)
9. Biological (bacterial bioslime, algae, mold, or fungi)

The nature and rapidity of fouling depends on a number of factors, such as the quality of the feed water and the system recovery rate. Typically, fouling is progressive, and if not controlled early, will impair the RO membrane element performance in a relatively short time. Cleaning should occur when the RO shows evidence of fouling, just prior to a long-term shutdown, or as a matter of scheduled routine maintenance. The elements shall be maintained in a clean or “nearly clean” condition to prevent excessive fouling by the foulants listed above. Some fouling is allowed as long as:

1. Normalized permeate flow decrease is less than 10%.
2. Normalized permeate quality decrease is less than 10%.
3. Normalized pressure drop, as measured between the feed and concentrate headers, increase is less than 15%.

Cleaning should be carried out before these values are exceeded to maintain the elements in a clean or “nearly clean” condition. Effective cleaning is evidenced by the return of the normalized parameters to their initial, Start-up, value. In the event you do not normalize your operating data, the above values still apply if you do not have major changes in critical operating parameters. The operating parameters that have to stay constant are permeate flow, permeate back-pressure, recovery, temperature, and feed TDS. If these operating parameters fluctuate, then it is highly recommended that you normalize the data to determine if fouling is occurring or if the RO is actually operating normally based on the change in a critical operating parameter.

Monitoring overall plant performance on a regular basis is an essential step in recognizing when membrane elements are becoming fouled. Performance is affected progressively and in varying degrees, depending on the nature of the foulants.

RO cleaning frequency due to fouling will vary by site. A rough rule of thumb as to an acceptable cleaning frequency is once every 3 to 12 months. If you have to clean more than once a month, you should be able to justify further capital expenditures for improved RO pre-treatment or a re-design of the system, contact Mar Cor. If the cleaning frequency is every one to three months, you may want to focus on improving the operation of your existing equipment but further capital expenditure may be harder to justify.

It is important to clean the membranes when they are only lightly fouled, not heavily fouled. Heavy fouling can impair the effectiveness of the cleaning chemical by impeding the penetration of the chemical deep into the foulant and in the flushing of the foulant out of the elements. If normalized membrane performance drops 30 to 50%, it may be impossible to fully restore the performance back to baseline conditions.

One RO design feature that is commonly overlooked in reducing RO cleaning frequency is the use of RO permeate water for flushing foulants from the system. Soaking the RO elements during standby with permeate can help dissolve scale and loosen precipitates, reducing the frequency of chemical cleaning.

What you clean for can vary from site to site depending on the foulant. Complicating the situation frequently is that more than one foulant can be present, which explains why cleanings frequently require a low pH and high pH cleaning regiment.

NOTE: *The membrane elements shall not be exposed to feed water containing oil, grease, or other foreign matter which proves to chemically or physically damage the integrity of the membrane.*

5.3 SANITIZATION OVERVIEW

The performance requirement of purified water equipment is generally based on measurement of conductivity and total organic carbon (TOC) in the product water. An additional requirement for the medical, biotechnical, and pharmaceutical industries is the monitoring of viable microbes and endotoxins limits. In order to meet these limits, specific sanitization processes are used to

disinfect and maintain the purified water equipment, allowing the product water to consistently meet requirements set by ISO 13959-11 standards. Therefore, from a regulatory perspective, the end user of the purified water equipment is obligated to ensure these processes are set in place in order to ensure these requirements are met. Today, it is common for end user protocols to stipulate very low, viable microbe and endotoxin levels in their product water.

In general, ambient operating temperatures and the absence of chlorine or other disinfectants in the water can create a microbe breeding environment on the surface of RO membranes. To attain better system control over microbe levels, end users have employed routine chemical sanitization processes where the cleaning frequency is often based on the bacterial level measured in the product water.

Chemicals used in RO membrane cleaning procedures include low and high pH-based cleaners for the primary purpose of removing hardness scale and organic deposits that may build up on the membrane surface. Some microbiological control is also achieved as a secondary benefit.

NOTE: *The RO system must be chemical cleaned or chemical sanitized by highly skilled personnel. Contact Mar Cor representative for assistance to perform appropriate cleaning procedures in the field.*

Bacteria counts taken from water streams are typically representative of the “loose” microbes in the stream and not necessarily of microbes that may exist in a biofilm within the water system. Biofilms are a complex organization of microbes formed by one of many common species, which adhere to a surface and cover themselves from the fluid stream in a somewhat protective polysaccharide boundary layer.

Biofilms form naturally but the growth may be mitigated with appropriate fluid velocity, smooth sanitary pipe designs, and periodic and appropriate sanitization processes. However biofilm has proven to be quite resilient to many chemical sanitization processes due to their protective boundary layer. Studies demonstrate biofilms may survive even after 60 minutes of exposure to chlorine-based sanitizing chemicals. Preventing the formation of biofilms, therefore, is a primary goal with respect to “regular and periodic” sanitization processes.

Although many chemical sanitization processes may be somewhat ineffective against biofilms, based on specific chemicals, exposure time and frequency of application, the sub-boundary layer bacterial colony cannot escape exposure to sanitizing temperatures during a hot water sanitization process.

Effective sanitization with hot water is accomplished through an appropriate combination of exposure time and temperature. While not officially established, the industry has “standardized” 80°C (176°F) for 30 to 60 minutes as the optimum sanitization target point.

A primary use for hot water sanitization is to inactivate viable microbes. However, it must be noted that endotoxin reduction is not achieved as a direct result of the hot water sanitization process. Based on the feed water source, system operating conditions and the end user’s preventative maintenance practices, traditional chemical cleaning processes may still be

required. Generally, hot water cleanings are performed more frequently than chemical cleanings, based on microbe counts and/or regular sanitization schedule.

5.4 CHEMICAL SELECTION AND USE OF CLEANING CHEMICALS

There are a number of factors involved in the selection of a suitable cleaning chemical (or chemicals) and proper cleaning protocol. The first time you have to perform a cleaning, it is recommended to contact the manufacturer of the equipment, the RO element manufacturer, or a RO specialty chemical and service supplier. Once the suspected foulant(s) are identified, one or more cleaning chemicals will be recommended. These cleaning chemical(s) can be generic or can be private-labeled proprietary chemicals. Typically, the generic chemicals can be of technical grades and are available from local chemical supply companies. The proprietary RO cleaning chemicals can be more expensive, but may be easier to use and one cannot rule out the advantage of the intellectual knowledge supplied by these companies. Some independent RO service companies can determine the proper chemicals and cleaning protocol for your situation by testing a fouled element pulled from your system at their facility.

It is not unusual to use a number of different cleaning chemicals in a specific sequence to achieve the optimum cleaning. Typically, a high pH cleaning is used first to remove foulants like oil or biological matter, followed by a low pH cleaning to remove foulants like mineral scale or metal oxides/hydroxides fouling. There are times that order of high and low pH cleaning solutions is reversed or one solution only is required to clean the membranes. Some cleaning solutions have detergents added to aid in the removal of heavy biological and organic debris, while others have a chelating agent like EDTA added to aid in the removal of colloidal material, organic and biological material, and sulfate scale. An important thing to remember is that the improper selection of a cleaning chemical, or the sequence of chemical introduction, can make the foulant worse.

Mar Cor recommends that the membrane system operator thoroughly investigate the signs of fouling before they select a cleaning chemical and a cleaning protocol. Some forms of fouling (iron deposits and scaling commonly associated with well waters) may require only a simple low pH cleaning. However, for most complex fouling phenomena, a multi-step solution is often required.

NOTE: *The RO system must be chemical cleaned or chemical sanitized by highly skilled personnel. Contact Mar Cor representative for assistance to perform appropriate cleaning procedures in the field.*

5.5 GENERAL PRECAUTIONS IN CLEANING CHEMICAL SELECTION AND USAGE

NOTE: *When using a proprietary chemical, make sure the chemical has been qualified for use with your membrane by the chemical supplier. The chemical supplier's*

instructions should not be in conflict with the manufacturer's recommended cleaning parameters and limits listed in this overview.

NOTE: *When using generic chemicals, make sure the chemical has been qualified for use with your membrane in this overview.*

NOTE: *Use the least harsh cleaning regiment to get the job done. This includes the cleaning parameters of pH, temperature, and contact time. This will optimize the useful life of the membrane.*

NOTE: *Clean at the recommended target temperatures to optimize cleaning efficiency and membrane life.*

NOTE: *Use the minimum amount of chemical contact time to optimize membrane life.*

NOTE: *Be prudent in the adjustment of pH at the low and high pH range to extend the useful life of the membrane. A "gentle" pH range is 4 to 10, while the harshest is 2 to 4 and 10 to 12.*

NOTE: *Oil and biologically-fouled membranes should not use a low pH clean-up first as the oil and biological matter will congeal.*

NOTE: *Cleaning and flushing flows should be in the same direction as the normal feed flow to avoid potential telescoping and element damage.*

NOTE: *When cleaning a multi-stage RO, the most effective cleaning is one stage at a time so cleaning flow velocities can be optimized and foulants from upstream stages do not have to pass through downstream stages.*

NOTE: *Flushing out detergents with higher pH permeate water can reduce foaming problems.*

NOTE: *Verify that proper disposal requirements for the cleaning solution are followed.*

NOTE: *If the system has been fouled biologically, you may want to consider the extra step of introducing a sanitizing biocide chemical before and after a successful cleaning. Biocides can be introduced before and immediately after cleaning, periodically (e.g. once a week), or continuously during service. User must be sure that the biocide is compatible with the membrane, does not create any health risks, is effective in controlling biological activity, and is not cost prohibitive.*

CAUTION: *For safety reasons, make sure all hoses and piping can handle the temperatures, pressures and pH's encountered during a cleaning.*

CAUTION: *For safety reasons, always add chemicals slowly to an agitated batch of make-up water.*

CAUTION: *For safety reasons, always wear safety glasses and protective gear when working with chemicals.*

CAUTION: *For safety reasons, do not mix acids with caustics. Thoroughly rinse any cleaning solution from the RO system before introducing the next solution.*

CAUTION: *When performing a cleaning/sanitization it is very important to use the proper safety apparatus. Refer to the Material Safety Data Sheet (MSDS) provided with the sanitant or cleaner before using it. Carefully follow the manufacturer's safety instructions on the labels affixed to the chemical containers.*

NOTE: *After each batch of solution is used, the local municipality may require it to be neutralized before discharging. Check the local municipality for waste discharge regulations.*

5.6 CLEANING/SANITIZATION FREQUENCY

It is recommended that a cleaning be performed every 6 month on the first pass membrane or anytime a pre-treatment upset has occurred. A full system (1st pass, 2nd pass, and loop) sanitization with MINNCARE® Disinfectant should also be performed every 6 months.

NOTE: *The RO system must be chemical cleaned or chemical sanitized by highly skilled personnel. Contact Mar Cor representative to provide assistance to perform appropriate cleaning in the field.*

5.7 RO CLEANING/SANITIZATION CHEMICALS

The following cleaners have been evaluated and found to be compatible with the membrane elements. In addition to the generic cleaners listed below, there are many chemical companies that market membrane cleaners. For best results, use Mar Cor's recommended brands for antiscalants and cleaners for your specific elements.

Type of Foulant	Chemical
Mineral scale and metal precipitates	HCl or citric acid, pH 2
Organics, silt, bacterial slime	Sodium hydroxide, pH 12.5
Bacteria	MINNCARE® Disinfectant, 1%

5.8 RO CLEANING/SANITIZATION CHEMICALS CONCENTRATIONS

The following cleaners have been evaluated and found to be compatible with the membrane element.

Chemical Used	Quantity When Diluted In 189.3 Litres (50 Gallons) Of RO Quality Water Or Better
12N Hydrochloric Acid	1 litre (0.25 gallons), pH 2.5

Sodium Hydroxide Pellets	0.2 Kg (0.42 lbs), pH 12.5
MINNCARE® Disinfectant	2 litres (0.5 gallons), 1%

Chemical Used	Residual Concentration After Flushing
12N Hydrochloric Acid	pH 6 to 8
Sodium Hydroxide Pellets	pH 6 to 8
MINNCARE® Disinfectant	MINNCARE® Residual Test Strips

CAUTION: *Read all chemical safety labels prior to use.*

CAUTION: *Pretreated quality water or better must be used when rinsing and preparing the cleaning/sanitization solutions. Failure to do so will result in damage to the membrane elements.*

CAUTION: *When performing a sanitization using the MINNCARE® disinfectant, an acid cleaning must be performed first. This will remove any metal build-up on the membranes and prevent the membranes from oxidization.*

NOTE: *When cleaning the system, ensure it is operating at high flow rates and low pressure. If high pressure >60 psig is used, contaminants can be pushed further inside the membranes.*

CAUTION: *Do Not Breathe chemicals. Perform mixing operations only in well ventilated areas. If spattered with chemicals, remove contaminated clothing and immediately flush the areas with water for 20 minutes, then notify your supervisor. Wear rubber gloves, a rubber apron and face shield when preparing solutions and moving hoses with chemicals still present.*

5.9 NEUTRALIZING CLEANING/SANITIZATION SOLUTIONS

After each batch of solution is used, the local municipality may require it to be neutralized for discharging. Each batch of solution will have a different concentration after each use. Therefore it will have to be calculated at each time. Refer to the accompanying product brochure of the chemical used for instructions on neutralizing the solution.

NOTE: *Check with the local municipality for waste discharge regulations.*

5.10 MINNCARE® DISINFECTANT USE

MINNCARE® Liquid Disinfectant may be used for reverse osmosis (RO) systems, which are compatible with diluted hydrogen peroxide solutions. The MINNCARE® Liquid Disinfectant has a shelf life of one year. MINNCARE® Liquid Disinfectant should not be stored at its 1% use dilution

since this may compromise its effective concentration. This product is not registered for use on kidney dialyzers and dialysis machines.

The RO manufacturer should be consulted prior to use of MINNCARE[®] Liquid Disinfectant to determine the temperature and pH range acceptable for the particular membranes. At a 100X dilution (1% concentration), the pH range for MINNCARE[®] Liquid Disinfectant is 3.0-3.5. RO water systems vary in design according to the particular needs of the user. Due to the variations in materials of construction at different facilities, some systems may have components made of materials of construction such as copper, brass, iron and other metals that are not compatible with long term exposure to MINNCARE[®] Liquid Disinfectant. Latex and Buna-N have also exhibited a decreased resistance to extended exposure to MINNCARE[®] Liquid Disinfectant solution. More typical system materials have been found to be compatible with MINNCARE[®] Liquid Disinfectant solution such as stainless steel, polypropylene, high-density polyethylene, polysulfone, PTFE, polycarbonate, neoprene, ABS, nylon, acrylic silicone, Plexiglas, ethylene propylene, VITON, FLUOREL, PVC and CPVC.

Due to the uniqueness of RO water systems in terms of design and construction as well as the quality of raw feed water, all or portions of the following procedure may be used in the sanitization process.

Biological or organic fouling of the membrane or other parts of the system should be removed with appropriate cleaner. It is important to follow the membrane manufacturer's recommended cleaning procedure. After cleaning, flush the system with RO permeate. Mineral deposits should be removed with an acidic cleaner prior to sanitizing of the membrane. Again, follow the membrane manufacturer's recommended cleaning procedure. After cleaning, flush the system with RO permeate. The presence of iron or other transition metals, in conjunction with the hydrogen peroxide in MINNCARE[®] Liquid Disinfectant, could cause membrane degradation. Prepare a 1% solution (1 part MINNCARE[®] Liquid Disinfectant to 99 parts water) by adding the MINNCARE[®] Liquid Disinfection solution to permeate water. Fill the entire water circuit to be sanitized with a 1% solution and allow the diluted solution to reach a minimum temperature of 20°C (68°F). Do not exceed the membrane manufacturer's recommended temperature. Recirculate the 1% MINNCARE[®] Liquid Disinfectant until the entire system is filled. Allow the elements to soak in the 1% MINNCARE[®] Liquid Disinfection solution for a minimum of 36 minutes at 20°C (68°F).

Rinse the RO system and check for residuals by following the directions on the MINNCARE[®] Residual Test Strips label. The residual test strip should indicate less than 2ppm. Rinse times will vary depending on the size of the RO system. Residual sanitizer that may enter the system due to chemical rebound can be eliminated by diverting product water to drain for a short period of time.

CAUTION: *When performing a sanitization using the MINNCARE[®] disinfectant, an acid cleaning must be performed first. This will remove any metal build-up on the membranes and prevent the membranes from oxidization.*

MINNCARE® Liquid Disinfectant may be used for reverse osmosis (RO) systems, which are compatible with diluted hydrogen peroxide solutions. The MINNCARE® Liquid Disinfectant has a shelf life of one year. MINNCARE® Liquid Disinfectant should not be stored at its 1% use dilution since this may compromise its effective concentration. This product is not registered for use on kidney dialyzers and dialysis machines.

The RO manufacturer should be consulted prior to use of MINNCARE® Liquid Disinfectant to determine the temperature and pH range acceptable for the particular membranes. At a 100X dilution (1% concentration), the pH range for MINNCARE® Liquid Disinfectant is 3.0-3.5. RO water systems vary in design according to the particular needs of the user. Due to the variations in materials of construction at different facilities, some systems may have components made of materials of construction such as copper, brass, iron and other metals that are not compatible with long term exposure to MINNCARE® Liquid Disinfectant. Latex and Buna-N have also exhibited a decreased resistance to extended exposure to MINNCARE® Liquid Disinfectant solution. More typical system materials have been found to be compatible with MINNCARE® Liquid Disinfectant solution such as stainless steel, polypropylene, high-density polyethylene, polysulfone, PTFE, polycarbonate, neoprene, ABS, nylon, acrylic silicone, Plexiglas, ethylene propylene, VITON, FLUOREL, PVC and CPVC.

Due to the uniqueness of RO water systems in terms of design and construction as well as the quality of raw feed water, all or portions of the following procedure may be used in the sanitization process.

Biological or organic fouling of the membrane or other parts of the system should be removed with appropriate cleaner. It is important to follow the membrane manufacturer's recommended cleaning procedure. After cleaning, flush the system with RO permeate. Mineral deposits should be removed with an acidic cleaner prior to sanitizing of the membrane. Again, follow the membrane manufacturer's recommended cleaning procedure. After cleaning, flush the system with RO permeate. The presence of iron or other transition metals, in conjunction with the hydrogen peroxide in MINNCARE® Liquid Disinfectant, could cause membrane degradation. Prepare a 1% solution (1 part MINNCARE® Liquid Disinfectant to 99 parts water) by adding the MINNCARE® Liquid Disinfection solution to permeate water. Fill the entire water circuit to be sanitized with a 1% solution and allow the diluted solution to reach a minimum temperature of 20°C (68°F). Do not exceed the membrane manufacturer's recommended temperature. Recirculate the 1% MINNCARE® Liquid Disinfectant until the entire system is filled. Allow the elements to soak in the 1% MINNCARE® Liquid Disinfection solution for a minimum of 36 minutes at 20°C (68°F).

Rinse the RO system and check for residuals by following the directions on the MINNCARE® Residual Test Strips label. The residual test strip should indicate less than 2ppm. Rinse times will vary depending on the size of the RO system. Residual sanitizer that may enter the system due to chemical rebound can be eliminated by diverting product water to drain for a short period of time.

CAUTION: *When performing a sanitization using the MINNCARE® disinfectant, an acid cleaning must be performed first. This will remove any metal build-up on the membranes and prevent the membranes from oxidization.*

5.10.1 MINNCARE® Test Strip Overview

MINNCARE® Test Strips provide quick results with easy-to-read indicators. Use the MINNCARE® 1% TS test strips as a pass/fail measurement for adequate concentration of sterilant after dilution. Then, verify residual levels after rinse-out using MINNCARE® Residual Test Strips.

The test strips verify adequate concentration of a 1% MINNCARE® solution during sanitization of your high purity water system. These convenient dip-and-read strips can be used at any test port in your system to provide an accurate pass/fail measurement of MINNCARE® solution. The Test Strips are easy-to-use and are conveniently packaged in 100 test strips per container. Each strip has a “reaction zone” pad located at one end of the test strip. When stored and used properly, the test strip pad can indicate residual MINNCARE® levels as low as 1 part per million (ppm). The reaction zones are set to provide six different reading levels at 100, 30, 10, 3, 1 and 0 (zero) parts per million.

CAUTION: *Store sealed vials in a cool dry environment. It is not necessary to refrigerate. Once open store at room temperature. DO NOT REFRIGERATE. Failure to store tests properly will result in inaccurate readings. Close vial tightly when not in use.*

NOTE: *Each container has an expiry date. Do not use after date printed on the container.*

5.10.2 Test Strip Instruction

MINNCARE® P/N 78339 Test Strips

1. Collect a 1ml or greater sample of the MINNCARE® solution.
2. Remove a test strip from the tube and immediately replace the lid.
3. Dip the test strip into the solution such that the reaction zone is properly wetted.
4. Remove the test strip and gently shake off excess liquid.
5. Within 10 seconds, note the entire reaction pad colour development.
 - 1% = Blue-grey or blue-black colour over the entire reaction pad within 2-5 seconds and which does not fade, indicates MINNCARE® levels of 1% or greater.
 - No colour development or immediate fading of blue-grey or blue-black indicates MINNCARE® levels of less than 1%.
 - 1% MINNCARE® = 400 ppm Peracetic Acid, test strips is accurate +/- 50 ppm.

NOTE: *The reaction pad may slowly change/develop colour after the initial reaction. This is normal and should be ignored for the purpose of test interpretation.*

MINNCARE® P/N 78338 RESIDUAL TEST STRIPS

1. Collect a sample of the rinse solution.
2. Remove a test strip from the tube and immediately replace the lid.
3. Dip the test strip into the solution to be tested for one second, such that the reaction zone is properly wetted.
4. Remove the test strip, gently shake off excess liquid, and compare the reaction zone with the colour scale between 5 and 10 seconds after exposure.

NOTE: *If test strip indicates dark blue to brown or green to brown, the concentrations are too high for the colour scale. Read the MINNCARE® disinfectant package insert instruction for additional information.*

5.11 SEMI-AUTOMATIC CHEMICAL CLEANING AND SANITIZATION

The chemical cleaning for the BioPure HX2® system is executed via the HMI display by manually going through step-by-step procedures and instructions. The operator executing the chemical cleaning operation will be prompted to manipulate controls, and start and stop the operation of automatic components on the system throughout the entire chemical cleaning operation.

The system must be in 'Off' mode to manually initiate a Chemical Clean. The process is semi-automatic requiring the user to manually open and close valves and add chemicals. This sequence is password protected requiring a Security level Supervisor or higher. Once accessed the user can choose what part of the system will be cleaned or sanitized:

- RO 1st Pass
- RO 2nd Pass
- Loop

NOTE: *The typical chemical clean frequency is once every six months to a maximum of every three months. If more frequent sanitizations are required it will affect the normal life span or calibration frequency of critical components on the system.*

The system is out of service during chemical clean.

Upon equipment selection and initiation of a chemical clean the user will have access to a screen with the following information:

- Start
- Tank Fill
- Line Rinse
- Add Chemical
- Mix (start/stop)
- Recycle (start/stop)

- Neutralize (start/stop)
- Tank Drain
- Flush (start/stop)
- RO Run To Drain (start/stop)
- RO Run to Loop Drain (start/stop)
- Complete
- Recycle Timer 0-120 minutes (default 30)
- Timer reset

The Chemical Clean process mode shall be displayed and the last chemical clean shall be recorded in the history.

XXXXXX Chemical Clean Completed
Started HH:MM:SS YYYY/MM/DD
Completed HH:MM:SS YYYY/MM/DD

5.11.1 Chemical Cleaning and Sanitization Preparation

Prior to performing a chemical cleaning or sanitization verify the following is completed.

1. Chemical quantities and types required are on hand and expiry dates are checked.
2. All chemical warnings and precautions labels are read and understood.
3. Appropriate personal safety equipment is on hand and worn during the execution of the procedures.
4. Flow rates, pressures and water quality readings are recorded to determine if the cleaning has improved the system performance.
5. Set the break tank level setpoint to the desired level. Add 37.9 litres (10 gallons) more to the setpoint to compensate for the line rinse sequence.
6. Set the cleaning recycle timer to the desired time.

5.12 1ST PASS CHEMICAL CLEANING / SANITIZATION PROCEDURE

Security level Supervisor or higher is required to perform the following procedure. The HMI has built-in delays between each icon push to allow the system to transition smoothly. The time delay is shown on the status screen.

1. Press the Mode Control icon.
2. Place the system in OFF mode.
3. Press the Clean Mode icon. Press the 1st Pass icon.
4. Press the Screen Select icon.
5. Press the Status icon.
6. Press the Tank Fill icon:
 - a. Wait for the flows to equalize.
 - b. Wait for the flush duration finish.

- c. 2nd Pass quality check.
7. Wait for tank to fill to 35 gal.
8. Press the Line Rinse icon:
 - a. Tank level drains from 35 gal to 25 gal.
9. Press the Next icon.
10. Remove the 6" cap on top of the break tank and add the desired chemical. Then close the cap when finished.

CAUTION: If performing a sanitization using the MINNCARE® disinfectant, an acid cleaning must be performed first. This will remove any metal build-up on the membranes and prevent oxidization on the membranes.

11. Press the Next icon.
12. Press the Chemical Mix followed by the Start Pump icon. After 5 minutes test the solution concentration at sample valve V3205:
 - a. If concentration is not within limits press the Stop Pump Icon and add more chemical.
 - b. Press the Start Pump icon to continue mixing and check the solution concentration again at sample valve V3205.
13. Allow the chemical mix timer to complete.
14. Press the Next icon.
15. Close the 2nd pass pump P3025 bypass valve V3260.
16. Press the Next icon.
17. Press the 1st Pass Recycle icon followed by the Start Pump icon. After 5 minutes test the solution concentration at sample valve V3280:
 - a. If concentration is not within limits press the Stop Pump icon and more chemical.
 - b. Press the Start Pump icon to continue recycling and check the solution concentration after 5-10 minutes at sample valve V3280.
18. The pump will stop once the cleaning counter elapses or the Pump Stop icon is pressed.
19. Press the Next icon.
20. Remove the 6" cap on top of the break tank and add chemical to the break tank to neutralize the solution. Close the cap when finished.
21. Press the Next icon.
22. Press the Neutralization icon. After 5 minutes test the solution concentration at sample valve V3205. Allow time run out:
 - a. If concentration is not within limits press the Stop Pump icon and add more chemical.
 - b. Press the Start Pump icon to continue neutralize and check the solution concentration again at sample valve V3205.
23. Press the Next icon to confirm neutralization complete.
24. Press the Next icon to start tank drain. Allow drain cycle to complete.
25. Press the Next icon followed by the Flush icon.
26. Press the Screen Select icon followed by the Overview icon. Monitor the product water quality AE3695. Once water quality begins to improve and stabilize, test for residual solution at sample valve V3280. If residual is not within limits allow the system to flush and re-check for residual in another 5 minutes.
27. Press the Screen Select icon.
28. Press the Status icon.
29. Press the Stop Pump icon once the residual is within limits.

30. Press the Next icon and open the 2nd pass pump P3025 bypass valve V3260.
31. Press the Next icon and wait for 5 minutes to allow system to settle. This will help dissipate any hide out pockets of solution.
32. Press the RO Run To Drain icon. Then wait for:
 - a. The flows balance out.
 - b. Flush duration finish.
 - c. RO 2nd pass quality check.
33. Test for residual solution at sample valve V3280. If residual is not within limits, allow the system to run to drain and re-check for residual in another 5 minutes. Continue to run to drain until no residual is detected at sample valve V3280.
34. Press the Clean Complete icon.
35. Place the system back into normal operation or perform another cleaning / sanitization.

NOTE: *After 15 minutes record the flow rates, pressures and water quality reading and determine if the cleaning has improved the system performance. If the performance has not improved, further cleanings / sanitizations could be required or membrane replacement.*

36. The pre-filters F3005 are to be replaced after a cleaning is performed. Refer to section 5.15.1 for instructions.

NOTE: *Once the system cleaning or sanitization is completed rinse the inside walls of the tank to flush away any residue using clean water.*

5.13 2nd PASS CHEMICAL CLEANING / SANITIZATION PROCEDURE

Security level Supervisor or higher is required to perform the following procedure. The HMI has built in delays between each icon push to allow the system to transition smoothly. The time delay is shown on the Status screen

1. Press the Mode Control icon.
2. Place the system in OFF mode.
3. Please the Clean Mode icon.
4. Press the 2nd Pass icon.
5. Press the Screen Select icon.
6. Press the Status icon.
7. Press the Tank Fill icon:
 - a. Wait for the flows to equalize.
 - b. Wait for the flush duration finish.
 - c. 2nd Pass quality check.
8. Wait for tank to fill to 35 gal.
9. Press the Line Rinse icon:
 - a. Tank level drains from 35 gal to 25 gal.
10. Press the Next icon.

11. Remove the 6" cap on top of the break tank and add the desired chemical. Then close the cap when finished.

CAUTION: *If performing a sanitization using the MINNCARE® disinfectant, an acid cleaning must be performed first. This will remove any metal build-up on the membranes and prevent oxidization on the membranes.*

12. Press the Next icon.
13. Press the Chemical Mix followed by the Start Pump icon. After 5 minutes test the solution concentration at sample valve V3205:
 - a. If concentration is not within limits press the Stop Pump Icon and add more chemical.
 - b. Press the Start Pump icon to continue mixing and check the solution concentration again at sample valve V3205.
14. Allow the chemical mix timer to complete.
15. Press the Next icon.
16. Close the 1st pass feed valve V3230.
17. Press the Next icon.
18. Acknowledge that the valve 3230 is closed by pressing the Next icon.
19. Press the 2nd Pass Recycle icon followed by the Start Pump icon. After 5 minutes test the solution concentration at sample valve V3280:
 - a. If concentration is not within limits press the Stop Pump icon and add more chemical.
 - b. Press the Start Pump icon to continue recycling and check the solution concentration after 5-10 minutes at sample valve V3280.
20. The pump will stop once the cleaning counter elapses or the Stop Pump icon is pressed.
21. Press the Next icon.
22. Remove the 6" cap on top of the break tank and add chemical to the break tank to neutralize the solution. Close the cap when finished.
23. Press the Next icon.
24. Press the Neutralization icon.
25. Press the Start Pump icon.
26. After 5 minutes test the solution concentration at sample valve V3205. Allow time run out:
 - a. If concentration is not within limits press the Stop Pump icon and add more chemical.
 - b. Press the Start Pump icon to continue neutralize and check the solution concentration again at sample valve V3205.
27. Press the Next icon to confirm neutralization complete.
28. Press the Next icon to start tank drain. Allow drain cycle to complete.
29. Press the Next icon followed by the Flush icon.
30. Press the Screen Select icon followed by the Overview icon. Monitor the product water quality AE3695. Once water quality begins to improve and stabilize, test for residual solution at sample valve V3280. If residual is not within limits allow the system to flush and re-check for residual in another 5 minutes.
31. Press the Screen Select icon.
32. Press the Status icon.
33. Press the Stop Pump icon once the residual is within limits.
34. Press the Next icon and open the 2nd pass pump P3025 bypass valve V3260.
35. Press the Next icon and wait for 5 minutes to allow system to settle. This will help dissipate any hide out pockets of solution.

36. Press the RO Run To Drain icon. Then wait for:
 - a. The flows balance out.
 - b. Flush duration finish.
 - c. RO 2nd pass quality check.
37. Test for residual solution at sample valve V3280. If residual is not within limits, allow the system to run to drain and re-check for residual in another 5 minutes. Continue to run to drain until no residual is detected at sample valve V3280.
38. Press the Clean Complete icon.

NOTE: *After 15 minutes record the flow rates, pressures and water quality reading and determine if the cleaning has improved the system performance. If the performance has not improved further cleanings / sanitizations could be required or membrane replacement.*

39. The pre-filters F3005 are to be replaced after a cleaning is performed. Refer to section 5.15.1 for instructions.
40. Press Screen Select icon.
41. Place the system back into normal operation or perform another cleaning / sanitization.

NOTE: *Once the system cleaning or sanitization is completed rinse the inside walls of the tank to flush away any residue using clean water.*

5.14 LOOP CHEMICAL CLEANING / SANITIZATION PROCEDURE

Security level Supervisor or higher is required to perform the following procedure. The HMI has built in delays between each icon push to allow the system to transition smoothly. The time delay is shown on the Status screen.

1. Place the system in 'Off' mode.
2. Press the Mode Control icon.
3. Press the Loop icon under the clean mode heading.
4. Press the Screen Select icon.
5. Press the Status icon.
6. Press the Tank Fill icon and allow the tank to fill:
 - a. Wait for the flows to equalize.
 - b. Wait for the flush duration finish.
 - c. 2nd Pass quality check.
7. Wait for tank to fill to 35 gal.
8. Press the Line Rinse icon.
9. Tank level drains from 35 gal to 25 gal.
10. Press the Next Icon.
11. Remove the 6" cap on top of the break tank and add the desired chemical. Then close the cap when finished.
12. Press the Next icon.
13. Press the Chemical Mix followed by the Start Pump icon. After 5 minutes test the solution concentration at sample valve V3205:
 - a. If concentration is not within limits press the Stop Pump Icon and add more chemical.

- b. Press the Start icon to continue mixing and check the solution concentration again at sample valve V3205.
14. Allow the chemical mix timer to complete.
15. Press the Next icon.
16. Press the Next icon.
17. Press the Loop Recycle followed by the Pump Start icon. After 10 minutes test the solution concentration at sample valve V3215:
 - a. If concentration is not within limits press the Stop Pump Icon and add more chemical.
 - b. Press the Start Pump icon to continue recycling and check the solution concentration again after 5-10 minutes at sample valve V3215.
 - c. Draw off a small amount of chemical solution at each point of use to make sure the solution makes contact with each point.
18. The pump will stop once the cleaning counter elapses or Stop Pump icon is pressed.
19. Press the Next icon.
20. Remove the 6" cap on top of the break tank and add the desired chemical. Then close the cap when finished.
21. Press Next icon.
22. Press the Loop Neutralization icon. After 10 minutes test the solution concentration at sample valve V3205:
 - a. If concentration is not within limits press the Stop Pump icon and more chemical.
 - b. Press the Start icon to continue neutralize and check the solution concentration again at sample valve V3205.
23. Press the Stop Pump icon if the solution concentration is within limits.
24. Press the Next icon.
25. Press the Next icon followed by the Tank Drain icon and allow the sequence to complete.
26. Press the Next icon followed by the Loop Flush icon.

NOTE: *There is residual Minncare® solution in the feed water line of the system, and the system will indicate a bad quality alert while the Minncare® solution is flushing to drain. It will take approximately 10-15 minute before the water quality drops below setpoint and diverts to flush the loop.*

27. After 15 minutes test for residual solution at sample valve V3215. Flush each point of use to drain and test for residual solution at those points. If residual is not within limits allow the system to flush and re-check for residual in another 5 minutes.
28. Press the Stop Pump icon once the residual is within limits.
29. Press the Clean Complete icon to finish the cleaning. At this point the system can be placed back in operation or another cleaning / sanitization can be completed.
30. The pre-filters F3005 are to be replaced after a cleaning is performed. Refer to section 5.15.1 for instructions.
31. Press the Screen Select icon.
32. Press the Mode Control icon.
33. Place the system in normal operation.

5.15 SERVICE / REPLACEMENT PROCEDURES FOR SYSTEM COMPONENTS

The BioPure HX2® system uses a number of electrical/electronic devices, manual and pneumatic control valves as identified below. Electrical shocks and high water pressure in the system can cause bodily injuries. Do not attempt to detach or unscrew any lines or valves in the system that are under pressure, and always be sure to place the system to OFF mode. If performing electrical work, switch off the main power disconnect and UPS before working on the system.

NOTE: *It is strongly recommended to have on hand spare O-ring, gaskets and clamps when servicing components on the system.*

NOTE: *It is recommended that any time the system piping has been open to atmosphere that a heat sanitization be performed.*

WARNING: *To prevent electrical shock, disconnect the voltage supply to the system before servicing.*

WARNING: *Use proper safety procedures when working in an electrical panel. Disconnect and/or lock out the power supply to the system and tag it out. The system must be properly grounded to prevent personal injury or damage to the system.*

WARNING: *Drain all the water out in the system and relieve pressure before servicing.*

CAUTION: *Do not energize a solenoid coil with the coil removed from the valve stem. Damage to the coil may occur.*

Below are components that can be found in the BioPure HX2® System. Each component or its subcomponent has a service life that requires replacement:

- 5.15.1 Pre-filter
- 5.15.2 Loop filter
- 5.15.3 Compressor
- 5.15.4 Conductivity Sensor Calibration

5.15.1 RO Pre-filter cartridge replacement Procedure

NOTE: *As with any water system, water spillage on the ground or on the equipment when servicing the equipment is expected. Make sure proper equipment in on hand to dry up excess water during maintenance.*

1. Press the Mode Control icon followed by the Off Mode icon.
2. Fully close manual valve (customer installed) to isolate feed water supply.
3. Open vent valve V3210 to depressurize filter chamber.
4. Open sample valve V3205 to drain off water in the housing.
5. Close valve V3205.
6. Unlock the V-clamp to remove filter top cover.

7. Remove the hold-down plate after turning the knob loose. The knob keeps the plate firmly against the filter cartridges. See [Figure 5.1](#).
8. Remove and discard old filter cartridges.
9. Inspect and replace o-ring if signs of cracking or degrading are found on the housing.
10. Wear clean gloves and install new filter cartridges.
11. Return the hold-down plate and make sure it is sitting flat on the cartridges.
12. Reinstall knob to keep hold-down plate firmly against the cartridges.
13. Install filter housing top cover and V-clamp.
14. Slowly open the manual feed water valve (customer installed).
15. Once a steady stream of water is seen exiting the vent valve V3210, close the valve.
16. Check for leaks and repair if any are found.




FIGURE 5-1

5.15.2 Loop filter Cartridge Replacement (Optional Component)

Security level Supervisor or higher is required to perform the following procedure.

NOTE: *As with any water system, water spillage on the ground or on the equipment when servicing the equipment is expected. Make sure proper equipment in on hand to dry up excess water during maintenance.*

1. Press the Mode Control icon followed by the Off Mode icon.
2. Fully close valves V3285 and V3300.
3. Open vent valve V3290 to depressurize the filter chamber.
4. Open sample valve V3295 to drain water from the housing.
5. Close valve V3295.
6. Unlock the V-clamp to remove filter top cover.
7. Remove the hold-down plate after turning the knob loose. The knob is for keeping the plate firmly against the filter cartridges. See [Figure 5.1](#).
8. Remove and discard old filters.
9. Inspect and replace o-ring if signs of cracking or degrading are found on the housing.

10. Wear clean gloves and install new filter cartridges.
11. Return the hold-down plate and make sure it is sitting flat on the cartridges.
12. Reinstall knob to keep hold-down plate firmly against the cartridges.
13. Install filter housing top cover and V-clamp.
14. Open valves V3285 and V3300.
15. Press the Maint Mode icon followed by the Loop Flush icon and allow the system to start-up.
16. Once a steady stream of water is seen exiting the vent valve V3290, close the valve.
17. Check for leaks and repair if any are found.
18. Press the Off Mode icon and then the Auto Mode icon followed by the Production icon.
19. Perform a loop heat sanitization. Refer to section 4.4


5.15.3 Compressor

The BioPure HX2 system is equipped with an air compressor that provides compressed air to the control valve and actuated valves. Periodic maintenance is required.


Purging the Compressor System of condensate:

This shall be performed once a month

1. Press the Mode Control icon followed by the Off Mode icon.
2. Locate the outlet filter F6030 and press the condensate drain push button valve V6260 to drain out excess condensate.



3. Turn the compressor switch located on top of the tank to the off position.
4. Locate the tank T6020 purge valve V6275. Connect a small tube to the outlet and direct the other end of the tube to a small plastic container.



5. While holding the tube end slowly open the purge valve to drain out the condensate from the tank.


WARNING: *If an oil/water mixture is observed coming out of the tank, it could possibly be a sign the pump is cycling too much. Check for air leaks in the lines, verify the tank is drained of condensate and make sure the correct oil is used. If the condition persists the compressor could be damaged and passing oil.*

6. Once drained close the purge valve.
7. Turn the compressor switch located on top of the tank to the on position and allow the compressor to fill the tank and stop. Check for air leaks.
8. Place the system in the mode it was in prior to servicing.

Relief Valve Test

Check every month

1. Locate the relief valve V6275 on the compressor tank.




2. Slowly pull the ring for 1 second to relieve some pressure.
3. Repeat two times.
4. Check that the valve does not leak. If so replace.

Compressor Inlet Filter Maintenance

Check every six months

1. Locate the black filter F6000 on the compressor motor.




2. Unscrew the filter off the compressor motor.
3. Inspect the inlet and outlet of the filter for signs of debris. If any is found replace the filter.
4. If no debris is found or a new filter is required screw it into the compressor motor.

Compressor Outlet Filter Maintenance

Check every six month and replace yearly

1. Press the Mode Control icon followed by the Off Mode icon.
2. Turn the compressor switch located on top of the tank to the off position.
3. Locate the tank T6020 purge valve V6275. Connect a small tube to the outlet and direct the other end of the tube to a small plastic container.




4. While holding the tube end slowly open the valve to drain out the condensate from the tank.

WARNING: *If an oil/water mixture is observed coming out of the tank, it could possibly be a sign the pump is cycling too much. Check for air leaks in the lines, verify the tank*

is drained of condensate and make sure the correct oil is used. If the condition persists the compressor could be damaged and passing oil.

5. Close the purge valve V6275.
6. Unscrew the filter bowl F6030. Be careful not to lose the o-ring.



7. Look at the filter for signs of contamination without touching it. Filter should be white. If signs of contamination unscrew the filter and replace with new one.
8. Screw in the filter bowl including the o-ring.
9. Turn the compressor switch located on top of the tank to the on position and allow the compressor to fill the tank and stop. Check for air leaks.
10. Place the system in the mode it was in prior to servicing.

Ref. Werther P15-TC Man Rev 1.

5.15.4 Conductivity Transmitter Calibration Procedure

The BioPure HX2 system is equipped with a conductivity transmitter and sensors to measure the feed water, 1st pass product and 2nd pass product temperature and conductivity. The transmitter is located inside the top control panel. The procedure below describes the calibration of the conductivity meter. For maximum measurement accuracy the transmitter should be calibrated once every year.

NOTE: *The transmitter must be calibrated by a Mar Cor Technician that is familiar with accessing and navigating the settings on the transmitter. General navigation and access instruction are not covered in this procedure. Refer to the transmitter operator manual.*

NOTE: *The sensors can be calibrated but only at the manufacturer. If sensor calibration is required it is recommended that spare sensors are purchased. Calibration certificates typically expire after one year so it is recommended that the spare sensor get calibrated close to the time when the sensor on the system calibration expires.*

NOTE: *Since the conductivity sensors used have a PT1000 RTD they do not need to be calibrated.*

The transmitter is password protected to prevent unauthorized access.

Factory password: 3658

Calibration solutions:

Certified NIST calibration solutions shall be required to perform the calibration. The solution chosen should be close to the current conductivity reading at the measurement point in the system during normal system operation or 65% or the full scale. Test solution shall be 25°C +/- 2°C (77°F +/- 3°F) when calibrating.

AE3635 Feed Conductivity, Expected Range is: 500 to 3,000 µS/cm.

AE3665 1st Pass Conductivity, Expected Range is: 1.0 to 50.0 µS/cm.

AE3695 2nd Pass Conductivity, Expected Range is: 0.5 to 10 µS/cm.

Other Equipment Required:

Qty: 3, 250ml plastic clean container

Qty: 1, 500ml deionized water

Qty: 1, 1 ½" tri-clamp stainless cap

Qty: 3, 1 ½" EPDM tri-clamp gaskets

Qty: 1, Step ladder


Qty: 1, Calibration Sticker

Qty: 1, Calibrated hand held
conductivity meter


Calibration Procedure:

NOTE: *As with any water system, there is always water spillage on the ground or on the equipment when servicing the equipment. Make sure proper equipment in on hand to dry up excess water during maintenance*


1. Press the Mode Control icon followed by the Off Mode Icon.
2. See below:
 - a. Close valve V3230 if calibrating feed conductivity sensor AE3635.
 - b. Close valve V3260 if calibrating 1st pass conductivity sensor AE3665.
 - c. No valve to be closed of calibrating 2nd pass conductivity sensor AE3695.
3. Open the top control panel and place the step ladder in front to view the transmitter in the top right corner.
4. Remove the tri-clamp from the sensor being calibrated.
5. Remove the sensor from the fitting and place on clean surface.
6. Install the 1 ½" tri-clamp cap, gasket and clamp on the open pipe connection.
7. Inspect the sensor for signs of damage.
8. Rinse one of the containers with deionized water and fill it with calibration solution. Just enough to submerge the 2" the sensor end.
9. Rinse the sensor with deionized water and give it shake to remove excess water.
10. Submerge the sensor tip in the calibration solution and move it around to get rid of the bubbles at the tip.
11. Access the conductivity calibration menu.




12. Enable the hold function so the transmitter will not alarm. This will have to be disabled after the calibration.



13. Set to transmitter to perform a manual calibration.




14. Verify the sensor is immersed and move around to make sure no bubble are trapped on the electrode. Allow the reading to stabilize.



The multiCELL displays the calibration result.

15. Enter the conductivity value of the calibration solution indicated on the bottle.
16. Press F4 (OK) to accept the calibration.
17. Disable hold function.



18. Remove the 1 ½" tri-clamp cap, gasket and clamp on the pipe connection.
19. Install the calibrated sensor, gasket and clamp in the open connection.
20. Dispose of calibration solution.
21. Fill in the calibration sticker with the date and date for next calibration along with your initials.
22. Repeat step 2 to 21 for the other two sensors.
23. See below:
 - a. Open valve V3230 if calibrating feed conductivity sensor AE3635.
 - b. Open valve V3260 if calibrating 1st pass conductivity sensor AE3665.
 - c. No valve to be closed of calibrating 2nd pass conductivity sensor AE3695.
24. Press the Mode Control followed by the Run to Drain Icon. Then allow the system to start-up.
25. Verify the water quality at each sample port using a calibrated hand held meter. If not with range repeat the calibration procedure.
26. Place the system in the mode it was in prior to calibration.

6.0 SYSTEM TROUBLESHOOTING

The following table refers to troubleshooting of BioPure HX2® System:

Tag#	Message	Possible Cause	Corrective Action
FALL 1807	Softener Interlock Alarm	Softener regeneration time set too close together.	Check and set timers so the regeneration times do not overlap.
		If only one softener is installed check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
		Softener setpoint incorrect.	Verify setpoint and correct if necessary.
FALL 1825	Carbon Interlock Alarm	Carbon backwash time set too close together.	Check and set timers so the regeneration times do not overlap.
		Check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
FALL 1845	Media Filter Interlock Alarm	Filter backwash time set too close together.	Check and set timers so the regeneration times do not overlap.
		If no filter or only one is installed check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
AOR 3615	Level Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
LALL 3615	Tank Level Low Low Alarm	Points of use valve left open or leaking.	Check all points of use and repair leaks if any found.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.

Tag#	Message	Possible Cause	Corrective Action
		Valves on system or piping connections are leaking.	Check system for leaks and repair if any are found.
		Level sensor is not working.	Visually check tank level and if good, check 4-20ma signal from sensor. If signal not within range replace sensor.
LAHH 3615	Tank Level High High Alarm	One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		Check valve is leaking.	Turn the system off and open tank. Check for signs that water is back feeding in the system. Determine source and repair.
		Level sensor is not working.	Visually check tank level and if good, check 4-20ma signal from sensor. If signal not within range replace sensor.
CAL 9874	Loop Return Flow Velocity Low Alert	Too much water is being consumed at points of use.	Check points of use consumption is within recommended design range.
		Leak in the loop or system.	Check for leaks and repair if any are found.
		Wrong pipe ID entered in HMI setpoint.	Verify loop internal diameter and enter it into the HMI.
AOR 3600	Pre-filter Inlet Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
AOR 3605	Pre-filter Outlet Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.

Tag#	Message	Possible Cause	Corrective Action
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
CAH 9850	Pre-filter Pressure Differential High Alert	Pre-filter cartridges are plugged.	Change filters.
		Pre-filter inlet or outlet pressure sensor is faulty.	Turn the system off and verify the pressure on each sensor is reading within limits. If one sensor is out by more than 5 psig replace sensor.
PALL 3605	Feed Pressure Low Low Alarm	Pre-filter cartridges are plugged.	Change filters.
		Facility feed water pressure is low.	Contact building services to correct issue.
		High pressure drops in pre-treatment.	Backwash or replace upstream filter cartridges.
		Pressure regulator has come out of adjustment or is damaged.	Adjust, repair or replace pressure regulator if required.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
TAHH 3610	Heater Element Temperature High High Alarm	A manual valve is closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		Flow restriction in the system.	Place the system in the mode it was in and visually verify there is flow going back to the tank. If little or no flow determine and remove the restriction.
AOR 3620		Terminal wire is loose.	Check wire terminals.

Tag#	Message	Possible Cause	Corrective Action
	Loop Return Flow Out Of Range Alarm	Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
AOR 3630	Loop Return Temperature Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
TALL 3630	Heat San Temperature Low Low Alarm	Heat SCR is not functioning correctly.	Check signal to SCR and replace if defective.
		Heater fuse is blown.	Check and replace.
		Low feed voltage.	Verify heater voltage and check with facility services.
		Low heater output power.	Check heater element resistance, if not within range replace heater.
		Loop not insulated or there is not enough insulation.	Inspect the loop insulation.
		Loop is too long.	Verify loop length and re-route to reduce the length.
AOR 3640	Feed Temperature Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
TAHH1 3640	Feed Temperature High High Alarm	Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Facility feed water temperature too high.	Check facility feed water temperature/quality. Contact building services to correct.
TAHH2 3640		A manual valve is closed.	Open manual valves that are supposed to be open.

Tag#	Message	Possible Cause	Corrective Action
	Heat Sanitization Temperature High High Alarm	One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		Flow restriction in the system.	Place the system in the mode it was in and visually verify there is flow going back to the tank. If little or no flow determine and remove the restriction.
KAHH 3640	Heat Sanitization Ramp Too Long Alarm	Timer setting is too low.	Increase timer setting.
		Heat SCR is not functioning correctly.	Check signal to SCR and replace if defective.
		Heater fuse is blown.	Check and replace.
		Low feed voltage.	Verify heater voltage and check with facility services.
		Low heater output power.	Check heater element resistance, if not within range replace heater.
		Loop not insulated or there is not enough insulation.	Inspect the loop insulation.
		Loop is too long.	Verify loop length and re-route to reduce the length.
AOR 3635	Feed Conductivity Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Feed water conductivity is out of sensor's range.	Check pre-treatment and flush to drain.
			Check softener regeneration timer. Repair softener and flush water to drain until feed water conductivity is within limits.
			Check facility feed water quality. Contact building services to correct.

Tag#	Message	Possible Cause	Corrective Action
AAH 3635	Feed Conductivity High Alert	Feed water conductivity is out of sensor's range.	Check pre-treatment and flush to drain.
			Check softener regeneration timer. Repair softener and flush water to drain until feed water conductivity is within limits.
			Check facility feed water quality. Contact building services to correct.
		Sensor is not reading correctly.	Verify water quality with hand held meter. If not correct perform sensor calibration.
UA 3015	1 st Pass Pump Drive Alarm	Look up error code in drive manual.	Correct issue based on drive error code.
AOR 3645	1 st Pass Pump Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
PAHH1 3645	1 st Pass Pump Pressure High High Alarm	A manual valve is closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		One or more automatic flow control valves are closed.	Check valve LED status versus the valve manual. Correct issue.
		HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Feed pressure spike.	Check feed water pressure regulator and repair or adjust if necessary.
		HMI pump pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor is damaged.	Replace sensor.

Tag#	Message	Possible Cause	Corrective Action
PAHH2 3645	1 st Pass Heat Sani Pressure High High Alarm	A manual valve is closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		One or more automatic flow control valves are closed.	Check valve LED status versus the valve manual. Correct issue.
		HMI valve pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor is damaged.	Replace sensor.
AOR 3650	1 st Pass Reject Pressure Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
AOR 3655	1 st Pass Reject Flow Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
AOR 3660	1 st Pass Reject Drain Flow Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.

Tag#	Message	Possible Cause	Corrective Action
FALL 3660	1st Pass Reject Drain Flow Low Low Alarm	A manual valve is closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		Reject to drain automatic flow control valve is not working properly.	Check valve LED status versus the valve manual. Correct issue.
		Reject to drain flow sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
CAH 9854	1st Pass Membrane Pressure Differential High Alert	Pre-treatment not functioning correctly.	Check water chemistry is within limits at RO pre-filter. Correct pre-treatment issue.
		Cleaning of the 1 st pass may be required if pressure differential has increased above 15% the initial pressure drop from 1 st time start-up or last cleaning.	Perform a low pH followed by a high pH cleaning.
		Reject automatic flow control valve is not working properly.	Check valve LED status versus the valve manual. Correct issue.
		1 st pass pump outlet or 1 st pass reject pressure sensor is faulty.	Check sensor operation and repair or replace.
		Adjustable pressure setpoint set too low.	Check setpoint and adjust if necessary. Setpoint should be 15% higher than the initial pressure differential at 1 st time start-up run to loop or after a cleaning has occurred.
AOR 3670	1st Pass Product Temperature Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.

Tag#	Message	Possible Cause	Corrective Action
AOR 3665	1st Pass Product Conductivity Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged. Feed water out of sensor's range.	Replace sensor. Check pre-treatment and flush to drain.
AAH 3665	1st Pass Product Conductivity High Alert	Cleaning of the 1 st pass required.	Perform a low pH followed by a high pH cleaning.
		Sensor is not reading correctly.	Verify water quality with hand held meter. If not correct perform sensor calibration.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Feed water out of sensor's range.	Check pre-treatment and flush to drain. Check facility feed water quality. Contact building services to correct.
UA 3025	2nd Pass Pump Drive Alarm	Look up error code in drive manual.	Correct issue based on drive error code.
AOR 3675	2nd Pass Pump Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm to see if the alarm goes away.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
PAHH1 3675	2nd Pass Pump Pressure High High Alarm	Manual valve closed.	Open manual valves that are supposed to be open.

Tag#	Message	Possible Cause	Corrective Action
		One or more automatic valves are in the wrong position.	Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve not opening or closing replace actuator. Check if LED on the corresponding air solenoid is on. If solenoid is supposed to be on replace solenoid.
		One or more automatic flow control valves are closed.	Check valve LED status versus the valve manual. Correct issue.
		HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		1 st pass pressure spike.	Check 1 st pass pump pressure setting or control and adjust if necessary.
		HMI pump pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor is damaged.	Replace sensor.
PAHH2 3675	2nd Pass Heat Sani Pressure High High Alarm	Manual valve closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve not opening or closing replace actuator. Check if LED on the corresponding air solenoid is on. If solenoid is supposed to be on replace solenoid.
		One or more automatic flow control valves are closed.	Check valve LED status versus the valve manual. Correct issue.
		HMI valve pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor is damaged.	Replace sensor.
		HMI pump pop-up configuration has been changed.	Verify pop-up configuration and reset.
AOR 3685	2nd Pass Recycle Pressure Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.
		Sensor is damaged.	Replace sensor.

Tag#	Message	Possible Cause	Corrective Action
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
CAH 9856	2nd Pass Membrane Pressure Differential High Alert	Pre-treatment not functioning correctly.	Check water chemistry is within limits at RO pre-filter. Correct pre-treatment issue.
		Cleaning of the 2 nd pass required. Pressure differential has increase above 15% the initial pressure drop from 1 st time start-up or last cleaning.	Perform a low pH followed by a high pH cleaning.
		2 nd pass recycle automatic flow control valve is not working properly.	Check valve LED status versus the valve manual. Correct issue.
		2 nd pass pump outlet or 2 nd pass reject pressure sensor faulty.	Check sensor operation and repair or replace.
		Adjustable pressure setpoint set too low.	Check setpoint and adjust if necessary. Setpoint should be 15% higher than the initial pressure differential at 1 st time start-up run to loop or after a cleaning has occurred.
AOR 3690	2nd Pass Recycle Flow Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
AOR 3700	2nd Pass Product Temperature Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
TALL 3700	Heat Sanitization Temperature Low Low Alarm	Heat SCR not functioning correctly.	Check signal to SCR and replace if defective.
		Heater fuse blown.	Check and replace.
		Low feed voltage.	Verify heater voltage and check with facility services.

Tag#	Message	Possible Cause	Corrective Action
		Low heater output power.	Check heater element resistance if not within range replace heater.
		Room temperature low.	Verify room temperature and increase.
		Fan or air duct blowing air across system.	Turn fan off or redirect air duct.
AOR 3695	2nd Pass Product Conductivity Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		1 st pass water out of sensor's range.	Check pre-treatment and flush to drain. Check facility feed water quality. Contact building services to correct.
AOR2 3695	2nd Pass Product Conductivity Out Of Range Alert For Longer Than 10 Minutes	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		1 st pass water out of sensor's range.	Check pre-treatment and flush to drain. Check facility feed water quality. Contact building services to correct.
AAH 3695	2nd Pass Product Conductivity High Alert	Cleaning of the 2 nd pass required.	Perform a low pH followed by a high pH cleaning.
		Sensor not reading correctly.	Verify water quality with hand held meter. If not correct perform calibration.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Feed or 1 st pass water out of sensor's range or poor quality.	Check pre-treatment and flush to drain. Check 1 st pass product.

Tag#	Message	Possible Cause	Corrective Action
			Check facility feed water or 1 st pass quality. Contact building services to correct or perform a low pH followed by a high pH cleaning.
		High CO ² level in feed water.	Check CO ² level. Install pH adjust system on feed or degas module on 1 st pass product line.
AAHH 3695	2nd Pass Product Conductivity High High Alarm	Cleaning of the 2 nd pass required.	Perform a low pH followed by a high pH cleaning.
		Sensor not reading correctly.	Verify water quality with hand held meter. If not correct perform calibration.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Feed or 1 st pass water out of sensor's range or poor quality.	Check pre-treatment and flush to drain. Check 1 st pass product
			Check facility feed water or 1 st pass quality. Contact building services to correct or perform a low pH followed by a high pH cleaning.
		High CO ² level in feed water.	Check CO ² level. Install pH adjust system on feed or degas module on 1 st pass product line.
AAHH2 3695	2nd Pass Product Conductivity High High Alarm For Longer Than 10 Minutes While Running To Drain	Cleaning of the 2 nd pass required.	Perform a low pH followed by a high pH cleaning.
		Sensor not reading correctly.	Verify water quality with hand held meter. If not correct perform calibration.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Feed or 1 st pass water out of sensor's range or poor quality.	Check pre-treatment and flush to drain. Check 1 st pass product
			Check facility feed water or 1 st pass quality. Contact building services to correct or perform a low pH followed by a high pH cleaning.

Tag#	Message	Possible Cause	Corrective Action
		High CO ² level in feed water.	Check CO ² level. Install pH adjust system on feed or degas module on 1 st pass product line.
AOR 3705	2nd Pass Product Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm to see if the alarm goes away.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
PAHH 3705	2nd Pass Product Pressure High High Alarm	A manual valve is closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		Product flow is too high for the loop inside diameter and or length.	Calculate loop pressure drop (including elbows, tees and elevation). If close to setpoint a large diameter loop will need to be installed or re-roughed to reduce the length.
		Feed water pressure too high.	Check feed water pressure and adjust, repair or replace regulator.
		HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Factory product flow setting too high.	Verify flow what flow rate was specified and change settings.
		2 nd pass pressure spike.	Check 2 nd pass pump pressure setting or control and adjust if necessary.
		HMI pump pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor reading incorrectly on HMI.	Verify sensor operation and change or calibrate if required.

Tag#	Message	Possible Cause	Corrective Action
		Sensor is damaged.	Replace sensor.
AOR 3710	2nd Pass Product Flow Out Of Range Alert	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
FALL 3710	2nd Pass Product Flow Low Low Alarm	A manual valve is closed.	Open manual valves that are supposed to be open.
		One or more automatic valves are in the wrong position.	Check if LED on the corresponding air solenoid is on. If off and solenoid is supposed to be on replace solenoid.
			Check automatic valve operation using the by-pass switch on the corresponding air solenoid. If valve is not opening or closing replace actuator.
		1 st or 2 nd pass pump control not working.	Verify pop-up configuration and reset.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Transmitter program is corrupt.	Verify and reprogram transmitter.
AOR 3715	Loop Filter Outlet Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Temperature shock occurred during heat sanitization.	Allow the system temperature to stabilize and reset alarm.
		Sensor is damaged.	Replace sensor.
		Filter option HMI selected and not installed.	Deselect from factory screen.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor reading incorrectly on HMI.	Check 4-20ma signal from sensor to PLC. If signal not within range replace sensor.
CAH 9852	Loop Filter Pressure Differential High Alert	Loop Filter cartridges are plugging.	Change filters.

Tag#	Message	Possible Cause	Corrective Action
		Filter inlet or outlet pressure sensor faulty.	Verify sensor operation and change or calibrate if required.
		Filter option HMI selected and not installed.	Deselect from factory screen.
AOR 6640	Compressed Air Pressure Out Of Range Alarm	Terminal wire is loose.	Check wire terminals.
		Cable wire is broken.	Check wire continuity and replace cable if necessary.
		Sensor is damaged.	Replace sensor.
		Sensor HMI pop-up configuration has been changed.	Verify pop-up configuration and reset.
		Sensor reading incorrectly on HMI.	Replace sensor.
PAL 6640	Compressed Air Pressure Low Alert	Air leak.	Check for leaks and repair.
		Compressor tank is full of water.	Purge compressor tank.
		Compressor outlet filter is plugged.	Check filter and replace.
		Compressor outlet filter is full of water.	Purge compressor filter.
		Compressor oil is low.	Check oil level and fill if necessary.
PALL1 6440	Compressed Air Pressure Low Low Alarm	Air leak.	Check for leaks and repair.
		Compressor tank is full of water.	Purge compressor tank.
		Compressor outlet filter is plugged.	Check filter and replace.
		Compressor outlet filter is full of water.	Purge compressor filter.
		Compressor oil is low.	Check oil level and fill if necessary.
PALL2 6440	Compressed Air Pressure Low Low Alarm	Air leak.	Check for leaks and repair.
		Compressor tank is full of water.	Purge compressor tank.
		Compressor outlet filter is plugged.	Check filter and replace.
		Compressor outlet filter is full of water.	Purge compressor filter.
		Compressor oil is low.	Check oil level and fill if necessary.
		Fuse blown.	Check and replace if necessary.
		Compressor pressure start /stop switch faulty.	Check switch and replace if faulty.
		Compressor faulty.	Check and replace if faulty.

Tag#	Message	Possible Cause	Corrective Action
CAL 9864	System Percent Rejection Low Alert	2 nd pass product water conductivity high.	Verify 2 nd pass product water. Perform a low pH followed by a high pH cleaning.
		Adjustable setpoint set too high.	Adjust setpoint to lower setting.
EAL 9800	Main Power Failure Alarm	Main disconnect fuses are blown.	Check fuse and replace.
		No power is going to panel main disconnect.	Check voltage at the main disconnect. If no voltage is detected contact facility services.
		Facility power outage.	Wait for power to be restored.
YSA 9805	Emergency Stop Alarm	Main panel red push button is pressed.	Consult with facility employees to verify who pushed the button. Reset the push button only if no emergency situation exists.
		Remotely wired red push button is pressed.	Consult with facility employees to verify who pushed the button. Reset the push button only if no emergency situation exists.
UA 9980	DI Alarm	DI tank is below quality setpoint.	Replace DI tank(s).
		If no DI used check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
UA 9985	DI Active Alarm	DI back up panel switch turned on.	Consult with facility employees why the switch is enabled. If DI back up is not required turn switch off.
		If no DI used check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
UA 9957	Water Leak Alarm	Water leak detector in alarm.	Verify if a water leak is active and fix leaks.
		If no leak detector installed check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
UA 9990	Auxiliary Alarm	Auxiliary device is in alarm.	Verify if device is in alarm and correct issue.
		If no auxiliary device installed check if panel quick connect jumper is removed.	Replace quick connect panel jumper.
	PLC Battery Alarm	Battery is low.	Replace battery.

7.0 SYSTEM STATUS MESSAGES

The following table refers to status messages on the screen of BioPure HX2® System.

TAG #	MESSAGE	MESSAGE TRIGGER
UA1800	Softener In Regeneration RO Is Off	Softener In Regeneration
	Softener In Regeneration RO Is On	Softener In Regeneration
UA1805	Softener In Regeneration RO Is Off	Softener In Regeneration
	Softener In Regeneration RO Is On	Softener In Regeneration
UA1810	Carbon In Backwash RO Is Off	Carbon In Backwash
	Carbon In Backwash RO Is On	Carbon In Backwash
UA1815	Carbon In Backwash RO Is Off	Carbon In Backwash
	Carbon In Backwash RO Is On	Carbon In Backwash
UA1820	Carbon In Backwash RO Is Off	Carbon In Backwash
	Carbon In Backwash RO Is On	Carbon In Backwash
UA1835	Filter In Backwash RO Is Off	Filter In Backwash
	Filter In Backwash RO Is On	Filter In Backwash
UA1840	Filter In Backwash RO Is Off	Filter In Backwash
	Filter In Backwash RO Is On	Filter In Backwash
UA1850	Pre-treatment Offline RO is Off	Pre-treatment Offline
	Pre-treatment Offline RO is On	Pre-treatment Offline
	1st Pass Heat Sanitization In Progress. Temperature XX.X °C/°F	1 st Pass Heat Sanitization in progress.
	2nd Pass Heat Sanitization In Progress. Temperature XX.X °C/°F	2 nd Pass Heat Sanitization in progress.
	Loop Heat Sanitization In Progress. Temperature XX.X °C/°F	Loop Heat Sanitization in progress.
	Loop POU Heat Sanitization In Progress. Temperature XX.X °C/°F	Loop POU Heat Sanitization in progress.
	1st Pass Heat Sanitization Completed. XXX Minutes at XX.X °C Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	1 st Pass Heat Sanitization Completed Successfully.
	2nd Pass Heat Sanitization Completed. XXX Minutes at XX.X °C Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	2 nd Pass Heat Sanitization Completed Successfully.

TAG #	MESSAGE	MESSAGE TRIGGER
	Loop Heat Sanitization Completed. XXX Minutes at XX.X °C Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	Loop Heat Sanitization Completed Successfully.
	Loop POU Heat Sanitization Completed. XXX Minutes at XX.X °C Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	Loop and Point Of Use Heat Sanitization Completed Successfully.
	1 st Pass Chemical Clean Complete. Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	1 st Pass Chemical Clean Completed Successfully.
	2 nd Pass Chemical Clean Complete. Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	2 nd Pass Chemical Clean Completed Successfully.
	Loop Chemical Clean Complete. Started HH:MM:SS DD/MM/YYYY Completed HH:MM:SS DD/MM/YYYY	Loop Chemical Clean Completed Successfully.
	1st Pass Heat Sanitization Aborted. Started HH:MM:SS DD/MM/YYYY Aborted HH:MM:SS DD/MM/YYYY	1 st Pass Heat Sani Aborted.
	2nd Pass Heat Sanitization Aborted. Started HH:MM:SS DD/MM/YYYY Aborted HH:MM:SS DD/MM/YYYY	2 nd Pass Heat Sani Aborted.
	Loop Heat Sanitization Aborted. Started HH:MM:SS DD/MM/YYYY Aborted HH:MM:SS DD/MM/YYYY	Loop Heat Sani Aborted.
	Loop POU Heat Sanitization Aborted. Started HH:MM:SS DD/MM/YYYY Aborted HH:MM:SS DD/MM/YYYY	Loop and Points Of Use Heat Sani Aborted.
	Heat Sani System Must Be In Auto Or Standby To Start Sani.	Heat Sani Was Triggered Automatically Or Manually And the System Was Not In Auto or Standby.
	PLC Battery Low. Time XX:XX, Date DD/MM/YYYY	PLC Sent A Message To The HMI Indicating The PLC Battery Is Low And Needs To Be Changed.
	Stop Alarm Disabled. Time XX:XX, Date DD/MM/YYYY	Alert Triggered When User disables Stop On The Setup Screen.
	Stop Alarms Enabled. Time XX:XX, Date DD/MM/YYYY	Alert triggered when operator enables stop alarms on the setup screen.

TAG #	MESSAGE	MESSAGE TRIGGER
	Heat Sani Paused Due To Pre-treatment Regeneration or Backwash. Time XX:XX, Date DD/MM/YYYY	A Pre-treatment Device Is In Regeneration Or Backwash. The Heat Sani Has Paused Until The Device Is Finished It's Cycle.
	Scheduled Heat Sani Could Not Be Completed.	Heat Sani Could Not Be Completed Due To An Alarm Or User Interruption.
	DST MANUAL PLUS ONE HOUR	User presses the "+ Plus 1 hour" push button on the HMI "Time and Date" screen.
	DST MANUAL MINUS ONE HOUR	User presses the "- Minus 1 hour" push button on the HMI "Time and Date" screen.
	DST AUTO PLUS ONE HOUR	The time is automatically changed if the "Automatic Daylight Savings time Enabled" is set to yes, and the current time is greater than the "Set DST Spring (+ Plus 1 hour) and less than "Set DST Fall (- Minus 1 hour).
	DST AUTO MINUS ONE HOUR	The time is automatically changed if the "Automatic Daylight Savings time Enabled" is set to yes, and the current time is less than the "Set DST Spring (+ Plus 1 hour) or greater than "Set DST Fall (- Minus 1 hour).
	WATER SAMPLE REQUIRED	15 minutes after entering modes 8, 11 or 13 and setpoint is enabled.

8.0 SPARE PARTS LIST

When ordering parts for the system please have the Biolab Serial # (located on right side of the main panel) and the tag number (on component) of the component you require.

NOTE: *Some components require programming prior to shipment and are identified with a * beside the tag number.*

BioPure HX2® Spare Parts List ⁴	
Part #	Part Description
AJT150	EC/FUSE/600V 150A CLASS J TD FUSE
AJT35	EC/FUSE/600V 35A CLASS J TD FUSE
AJT4	EC/FUSE/600V 4A CLASS J TD FUSE
AJT40	EC/FUSE/600V 40A CLASS J TD FUSE
AJT60	EC/FUSE/600V 60A CLASS J TD FUSE
ATDR1	EC/FUSE/600V 1A CC TD FUSE
ATDR15	EC/FUSE/600V 15A CC TD FUSE
ATDR20	EC/FUSE/600V 20A CC TD FUSE
ATDR25	EC/FUSE/600V 25A CC TD FUSE
ATDR4	EC/FUSE/600V 4A CC TD FUSE
ATDR7	EC/FUSE/600V 7A CC TD FUSE
ATDR8	EC/FUSE/600V 8A CC TD FUSE
DFJ35	EC/FUSE/460V 35A DFJ CLASS J FUSE
DFJ70	EC/FUSE/208V 70A DFJ CLASS J FUSE
GMD1	EC/FUSE/250V 1A GGA TD FUSE 5/PK
GMD2	EC/FUSE/250V 2A GGA TD FUSE 5/PK
GMD3	EC/FUSE/250V 3A GGA TD FUSE 5/PK
GMD4	EC/FUSE/250V 4A GGA TD FUSE 5/PK
GMD500	EC/FUSE/250V 0.5A GGA TD FUSE 5/PK
LPCC1	EC/FUSE/TIME DELAY CLASS CC, 1A
LPCC20	EC/FUSE/TIME DELAY CLASS CC, 20A
LPCC9	EC/FUSE/TIME DELAY CLASS CC, 9A
LPJ30SP	EC/FUSE/TIME DELAY CLASS J 30A
LPJ40SP	EC/FUSE/TIME DELAY CLASS J 40A
LPJ90SP	EC/FUSE/TIME DELAY CLASS J 4A
LPJ90SP	EC/FUSE/TIME DELAY CLASS J 90A
98887181	CRN10-10 PUMP WITH 7.5 HP 208V 460 VAC MOTOR
99051987	CRN3-25 PUMP WITH 5.0 HP 208V 460 VAC MOTOR
96315812	CRN10-10 PUMP WITH 7.5 HP 575 VAC MOTOR
96315813	CRN3-25 PUMP WITH 5.0 HP 575 VAC MOTOR
96501053	CRN10-18 PUMP WITH 10 HP 380-415 VAC MOTOR
96513493	CRN5-29 PUMP WITH 5.5 HP 380-415 VAC MOTOR
7000302A020SP	CARTRIDGE/PLT/20IN/0.2MICRON/222 FIN
FM00120X	20" 1 MICRON MICROFIBER FILTER DOE
SANROHS28	8" X 40" HOT MEMBRANE THIN FILM

⁴ Not all parts will be used on all systems, contact your local Mar Cor Sales representation for clarification regarding spare parts and your specific model

BioPure HX2® Spare Parts List ⁴	
Part #	Part Description
466040	GLOBE VALVE 316L/3/4" FAIL OPEN TC ENDS
556575	TYPE 8012/FLOW SENSOR/1.5"/SS/20 RA, MT
556613	FLOW SENSOR/TYPE 8012/3/4"TC/SS BODY/EP
557065	FLOW SENSOR/TYPE 8012/4-20MA 1" AND UP
560766	TYPE 8012/FLOW SENSOR/1"/SS/20 RA, MT
563419	METER/CON,RES/8619/3 CONDO INPUTS, 2 OUT
992180	1-1/2" TC SS RTD TEMP TRANSMITTER 0-200
7093122	S-10 PRESS TRANSMITTER 0-100 WC 1/4" MNPT
13299701	1 1/2" TRI-CLAMP 0-100PSIG SANITARY PRESSURE TRANSM, M12, CAL
13333127	1/4" MNPT, 316 SS 0-100 PSIG, 4-20MA PRESSURE SENSOR, 100C
13333127	A-10 PRESS TRANSMITTER 0-100 PSIG 1/4" MNPT
13335384	1/4" MNPT, 316 SS 0-500 PSIG, 4-20MA PRESSURE SENSOR, 100C, A-10
13335384	A-10 PRESS TRANSMITTER 0-500 PSIG 1/4" MNPT
14020300	1-1/2" TC SA-11 PRESS TRANSMITTER 0-100PSI 1/4"
92709095	CONDUCTIVITY PROBE, 1.5"TC 316SS, 0.1K
92709773	COND PROBE, 1.5" TC 316SS, 0.01K 20'CABLE
92709775	CONDUCTIVITY SENSOR W/ATC 1.5" TC, K=1 C
TWR5202100	3/4" MNPT TEMPERATURE SWITCH NC 100 C
CXS620F332R001	IMMERSION HEATER 208V 60HZ 20KW CERT
CXS620F332R002	IMMERSION HEATER 460V 60HZ 20KW CERT
CXS620F332R005	IMMERSION HEATER 380V 50HZ 20KW NO CERT
CXS620F332R007	IMMERSION HEATER 415V 50HZ 20KW NO CERT
00186260	TYPE 8640/3/2WAY/SOL VALV/24VDC
00241627	VALVE BLOCK 8640 11 BANK 10 VALVE 6524
00278773	VAL / ANG SEAT / 1/2" OD TUBE /DN15 / 4-20MA CON
00278774	VAL / ANG SEAT / 1-1/2" OD TUBE /DN40 / 4-20MA CON
00278775	VAL / ANG SEAT / 1" OD TUBE / 4-20MA CON
00456443	VALVE/SV/1/2" FNPT/6213/NC/24VDC
00458732	6213 SOL VALVE/2/2 WAY/1/4" FNPT/NC/24VAC
458776	VALVE/6213/316SS/1/4" FPNT/NC/24VDC
3WAY112CCWSR	3 WAY VALVE 1.5" TC WITH ACTUATOR
1701502271154J	3WAY112CCWSR SEAL REPAIR KIT
3WAY1CCWSR	3 WAY VALVE 1" TC WITH ACTUATOR
3WAY1CWSR	3 WAY VALVE 1" TC WITH ACTUATOR
1701002271154J	3WAY1CCWSR / 3WAY1CWST SEAL REPAIR KIT
BV112FC	BALL VALVE 1.5" TC FC WITH ACTUATOR
14015022703578	1.5" BALL VALVE SEAL KIT
BV1FC	BALL VALVE 1" TC NC WITH ACTUATOR
MBV1	MANUAL BALL VALVE 1" TC
MBV112	MANUAL BALL VALVE 1.5" TC
14010022703578	1" BALL VALVE SEAL KIT
OFSV90F1	VALVE/PV,90/316/1"-1.5" TC/90/25RA/EP
OFSV90F34	VALVE/PV,90/316/3/4" TC/90/25RA/EP/COC/MT
OFSV90RK	VALVE REPAIR KIT, 90 DEGREE, REPAIRS (2)
OFSVILF	1/2"TCx3/8"HB SANITARY SAMPLE VALVE 316L

BioPure HX2® Spare Parts List ⁴	
Part #	Part Description
15280	AIR COMPRESSOR/SYNTHETIC OIL ⁵
C0164N6	AIR COMPRESSOR REGULATOR FILTER ELEMENT
C0222	MOTOR/COMPRESSOR/AIR/220VAC/60HZ/P50
C1794	AIR COMPRESSOR INTAKE FILTER
M03092SM2	AIR COMPRESSOR/PT50TCAL/2.1CFM 120PSI MA
114.77EGALV	1-1/4" GALVANIZED VICTAULIC CLAMP EPDM
2.77EGALV	2" GALVANIZED VICTAULIC CLAMP EPDM
R13HC112	1" & 1-1/2" 304 SS SINGLE PIN CLAMP
R13HC2	2" 304 SS SINGLE PIN CLAMP
R13HC34	1/2" & 3/4" 304 SS SINGLE PIN CLAMP
R13HP112	1" & 1-1/2" 304 SS 2 PIN HP CLAMP
40MPE100	1" EPDM TRI CLAMP GASKET
40MPE150	1-1/2" EPDM TRI CLAMP GASKET
40MPE200	2" EPDM TRI CLAMP GASKET
42MPGTS050	GASKET/TUFFSTEEL/1/2" TC
42MPGTS075	3/4" TUFFSTEEL TRI CLAMP GASKET
854030E	3" FLANGE GASKET 150# EPDM
780779	SILENCER, PE, threaded G 1/8
901108	1/4" TUBING X G1/8 ANGLED PNEUMATIC CONNECTOR
EVC528	M12 5 PIN X FLYING LEAD CABLE 10 METER
BPHX2NURSESTN	NURSES STATION KIT FOR HX 2

⁵ Available in Continental America only, if required outside of this area please contact your local Mar Cor representative.

9.0 GENERAL LIMITED WARRANTY

Mar Cor Purification and Biolab Equipment Ltd. ("Mar Cor®") warrants all products manufactured by Mar Cor® to be free from defects in material and workmanship under normal use for a period of one year from date of shipment. Warranty is only provided when the equipment is properly maintained in accordance with Mar Cor® instructions. Our obligation shall be limited to the repair or replacement at our option of any parts we deem defective. Consequential damages due to the failure of any part or system are not provided. Mar Cor® provides no other warranty, expressed or implied. There is no warranty of merchantability or fitness for a particular purpose. Other manufacturer's equipment offered, as part of any proposal shall carry the manufacturer's warranty. On site, labour is covered for the first 90 days. Travel time and expenses are not covered.

9.1 CONTACT INFORMATION

For further information regarding the operation or maintenance of the Mar Cor® water system and for ordering, please contact your local Mar Cor Purification office. When ordering replacement parts, please provide the Biolab Serial # located on the right side of the control panel.

CANADA

Manufactured by: Biolab Equipment Ltd o/a Mar Cor Purification.
3250 Harvester Rd - Unit 6
Burlington, ON L7N 3W9
Tel: (905) 639-7025
Fax: (905) 639-0425

Toll Free: (800) 268-5035

UNITED STATES OF AMERICA

Distributed By: Mar Cor Purification
4450 Township Line Rd
Skippack, PA 19474-1429
Tel: (484) 991-0220
Fax: (484) 991-0230
info@mcpur.com

Toll Free: (800) 346-0365

WORLD WIDE WEB

www.mcpur.com

This Page Has Intentionally Been Left Blank