

Mar Cor Purification  
introduces the

# Millenium HX<sup>®</sup>



Portable Reverse Osmosis System with **Heat Disinfection**.  
For Use in Home and Acute Settings.



# Millenium HX Water Purification System

## The Next Generation of Portable Reverse Osmosis

Mar Cor Purification offers the next generation portable RO unit, the Millenium HX. This compact automatic RO represents a quantum leap in reliable water purification production with the added feature of automatic hot water disinfection.

The Millenium HX is designed to meet AAMI grade water standards according to ANSI/AAMI 13959 and holds an FDA cleared claim for use in home and acute dialysis settings. The automated heat disinfect process allows for greater bacteria control, reducing the need for chemicals, which leads to increased peace of mind and less downtime.

### What Makes the Millenium HX Different ?

- **Complete System Disinfection** - Unit has the ability to disinfect the RO and the product water supply tubing to the dialysis machine.
- **Quiet Operation** - Low noise level for regular and nocturnal treatments.
- **Home Use** - The only portable Reverse Osmosis System holds an FDA cleared claim for home use.
- **Internet Connectable** - Data logging of important events along with password protected set-point & interactive troubleshooting screens.
- **Safety** - Programmable auto-flushing and continuous product loop eliminates stagnation and ensures system cleanliness and readiness.
- **User Friendly** - Color touch screen LCD control panel offers display of operational functions & key parameters.
- **Isolated Disinfection** - Unit can complete a hot water disinfection disconnected from water source.



Mar Cor Purification  
4450 Township Line Road  
Skipack, PA 19474-1429  
Tel: (484) 991-0220  
Toll Free: (800) 346-0365  
Fax: (484) 991-0230

Mar Cor Purification  
14550 28th Avenue North  
Plymouth, MN 55447  
Tel: (484) 991-0220  
Toll Free: (800) 633-3080  
Fax: (763) 210-3868

Mar Cor Purification  
3250 Harvester Road - Unit 6  
Burlington, ON L7N 3W9  
Tel: (905) 639-7025  
Toll Free: (800) 268-5035  
Fax: (905) 639-0425

Mar Cor Purification  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: (+31) 45 5471 471  
Fax: (+31) 45 5429 695

Mar Cor Purification  
1A International  
Business Park, #05-01  
Singapore 609933  
Tel: (+65) 6227 9698  
Fax: (+65) 6225 6848