

PRODUCT BULLETIN #3240

Date: September 11, 2023

From: Michael Verguldi, Sr. Product Manager

Subject: **Notification of End of Life for Posiclear® Filter**

Effectivity: End of Sales: September 11, 2023
End of Shipment: October 31, 2023

Evoqua Water Technologies, which acquired Mar Cor® Water Purification Systems and Services in January 2022, is providing notification of the discontinuance of the Posiclear Filter line. The Posiclear Filter is used exclusively as a final filter to remove endotoxins and bacteria as part of the CWP 100 WRO H/S and BioPure HX2® water treatment systems for hemodialysis.

With the discontinuation of the Posiclear Filters, Evoqua has developed the RenaPure® Filter a direct replacement for new and existing CWP 100 WRO H/S and BioPure HX2 water treatment systems for hemodialysis. It is designed to fit into existing housing configurations and remains an effective safeguard against potential endotoxins and bacteria. This final filter is available immediately and the corresponding part numbers are listed in the chart below. Please direct all questions to your local sales representative or to info@mcpur.com.

Evoqua P/N	Mar Cor P/N	Description	Replacement P/N
W2T913258	WT7002502A020SPD	Filter, Posiclear, 20", 222, Flat	W3T590100
W2T913257	WT7000302A020SPD	Filter, Posiclear, 20", 222, Fin	W3T590099

The Posiclear Filters will no longer be offered for sale as of September 11, 2023. However, due to demand, supply shortages, and manufacturing constraints, these dates are subject to change and an existing order may need to be modified.

Please direct all questions to your local sales representative or to info@mcpur.com.

BioPure HX2, Mar Cor and RenaPure are trademarks of Evoqua Water Technologies LLC, its subsidiaries or affiliates in some countries. Posiclear is a trademark of Gambro Lundia in some countries and used under license.