

ROVER ALX[®]

Dialysis Water Transportation System



MAR COR

A BETTER WAY TO MOVE WATER.

The ROVER ALX dialysis water transport system is our smallest transport system footprint designed to fit the tightest spaces without sacrificing safety. Unparalleled in today's market, it streamlines the movement of water dialysis treatment equipment and replaces many outdated, risk-prone carts currently used in U.S. hospital settings, nursing care facilities as well as home care spaces.

RELIABLE AND EFFICIENT.

The ROVER ALX transport system is ready when your patient needs it. This compact, open-frame transport system can carry everything required to perform an effective hemodialysis patient treatment and still fit in a crowded patient care area. The ROVER ALX transport system dramatically improves the process of moving and storing a hemodialysis water treatment system.

COMPACT ALL-IN-ONE WATER TREATMENT DELIVERY SYSTEM

An all-in-one, compact water transport system that holds a portable Reverse Osmosis (RO) device and the necessary pretreatment components to deliver proven hemodialysis water treatment technology to patients – at an affordable price point. The ROVER® ALX transport system is built of corrosion-resistant materials that can withstand the daily demands of a wet environment and eliminates the concern of metal deterioration and rust often seen on other carts. The monitoring station and sample ports are clearly labeled and perfectly appointed for easy access. Additionally, the transport system has completed and passed stability and threshold testing to ensure a safe experience for the user.

STANDARD FEATURES & BENEFITS.

- Rust-proof and corrosion resistant design
- Open concept for easy access and monitoring
- Dual wheel casters for a smooth ride with an integrated braking mechanism
- Sample ports with removable drip tray align with integrated pressure gauges
- Easy to clean and disinfect surfaces to help avoid contamination and lower patient risk
- Complies with IEC Tip and Threshold Testing Standards
- Complete with a handy start-up and shut-down reference guide



UNIVERSAL DESIGN

Compatible with a variety of portable hemodialysis RO systems, including: **EON®, MILLENIUM HX®, WRO 300 H, WRO 300**



SMALL FOOTPRINT

The smaller size of the ROVER ALX makes it a great solution for acute treatments in hospital settings, nursing care facilities or home care spaces.



USER FRIENDLY

Featuring a sturdy bracket for cord and hose management, as well as easily accessible sample ports with drip tray.



CORROSION RESISTANT

Constructed of Aluminum, the ROVER ALX system is rust proof and easy to clean.



Visit mcpur.com for more information or call 1-800-633-3080

EON, Millenium HX, and ROVER ALX are trademarks of Evoqua Water Technologies LLC, its subsidiaries or affiliates in some countries.

All information presented herein is believed reliable and in accordance with accepted engineering practices. Evoqua makes no warranties as to the completeness of this information. Users are responsible for evaluating individual product suitability for specific applications. Evoqua assumes no liability whatsoever for any special, indirect or consequential damages arising from the sale, resale or misuse of its products.

© 2023 Evoqua Water Technologies LLC Subject to change without notice P/N:W3T595391 Rev.A