

### PRODUCT BULLETIN #3248

Date: March 29, 2024

From: Michael Verguldi, Sr. Product Manager

Subject: **Availability of Minncare® Cold Sterilant and Minncare HD Disinfectant in Canada**

Effectivity: Immediately

This Product Bulletin is to inform the marketplace that Minncare Cold Sterilant and associated test strips are no longer available having been discontinued by STERIS Life Sciences.

This Product Bulletin is also to inform the marketplace that Minncare HD Disinfectant and associated test strips are only available to the hemodialysis marketplace and are not available for sale outside the United States unless prior approval has been obtained. This encompasses the entire product line, including various sizes as well as test strips.

Minncare HD Disinfectant is a clear liquid disinfectant solution that consists of a stabilized mixture of hydrogen peroxide, peracetic acid and acetic acid. Minncare HD Disinfectant is intended for hemodialysis water purification system disinfection when diluted to a 1% concentration (1 part Minncare HD to 99 parts water) and applied for a minimum contact time of 36 minutes at 20°C.

If operating a Mar Cor® Water Purification System in Canada and a chemical disinfection process is required or becomes necessary, be advised of the following recommendation:

<u>Country</u>	<u>Product</u>	<u>Details</u>
Canada	PROCLEAN PC 110	5-Gallon pail; P/N W2T940939
	PROCLEAN PC 110TS-MED	50-Test Strip Pack; P/N W2T941257
	PROCLEAN PC 110TS-HI	50-Test Strip Pack; P/N W2T941258

Please note that the chemistry of the afore mentioned product may be different from Minncare products, making specifications, usage guidelines, and other relevant details different. Likewise, Minncare test strips are NOT to be used with any other chemical disinfectant product. Please be sure to consult with the manufacturer's shipping and handling instructions, SDS, and recommended sanitization concentration levels suitable for both RO membranes as well as for process piping and components.

Mar Cor is a trademark of Evoqua Water Technologies, LLC its subsidiaries, and affiliates in some countries.  
Minncare is a trademark of Medivators, used under license.