


SDS/PSDS PLASTIC PUMP STAND

This Tech Note identifies a product improvement as it relates to the metal pump stand (i.e., base) that is directly mounted to the mix pump, specifically used on the Solution Delivery System (SDS) and the Pressurized Solution Delivery System (PSDS) unit(s).

We have produced a plastic base for the mix pump/motor set-up. Units built as of September 1, 2016, will have the new plastic pump base installed as standard. Clinics who prefer or require a replacement pump stand on an existing SDS or PSDS unit due to corrosion or for aesthetic reasons, can order a pump upgrade kit, P/N W3T579009. Included in the kit is a set of instructions as well as all of the necessary hardware.



This is a direct replacement; no modifications are required. Please note that a new pump/motor is NOT included in the kit. If you have any questions, please contact Technical Service at 888-595-0666.