

# WRO 300 H

PORTABLE REVERSE  
OSMOSIS SYSTEM WITH  
AUTOMATED HOT WATER  
DISINFECTION



## MARCOR

### CLEAN BY DESIGN

The WRO 300 H is designed to overcome the problem of microbiological build up in portable dialysis reverse osmosis systems that may have extended idle periods. By integrating an automated heat disinfection system into the unit, the WRO 300 H system is engineered and manufactured to be clean by design.

The WRO 300 H water purification RO unit is designed for single hemodialysis patient operation. It satisfies the stringent demands of acute providers, while addressing the typical concerns of the home patient. This device meets AAMI grade water standards with the added capability of efficiently disinfecting without the use of chemicals, all in a validated format.

# WRO 300 H WATER PURIFICATION SYSTEM

## WHAT MAKES THE WRO 300 H DIFFERENT?

- **Ease of Use**  
Automated procedures to minimize caregiver involvement
- **Quiet Operation**  
Low noise level for regular and nocturnal treatments
- **Low Maintenance**  
Validated, Automated Hot Water or Chemical disinfection for bacteria control
- **Safety**  
Programmable auto-flushing and continuous product loop eliminates stagnation and ensures system cleanliness and readiness
- **Reliability**  
Consistent production of quality water on demand
- **Service & Support**  
Technical assistance, manuals, training DVD for help when needed



Visit [mcpur.com](http://mcpur.com) for more information or call 1-800-633-3080

All information presented herein is believed reliable and in accordance with accepted engineering practices. Evoqua makes no warranties as to the completeness of this information. Users are responsible for evaluating individual product suitability for specific applications. Evoqua assumes no liability whatsoever for any special, indirect or consequential damages arising from the sale, resale or misuse of its products.