

### DISCONTINUED ELEMENT INFORMATION (MEDICAL DEVICES)

The following elements have been discontinued by the manufacturer. We have optional replacements that will require purchase of a different style element and/or a complete ROMA (Reverse Osmosis Membrane Assembly). This generally includes the housing, element, end caps, and fittings. In some cases, additional fittings and tubing will be required (i.e., 700 Series RO units).

#### CAUTION

If switching from a CA (Cellulose Acetate) material to the more common TFC (Thin Film Composite) material, please note that chlorine removal must be performed before the RO element or damage to the element(s) will occur.

Discontinued Part Number	Replacement Part Number	Description of Change
ME40179A	W3T579381	Tape wrapped element with extended brine seal to fit in full size housings.
ME40599A	3024436 - No longer available	Complete ROMA and required parts. Required on first converting from old style housing. For element replacement after initial installation, element ME40599 (no longer available) must be purchased.
ME40091A and ME40510A	W3T578353 ***Additional fittings may be required depending on RO, consult Technical Support.	Converts from CA element material (ME40091A only) and housing(s) to current TFC element and housing(s). For element replacement after initial installation, element W2T913222 must be purchased.

### DISCONTINUED ELEMENT INFORMATION (NON-MEDICAL DEVICES)

The following elements have been discontinued by the manufacturer with no replacement. Therefore, these are discontinued with no direct replacement available. Please contact Tech Support for assistance in converting to a potential compatible element or to a newer device.

#### Discontinued Part Numbers

475-02-001	112 HR(CA) element	475-02-004	112 HR(CA) element
475-02-051	192 HR element	475-02-054	192 HR(PA) element
475-04-054	411 HF(PA) element	475-04-154	416 HR element
475-04-151	416 HR(PA) element	475-04-001	704 HR(CA) element