

**WRO 300 / 300 H PART CHANGES**

This is an update to what was previously WATER,SP,10091-12/2009.

There have been changes made to the parts listed below:

1. The SS Y-connector for the WRO product line has in the past been prone to leaks. The valve part of the connector was re-designed to help prevent these failures. The valve part and O-ring between the valve and body are now available as spare parts.
  - a. Valve connector part number is W2T914795.
  - b. Connector O-ring part number is W2T914567.
2. The float switches for the inlet tank are now only available to order in quantity of one. In the past these parts were shipped two per package.
  - a. The new part number is W3T572520.
  - b. The old part number WTW508013001 (package of two) is no longer available.

