

BIOPURE HX2® DRIVE MOTOR REPLACEMENT

As part of our commitment to continual product improvement, we have updated the drive motors on the BioPure HX2 reverse osmosis machine. These new drive motors are equally reliable and will provide improved efficiency while still delivering proper operation. This change starts with serial number BM16-067. However, the new drive motors are equivalent and can be used on existing BioPure HX2 devices as direct replacement on those machines that have the older drive motors.

Ordering Information

For ordering purposes, the following table identifies the replacement part numbers as required:

NEW P/N	OLD P/N	NEW DESCRIPTION	VOLTAGE
W2T921212	3032838	DRIVE,MOTOR,208VAC,3.7KW,17.5A,3PH,5HP	208V/60HZ
W2T921211	3032839	DRIVE,MOTOR,208VAC,5.5KW,24A,3PH,7.5HP	208V/60HZ
W2T921210	3032840	DRIVE,MOTOR,460VAC,3.7KW,8.7A,3PH,5HP	460V/60HZ
W2T921209	3032841	DRIVE,MOTOR,460VAC,5.5KW,13A,3PH,7.5HP	460V/60HZ
W2T921207	3032844	DRIVE,MOTOR,575VAC,4KW,6.6A,3PH,5HP	575V/60HZ
W2T921208	3032845	DRIVE,MOTOR,575VAC,5.5KW,9.9A,3PH,7.5HP	575V/60HZ

BioPure HX2 is a trademark of Evoqua Water Technologies LLC, its subsidiaries or affiliates in some countries.