

BIOPURE HX2® PREFILTER

Due to changes in the maximum temperature rating of the existing prefilter used with the BioPure HX2 RO System, part number W2T912567 (275-10-596), we are now referencing a different depth filter cartridge. The replacement depth filter cartridge is part number W2T912633. The specifications for this filter are as follows:

Description

String wound polypropylene filter cartridge

Materials of Construction

Filter Media:	Polypropylene
Core Type:	Polypropylene
Nominal Length:	20 inches
Outside Diameter:	2.5 inches
Inside Diameter:	1 inch
End Style:	Double Open End, no gasket

Performance

Micron Rating:	1 micron
Max. Temperature Rating:	82°C (180°F)
Max. Operation Differential:	60 psid (4 bar) at 21°C (70°F)
Clean Flow Rate:	2 psid drop @ 5 gpm per 20" length

The replacement depth filter is equivalent in fit, form and function. The only change is the higher maximum operating temperature capacity which is normally only relevant during heat sanitization mode of operation. All other BioPure HX2 system operational parameters are not impacted by this replacement prefilter. System monitoring and replacement frequency should remain the same.

BioPure HX2 is a trademark of Evoqua Water Technologies, LLC., its subsidiaries, and affiliates in some countries.